

## ProCare Dent 40

### Rinsing agent

#### Use:

Liquid concentrate. For the automated reprocessing of dental instruments including transmission instruments in Miele Professional washer-disinfectors.

#### Quality features:

*ProCare Dent 40* complies with European guidelines for medical devices. **CE MD**

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

#### Product properties:


- Improves and accelerates the drying of dental instruments in machine-based reprocessing methods.
- Neutralises alkaline carryover and alkalinity.

#### Use and dispensing:

*ProCare Dent 40* is formulated for use in processes in Miele Professional washer-disinfectors.

- Suitable dosing dispensers should be used.
- Application concentration: 0.0375 fl oz/gal / 0.3 ml/l via door dispenser for rinsing agents
- Programme selection: Vario TD
- Minimum dispensing temperature:  $\geq 158$  °F / 70 °C

#### Programme sequence:



#### Vario TD process:

- 1 Pre-wash**  
Cold water (softened), no heating
- 2 Main wash**  
Cold / hot water (softened), with heating, with process chemicals
- 3 Neutralisation**  
Cold / hot water (softened), no heating, with process chemicals
- 4 Intermediate rinse**  
Cold / hot water (softened), no heating
- 5 Thermal disinfection**  
at  $> 194$ °F / 90°C with 5-minute holding time in the final rinse with *ProCare Dent 40*
- 6 Drying**

Miele recommends the use of fully demineralised water in the final rinse in order to protect the load and the washer-disinfectors against stains and deposits and to protect anodised aluminium.

#### Please note!

- Do not mix with other products.
- Before switching products, rinse the dispensing system, including intake hoses, with water.
- Reprocessing should comply with all ordinances pursuant to medical device regulations and should be performed with appropriate validated processes.
- Observe recommendations on reprocessing issued by medical device manufacturers with respect to the provisions of DIN EN ISO 17664 (CAN/CSA-Z17664).
- Operating instructions for Miele Professional washer-disinfectors must be observed.
- For commercial applications only.

## ProCare Dent 40

### Rinsing agent


**Technical data:**

- Density (at 68°F / 20 °C): approx. 1.1 g/cm<sup>3</sup>
- pH value (determined using fully demineralised water, at 68°F / 20 °C) 0.3 ml/l: approx. 3.8

**Ingredients according to Regulation (EC) No 648/2004 on detergents:**

- < 5 % phosphonates
- 15 % to 30 % non-ionic tensides
- Also preservatives (methylchloroisothiazolinone, methylisothiazolinone)

**Storage information:**

- Storage temperature: 32 °F / 0 °C to 86 °F / 30 °C
- Can be used for 3 years under proper storage conditions.
- Use-by date: See label next to hourglass symbol 

**Warning and Safety notes:**

For safety and environmental information, cf. safety data sheets. These are available from the Miele websites listed below.

Dispose of containers only when empty and sealed. Disposal of product residue: refer to safety data sheet.

All information contained herein is current and correct to our knowledge. This does not however release users from obligations to perform their own trials and experiments. This cannot be deemed as a legal warranty of product characteristics or properties.

**Packaging unit:**

Designation	Packaging unit
<i>ProCare Dent 40</i>	1 l (33.8 fl oz (1.05 qt)) bottle

**Contact information:****Deutschland**

Miele Vertriebsgesellschaft Deutschland KG  
Carl-Miele-Straße 29  
33332 Gütersloh  
Tel.: +49 800 22 44 644  
vertrieb@miele-professional.de  
www.miele.de/professional

**Canada**

Miele Limited  
161 Four Valley Drive  
VAUGHAN  
Ontario L4K 4V8  
KANADA  
Tel. +1 905 532-2270  
customercare@miele.ca  
www.miele.ca

**USA**

Miele Inc.  
9 Independence Way  
Princeton, NJ 08540  
Tel.: 800-991-9380  
www.miele-professional.com

**Germany**

Chemische Fabrik  
Dr. Weigert GmbH & Co. KG  
Mühlhagen 85  
20539 Hamburg, Germany  
www.drweigert.com