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
PG 8504

Legend:

Circled, bold-type abbreviations:  
Connection required

Abbreviations surrounded by broken circle:  
Connection optional or required, depending on model

<b>(EL)</b>	Electrical connection	Voltage <b>(default)</b>	3 AC 208V 60 Hz	3 AC 208V 60 Hz
		Rated load	6.1 kW	6.1 kW
<b>(EL)</b>	Electrical connection	Fuse rating	3 x 20 A	3 x 20 A
		Minimum cross-section	4 x AWG 12	4 x 3.3 mm <sup>2</sup>
		Length of electrical cord (UL 921) with plug	approx. 5'9"	approx. 1.8 m
		Electrical connection	NEMA L15-20	NEMA L15-20
		Voltage <b>(convertible)</b>	AC 208V 60 Hz	AC 208V 60 Hz
<b>(EL)</b>	Electrical connection	Rated load	6.1 kW	6.1 kW
		Fuse rating	2 x 30 A	2 x 30 A
		Minimum cross-section	3 x AWG 10	3 x 5.2 mm <sup>2</sup>
		Length of electrical cord (UL 921) with plug	approx. 5'9"	approx. 1.8 m
		Electrical connection	NEMA L6-30	NEMA L6-30
<p>The machine is supplied with electrical cord and plug. The receptacle must be accessible after machine installation. For hard-wired machines, connection should be made via a suitable mains switch with all-pole isolation. The contact opening between all open contacts must be at least 3 mm wide and the mains switch must be lockable in the open position. If necessary, an equipotential bond with good contact connection must be provided in accordance with all appropriate national and local regulations.</p> <p>Splash protection class IP 21</p>				
<b>(PA)</b>	Equipotential bonding and earthing	<p>There is a screw connection point at the back of the machine, to which equipotential bonding can be connected if required.</p> <p><b>Connect equipotential bonding and protective earth conductor!</b></p>		
<b>(KW)</b> <b>(WW)</b>	Cold Water / Hot Water Connection	Max. temperature	150 °F	65 °C
		Max. water hardness	73 gpg	70 °dH
		Recommended flow pressure	29 psi	200 kPa
		Min. flow pressure	5.8 psi	40 kPa
		Max. pressure	145 psi	1,000 kPa
		Throughput	2 gal/min	7.5 l/min
		Connection thread - on site (flat sealing)	3/4" Male Garden Hose Thread	3/4" Male Garden Hose Thread
		Length of connection hose (supplied)	approx. 5' 7"	approx. 1.7 m

	Demineralised water (DI)	Max. temperature Recommended flow pressure Min. flow pressure Max. pressure  Connection thread - on site (flat sealing)  Throughput Length of connection hose (supplied)    <b>Model with feed pump for unpressurised demineralised water (ADP)</b> Minimum flow pressure (DI unpressurized) Maximum pressure (DI unpressurized) Machine connection nipple (outer dia. x l)  Hose from supply container to machine to be fitted on site.	150 °F 20.3 psi 4.4 psi 145 psi  3/4" Male Garden Hose Thread 2 gal/min approx. 5' 7"    1.3 psi 8.7 psi 1/4" x 3/16"	65 °C 140 kPa 30 kPa 1,000 kPa  3/4" Male Garden Hose Thread 7.5 l/min approx. 1.7 m    8.5 kPa 60 kPa 6 mm x 30 mm
	Drain Water	Drainage temperature Length of drain hose Drain hose, max. drainage length Drain pump head height from floor level, max. Drain hose delivery, max. On-site hose connector (external diameter x length)	200 °F approx. 4' 6" approx. 13' 1" approx. 3' 3" 4.2 gal/min 7/8" x 1 3/16"	93 °C approx. 1.4 m approx. 4.0 m approx. 1 m 16 l/min 22 mm x 30 mm
	External dispensing (optional)	Connection for up to two dispensing units for liquid agents on the back panel.		

	Machine feet	Height-adjustable Foot diameter Entry point into housing	0-2 3/8" 1 3/8" 8 M	0 – 60 mm 35 mm 8 M
	Heat Dissipation to Room	Heat transfer by emission Load on removal	1,194 btu/h 1,365 btu/h	0.35 kWh/h 0.40 kWh/h
	Machine data	Height Height without lid Width Depth Door height Net weight Dynamic floor load Width including transport pallet Depth including transport pallet Height including transport pallet	32 7/8" 32 1/4" 23 5/8" 23 5/8" 24 1/2" 163 lbs 1,200 N 26 3/8" 29 3/8" 36 1/4"	835 mm 820 mm 598 mm 598 mm 622 mm 74 kg 1,200 N 670 mm 740 mm 920 mm
	Set-up Conditions	Ambient temperature Max. humidity Max. altitude above sea level	40-105 °F 80% 4,921'	5-40 °C 80% 1,500 m
<p>It is essential to observe the installation instructions when installing the machine.</p> <p>In order to reduce the risk of steam damage, the area around the machine should be limited to furniture and fittings that are designed for use in commercial environments.</p> <p>Installation should only be performed by Miele Service or authorized installers in accordance with valid regulations, relevant standards and health and safety codes!</p> <p>Subject to change!</p>				

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