Installationsplan / Installation plan

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
PG 8536 CD

Legend:

Highlighted circled abbreviations:
Connection required

CDN-GB

Broken line circled abbreviations:
Connection optional or required depending on model version

EL	Electrical connection	1. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug		V/Hz kW A mm ² m	3N AC 400/50 9.7 3 × 16 5 × 2.5 1.7
	Country variations:				
		2. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug		V/Hz kW A mm ² m	3N AC 400/50 7.3 3 × 16 5 × 2.5 1.7
		Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug	Convertible	V/Hz kW A mm ² m	3 AC 230/50 7.3 3 × 20 5 × 2.5 1.7
		3. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug		V/Hz kW A mm ² m	3 AC 230/50 7.3 3 × 20 5 × 2.5 1.7
		Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug	Convertible	V/Hz kW A mm ² m	3 N AC 400/50 7.3 3 × 16 5 × 2.5 1.7


	4. Standard voltage (at delivery) Connected load Fuse rating (integrated in plug) Min. cross-sectional area of supply lead Connection cable length with plug		V/Hz kW A mm ² m	3N AC 400/50 9.2 3 x 15/16 5 x 2.5 1.7
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug	Convertible	V/Hz kW A mm ² m	AC 230/50 6.7 1 x 30 3 x 4 2.1
	5. Standard voltage (at delivery) Connected load Fuse rating (integrated in plug) Min. cross-sectional area of supply lead Connection cable length with plug		V/Hz kW A mm ² m	AC 230/50 6.7 1 x 30 3 x 4 2.1
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug	Convertible	V/Hz kW A mm ² m	3N AC 400/50 9.2 3 x 15/16 5 x 2.5 1.7
	6. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length with plug		V/Hz kW A mm ² m	3N AC 400 - 415/50 9.2 - 9.9 3 x 15 5 x 2.5 2.1
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length with plug	Convertible	V/Hz kW A mm ² m	AC 230 - 240/50 6.7 - 7.2 1 x 30 3 x 4 2.1
	7. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length with plug		V/Hz kW A mm ² m	AC 230 - 240/50 6.7 - 7.2 1 x 30 3 x 4 2.1
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length with plug	Convertible	V/Hz kW A mm ² m	3N AC 400 - 415/50 9.2 - 9.9 3 x 15 5 x 2.5 2.1
	8. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length without plug		V/Hz kW A AWG m	3 AC 208/60 6.0 3 x 20 4 x 12 1.9
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length without plug	Convertible	V/Hz kW A AWG m	AC 208/60 6.0 2 x 30 3 x 10 1.9

	9. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length without plug		V/Hz kW A AWG m	AC 208/60 6.0 2 × 30 3 × 10 1.9
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length without plug	Convertible	V/Hz kW A AWG m	3 AC 208/60 6.0 3 × 20 4 × 12 1.9
	10. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug		V/Hz kW A mm ² m	3 AC 200/50 5.7 3 × 20 5 × 2.5 2.3
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug	Convertible	V/Hz kW A mm ² m	AC 200/50 5.7 2 × 30 3 × 4 2.1
	11. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug		V/Hz kW A mm ² m	3 AC 200/60 5.7 3 × 20 5 × 2.5 2.3
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug	Convertible	V/Hz kW A mm ² m	AC 200/60 5.7 2 × 30 3 × 4 2.1
	12. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug		V/Hz kW A mm ² m	3N AC 380/60 6.7 3 × 16 5 × 2.5 2.3
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug	Convertible	V/Hz kW A mm ² m	3 AC 220/60 6.7 3 × 20 5 × 2.5 2.3
	13. Standard voltage (at delivery) Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug		V/Hz kW A mm ² m	3 AC 220/60 6.7 3 × 20 5 × 2.5 2.3
	Voltage Connected load Fuse rating Min. cross-sectional area of supply lead Connection cable length (H05(07)RN-F) without plug	Convertible	V/Hz kW A mm ² m	3N AC 380/60 6.7 3 × 16 5 × 2.5 2.3

		<p>It is recommended to make electrical connection via a plug and socket so that electrical safety checks can be carried out easily. The socket must be accessible after machine installation.</p> <p>With a fixed installation, this must be via a main switch to be provided on site, which must totally isolate the machine from the mains with a contact gap of at least 3 mm.</p> <p>The fitting of an earth leakage circuit breaker is strongly recommended.</p> <p>If necessary, an equipotential bond with good contact connection must be provided in accordance with all appropriate national and local regulations.</p>		
(KW)	Cold water	<p>Two cold water connections are required with machines with steam condenser (DK). If necessary connect the second cold water hose via the Y-piece provided</p> <p>Max. temperature Max. water hardness Min. flow pressure Max. pressure Flow rate On-site threaded union</p> <p>Connection hose length - Disinfector (supplied) Connection hose length - Steam condenser (supplied with DK variant)</p> <p>Machine is DVGW certified. Direct connection to mains water supply without special safety stopcock permitted.</p> <p>Water connection Australia: To prevent back-syphonage, the water connection must be via an approved Dual-Check-Valve supplied with the machine.</p>	<p>°C °dH kPa kPa l/min inch mm mm</p>	<p>20 60 250 1000 8 ¾ male thread (USA: 11.5 NH) 1500 1500</p>
(WW)	Hot water	<p>Max. temperature Max. water hardness Min. flow pressure Max. pressure Flow rate On-site threaded union</p> <p>Connection hose length (supplied)</p> <p>Machine is DVGW certified. Direct connection to mains water supply without special safety stopcock permitted.</p> <p>Water connection Australia: To prevent back-syphonage, the water connection must be via an approved Dual-Check-Valve supplied with the machine.</p>	<p>°C °dH kPa kPa l/min inch mm</p>	<p>70 60 250 1000 8 ¾ male thread (USA: 11.5 NH) 1500</p>
(VE)	Demineralised water	<p>Standard version: Min. flow pressure Min. flow pressure with extended water intake Max. pressure Flow rate</p> <p>Variant with demineralised water pump (aqua destillata pump) (ADP): Min. flow pressure Max. pressure</p> <p>On-site threaded union</p> <p>Connection hose length (supplied)</p> <p>Water connection Australia: To prevent back-syphonage, the water connection must be via an approved Dual-Check-Valve supplied with the machine.</p>	<p>kPa kPa kPa l/min kPa kPa inch mm</p>	<p>250 100 1000 8 5 30 ¾ male thread (USA: 11.5 NH) 1500</p>

AW	Drain water	Two drain hoses with DK machine variants		
		Drain hose (inner dia. x thickness x l) Drain pump (supplied) Steam condenser (supplied with DK variant) Max. drain pump head height from bottom edge of machine Max. brief flow rate per hose Hose connector on site for drain hose (outer dia. x l) Drain pump Steam condenser	mm mm m l/min mm mm	22 x 6 x 1500 22 x 6 x 1500 1.0 50 22 x 30 22 x 30
AL	Cabinet vent (optional)	Required on variants without DK (steam condenser) Vent connection, 2 (oval), for hose diameter (Connection hoses not supplied as standard) Use genuine Miele connection hoses.	mm	25
		Temperature - Average / Brief max. Relative air humidity - Average / Brief max.	°C/°C %/%	70/95 80/100
NL		I. Connection to an external vent system with fan: Flow rate per machine Vent coupling (with open air inlet) to be provided on site.	m³/h	70
		II. Vent direct to atmosphere without fan: Flow rate - Drying unit Max. permitted duct pressure loss Do not combine vent ducts! Ducts must be laid separately!	m³/h Pa	55 10
		Condensate must not flow back to the machine! The vent duct must first rise vertically and then slope downwards away from the machine. A condensate trap should be provided at the lowest point.		
		An adequate supply of fresh air must be available at the installation site to replace that extracted by the air vent.		
	Machine data	Height without lid Height with lid (optional) Width Depth Depth with removable side panels (standard delivery) Net weight Floor load in operation Min. installation location access width including transport pallet Min. installation location access height including transport pallet	mm mm mm mm mm kg N mm mm	820 850 900 600 700 120 2000 800 1100

When using plinths (UC 30-90/60-78, UC 30-90/70-78) these must be secured to the floor. It is essential to fasten the machine to the plinth.
The electrical, water, drain and exhaust connections must be made in accordance with all appropriate local and national safety regulations and applicable standards.
During installation it is essential to follow the installation instructions provided! Modification rights reserved! Dimensions in mm.

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