A 312

cs	Návod k obsluze Injektorový modul
el	Οδηγίες χρήσης μονάδα με εγχυτήρες
en	Operating instructions Injector module
hu	Használati utasítás Injektormodul
ja	取扱説明書 インジェクターモジュール
ko	사용 설명서 인젝터 모듈
pl	Instrukcja użytkowania Moduł iniekcyjny
ro	Instrucțiuni de utilizare Modul de injectie
ru	Инструкция по эксплуатации Инжекторный модуль
sk	Návod na použitie Injektorový modul
tr	Kullanım Kılavuzu Enjektör Modülü
zh	操作说明书 注射器模块

cs 4
el 14
en 24
hu 34
ja 44
ko 53
pl 63
ro 73
ru 83
sk 93
tr103
zh113

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Notes about these instructions	24
Queries and technical problems.....	24
Definition of terms	24
Appropriate use	25
Items supplied	26
Scope of delivery A 312	26
Disposal of the packing material	26
Warning and Safety instructions	27
Installation	28
Tools required:.....	28
Injector nozzles.....	28
Shortening the plastic hose.....	29
Areas of application	30
Injector modules.....	30
Inserting an injector module	30
Removing an injector module.....	30
Preparing the load	31
Before loading the machine and before starting a programme	32

Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

Queries and technical problems

In the event of queries or technical problems, please contact Miele. Contact details can be found at the end of the operating instructions for your cleaning machine or at www.miele.com/professional.

Definition of terms

- | | |
|-------------------------|---|
| Cleaning machine | In these operating instructions, the washer-disinfector is referred to as “the cleaning machine”. |
| Load items | The term “load items” is used wherever the items to be processed are not defined in any further detail. |
| Wash water | The term “wash water” is used for the mixture of water and process chemicals. |

This module can be used to reprocess machine-reprocessable laboratory glassware and utensils in a Miele washer-disinfector for laboratory glassware and utensils. Follow the operating instructions for your washer-disinfector as well as the instructions of the glassware and utensil manufacturer on how to reprocess their items by machine.

The A 312 injector module is designed for reprocessing glass capillary viscometers with the following dimensions:


- | | |
|--|---------------------------|
| - Height | Max. 310 mm |
| - Capillary tube inner diameter | Min. \varnothing 3.0 mm |
| - Filling tube inner diameter | Min. \varnothing 5.0 mm |
| - Centre of capillary tube to centre of filling tube | Approx. 20-25 mm |

The module has 8 nozzles for reprocessing 8 capillary viscometers.

Up to 2 injector modules can be used in one A 150 lower basket. A maximum of 16 capillary viscometers can be reprocessed per cycle.

No upper baskets can be used when the A 312 injector module is used.

Scope of delivery A 312

- A 312 injector module
 - Height 73 mm, width 133 mm, depth 475 mm
 - Max. vertical clearance (from injector nozzle supports) 310 mm

Items supplied

- A 867, injector nozzle with 2 irrigation connectors for viscometer; included: 1 nozzle with plastic support, length 110 mm, Ø 2.5 mm and 1 plastic hose, length ex works approx. 300 mm, Ø 3 mm

Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in landfill sites.

en - Warning and Safety instructions

To avoid the risk of accidents and damage to this module please read these instructions carefully before using it for the first time. Keep these instructions in a safe place where they are accessible to users at all times.

 Please also read the operating instructions for your washer-disinfector and pay particular attention to the Warning and Safety instructions.


- ▶ The module is approved solely for the applications specified in the “Areas of application” chapter of these operating instructions. Components such as nozzles may only be replaced with Miele components or genuine original spare parts.
- ▶ Inspect all mobile units, baskets, modules and inserts daily as described in the “Maintenance” section in the Operating instructions for the washer-disinfector.
- ▶ Only items which have been declared by their manufacturer as suitable for machine reprocessing may be processed. The manufacturer's specific reprocessing instructions must be observed.
- ▶ New load carriers must be cleaned in the washer-disinfector without a load prior to first use.
- ▶ Broken glass can result in serious injury when loading or unloading. Broken glass items must not be processed in the machine.
- ▶ Any hollow items that have become disconnected from their injector nozzles, adapters or holders during processing must be re-processed.
- ▶ Only place empty, unloaded modules in baskets. Check that they are correctly engaged before loading them. Modules must be completely emptied before taking them out. Placing loaded modules in baskets or removing them can cause damage to the wash load and cause injury should glassware be broken.
- ▶ The reprocessing result must be checked as necessary using a suitable test method and not just by a visual check.

Tools required:

- 9 mm Spanner (WAF 9)
- 13 mm spanner (WAF 13)

Injector nozzles

The nozzles of the A 312 module consist of a Y-piece with 2 irrigation connections. The supplied nozzles must be attached to the branches of the Y-piece.


- Screw the injector nozzle with the plastic support onto the longer branch of the Y-piece.
- Screw the nozzle with the plastic hose into the shorter branch of the Y-piece.
- Use a spanner to tighten both nozzles.
- Screw the Y-piece into one of the module's holders and use a spanner to tighten this too.
- Repeat the process until nozzles have been attached to all of the module's threaded holders.

Shortening the plastic hose

The plastic hose of the A 867 nozzle must be cut to the height of the viscometers.

- To do this, place the viscometer with the capillary tube on the nozzle with the plastic support. The viscometer must rest on the plastic support.
- Turn the viscometer until the filling tube is alongside the plastic hose.

The plastic hose must end above the container (reservoir, bulb) in the connection piece between the capillary tube and filling tube. For re-processing, it is important that the hose is cut to a point.


- With sharp sturdy scissors, carefully cut the plastic hose at an angle of approx. 45° in the direction of the nozzle to form a point at the end of the hose. The lower end of the hose must be pointing towards the nozzle.

Injector modules

⚠ Only place empty, unloaded modules in baskets. Check that they are correctly engaged before loading them.
Modules must be completely emptied before taking them out.
Placing loaded modules in baskets or removing them can cause damage to the wash load and cause injury should glassware be broken.

Inserting an injector module

This injector module is only suitable for use in lower basket A 150. The injector nozzles must point upwards.

- Carefully slide the injector module in and connect it to the water connection point.
- Then press down on the module to secure it on the rails in the basket.

Removing an injector module

⚠ Do not hold the basket by the injector module when removing it from the cleaning machine.
The module could come loose or the water connection could be damaged.


- Release the module from the locking mechanism and pull it out of the water connection at a flat angle. If the module is removed at too steep an angle, the lever action may damage the connection point and the connector.

Preparing the load

The A 312 injector module is designed for reprocessing glass capillary viscometers with the following dimensions:

- | | |
|--|---------------------------|
| - Height | Max. 310 mm |
| - Capillary tube inner diameter | Min. \varnothing 3.0 mm |
| - Filling tube inner diameter | Min. \varnothing 5.0 mm |
| - Centre of capillary tube to centre of filling tube | Approx. 20-25 mm |

- Check that the viscometers and the capillaries in particular are free of blockages.

- Place the viscometers on the A 867 injector nozzles.
 - The capillary tubes must be connected to the injector nozzles with the plastic support.
 - The filling tubes must be connected to the nozzles with the plastic hoses.
 - If there are any ventilation tubes, these do not require extra cleaning.

Make sure that the ends of the plastic hoses do not touch the glass walls of the viscometers.

If they do, the vessels will not be cleaned satisfactorily.

- The height of the plastic supports can be adjusted. Move the plastic supports, if necessary. For example, if the viscometers have tubes with different lengths, the supports may need to be adjusted to compensate for the length difference. Or if the tubes are particularly close together, the supports may need to be moved to prevent the tubes from being bent apart.

Before loading the machine and before starting a programme

- Check that fittings such as irrigation sleeves and injector nozzles are securely screwed into position.

⚠ Make sure that all screw connectors are fitted with injector nozzles, adapters, irrigation sleeves or blind stoppers to ensure that all fittings in use are supplied with sufficient standardised pressure.

Damaged fittings such as injector nozzles, adapters and irrigation sleeves must not be used.

Fittings not equipped with wash items do not have to be replaced blind stoppers.

- Check that the baskets and modules are correctly connected to the water supply.

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