A 313

da	Brugsanvisning Injektormodul
de	Gebrauchsanweisung Injektormodul
en	Operating instructions Injector module
es	Instrucciones de manejo Módulo inyector
fi	Käyttöohje Suorasuihkumoduuli
fr	Mode d'emploi Module à injection
hr	Upute za uporabu Injektorski modul
it	Istruzioni d'uso Modulo iniettore
nl	Gebruiksaanwijzing injectormodule
no	Bruksanvisning injektormodul
pt	Instruções de utilização Módulo injetor
sl	Navodila za uporabo Injektorski modul
sv	Bruksanvisning injektormodul

da	4
de	18
en	32
es	46
fi	60
fr	74
hr	88
it	102
nl	116
no	130
pt	144
sl	158
sv	172

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Notes about these instructions	32
Queries and technical problems	32
Appropriate use	33
Items supplied	34
Scope of delivery A 313	34
Disposal of the packing material	34
Optional components	35
Warnings and safety notes	36
Areas of application	37
Installation	37
Tools required	37
Loading dimensions	37
Adjusting the height of the frame	37
Inserting retaining bars	38
Flasks and round flasks	39
Installing A 846 for flasks and round flasks	39
Preparing the load	40
Erlenmeyer flasks	41
Installing A 847 for Erlenmeyer flasks	41
Preparing the load	41
Measuring flasks	42
Installing A 848 for measuring flasks	42
Preparing the load	42
Injector modules	43
Inserting the module	43
Removing an injector module	43
To check before loading the machine and before starting a programme	44
Knocking noises during reprocessing	44

Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

Queries and technical problems

In the event of queries or technical problems, please contact Miele. Contact details can be found at the end of the operating instructions for your cleaning machine or at www.miele.com/professional.

This module can be used to reprocess machine-reprocessible laboratory glassware and utensils in a Miele washer-disinfector for laboratory glassware and utensils. Follow the operating instructions for your washer-disinfector as well as the instructions of the glassware and utensil manufacturer on how to reprocess their items by machine.

The A 313 module is designed for reprocessing large-volume laboratory glassware with a maximum diameter of 232 mm. The maximum vertical clearance is determined by the type of components used.

The module is supplied without injector nozzles. Miele offers various component sets which can be fitted based on the relevant requirements. The sets consist of injector nozzles and holders which are designed for reprocessing different types of large-volume laboratory glassware.

Depending on the height and diameter of the load items, the module is suitable for laboratory glassware with the following volumes:

- Laboratory flasks from 2 l to 10 l (max. height 410 mm)
- Round flasks from 2 l to 6 l
- Erlenmeyer flasks from 2 l to 5 l
- Measuring flasks from 2 l to 5 l (max. height 470 mm)

This module is suitable for use in the following:

- Lower basket A 150
- APLW 106 basket

Because of its height, this module can only be used in a lower basket. It is not possible to use an upper basket at the same time.

Up to 2 injector modules can be placed in one lower basket.

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessable laboratory glassware and utensils are referred to as “items” if they are not more closely defined.

Scope of delivery A 313

- A 313 injector module
 - Height 221 mm, width 254 mm, depth 495 mm
 - The max. vertical clearance is determined by the type of accessories used.

Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in land-fill sites.

Optional components

- A 846, injector nozzle, length 127 mm, Ø 10 mm, support and 2 rods
This combination is particularly suitable for laboratory flasks and round flasks.
 - A 847, injector nozzle, length 127 mm, Ø 10 mm and support
Ø 142 mm
This combination is particularly suitable for Erlenmeyer flasks.
 - A 848, injector nozzle, length 350 mm, Ø 10 mm and support
Ø 142 mm
This combination is particularly suitable for measuring flasks.
- Further components are available from Miele as optional extras.

en - Warnings and safety notes

To avoid the risk of accidents and damage to this module please read these instructions carefully before using it for the first time. Keep these instructions in a safe place where they are accessible to users at all times.

 Please also read the operating instructions for your dishwasher or the washer-disinfector and pay particular attention to the Warning and Safety instructions.

- ▶ The module is approved solely for the applications specified in the “Areas of application” chapter of these operating instructions. Components such as nozzles may only be replaced with Miele components or genuine original spare parts.
- ▶ New load carriers must be cleaned in the washer-disinfector without a load prior to first use.
- ▶ Inspect all mobile units, baskets, modules and inserts daily as described in the “Maintenance” section in the Operating instructions for the washer-disinfector.
- ▶ Only items which have been declared by their manufacturer as suitable for machine reprocessing may be processed. The manufacturer's specific reprocessing instructions must be observed.
- ▶ Broken glass can result in serious injury when loading or unloading. Broken glass items must not be processed in the machine.
- ▶ Any hollow items that have become disconnected from their injector nozzles, adapters or holders during processing must be re-processed.
- ▶ Only place empty, unloaded modules in baskets. Check that they are correctly engaged before loading them. Modules must be completely emptied before taking them out. Placing loaded modules in baskets or removing them can cause damage to the wash load and cause injury should glassware be broken.
- ▶ The reprocessing result must be checked as necessary using a suitable test method and not just by a visual check.


Installation

Tools required

- 17 mm spanner (WAF 17)
- TORX T20 screwdriver

Loading dimensions


The maximum diameter of the items depends on the position on the module and the number of nozzles used.

- ① Ø max. 232 mm
- ② Ø max. 236 mm
- ③ Ø max. 146 mm

Adjusting the height of the frame

The height of the frame can be adjusted to different levels. The height of the frame is determined by the nozzles used and the items to be re-processed.

- Use a TORX T 20 screwdriver to remove the 4 screws on the struts of the module ①.
- Move the frame to the required height ② and secure back in place using the 4 screws ①.


en - Areas of application

Inserting retaining bars

The component sets for reprocessing large-volume laboratory glassware include rods or supports which are used to support the items during reprocessing. These must be secured in place with the aid of the module's retaining bars.

The supports are matched to the shape of the items used. They must be secured in place using the retaining bars.

The rods can be used as required and are combined with the retaining bars to support the load items. Once the retaining bars have been inserted, the distances between the rods and the locking pins must be roughly equal to the diameter of the load items to ensure that the items move as little as possible during reprocessing.


- Place the rods and/or supports on the top of the frame.
- Guide one end of the retaining bar through the relevant opening on the frame up to the corresponding opening on the opposite side ① and push the bar through the opening.
- Turn the retaining bar until the locking pins at the ends of the bar are pointing downwards ②.

The rods and/or supports are now secured in place.


Installing A 846 for flasks and round flasks

Flasks and round flasks

- Adjust the height of the frame.
 - Flasks:
The height of the frame should be adjusted so that it comes to approximately halfway up the flasks.
 - Round flasks:
The height of the frame should be adjusted so that it is approximately at the point of the largest diameter on the round flasks.
- Remove one of the knurled nuts out of the module's injector rails.

- Place the support for the A 846 on the threaded fitting and screw the injector nozzle in ①.
- Use a 17 mm spanner (WAF 17) to tighten the injector nozzle ②.

- Place the rods on the frame above the injector nozzle.
- Secure the rods on the right and left using the module's retaining bars.

The distances between the rods and the retaining bars must be roughly equal to the diameter of the load items.

en - Areas of application

Preparing the load

Flasks

The laboratory flasks rest on 4 points of the support; the neck of the flask should not touch the water supply pipe.

The rods inserted in the frame prevent the flasks from touching each other if they are moved by the force of water during reprocessing.

2 l laboratory flasks are held securely on the support and do not need any additional supports at the sides. Rods are therefore not required.

Round flasks


The round flasks rest on 4 points of the support; their necks should not touch the water supply pipe.

The rods inserted in the frame prevent the round flasks from touching each other if they are moved by the force of water during reprocessing.

Erlenmeyer flasks

Installing A 847 for Erlenmeyer flasks

- Adjust the height of the frame.
The height of the frame should be adjusted so that the items do not touch the module's injector rail or the injector nozzle's screw thread.

- Remove a knurled nut from the module's injector rail and replace it with the injector nozzle (1).
- Use a 17 mm spanner (WAF 17) to tighten the injector nozzle (2).
- Place the support on the frame above the injector nozzle (1).
- Secure the support on the right and left using the module's retaining bars (3).

Preparing the load


The Erlenmeyer flasks rest on the support in the module frame. The cone on the bottom of the injector nozzle should be inside the neck of the item. This protects the rim from becoming chipped.

Measuring flasks


Installing A 848 for measuring flasks

- Adjust the height of the frame.
The height of the frame should be adjusted so that the items do not touch the module's injector rail or the injector nozzle's screw thread.

- Remove a knurled nut from the module's injector rail and replace it with the injector nozzle ①.
- Use a 17 mm spanner (WAF 17) to tighten the injector nozzle ②.
- Place the support on the frame above the injector nozzle ①.
- Secure the support on the right and left using the module's retaining bars ③.

Preparing the load


The measuring flasks rest on the support in the module frame.
The cone on the bottom of the injector nozzle should be inside the neck of the item. This protects the rim from becoming chipped.

Injector modules

⚠ Only place empty, unloaded modules in baskets. Check that they are correctly engaged before loading them. Modules must be completely emptied before taking them out. Placing loaded modules in baskets or removing them can cause damage to the wash load and cause injury should glassware be broken.

Inserting the module

Ensure that the module is correctly orientated when inserting it. For example, the nozzles must point upwards.

Example of a flat insertion angle

- Insert the module into the water connection with the connector at a flat angle.
- Then engage the module by pushing it down onto the struts of the basket.

Removing an injector module

⚠ Do not hold the basket by the injector module when removing it from the cleaning machine. The module could come loose or the water connection could be damaged.

- Release the module from the locking mechanism and pull it out of the water connection at a flat angle. If the module is removed at too steep an angle, the lever action may damage the connection point and the connector.

The supports and rods are not designed to bear the weight of the module.

Therefore, only lift the module by the handle or frame.

To check before loading the machine and before starting a programme

- Check that fittings such as irrigation sleeves and injector nozzles are securely screwed into position.

⚠ Make sure that all screw connectors are fitted with injector nozzles, adapters, irrigation sleeves or blind stoppers to ensure that all fittings in use are supplied with sufficient standardised pressure. Damaged fittings such as injector nozzles, adapters and irrigation sleeves must not be used.

Fittings not equipped with wash items do not have to be replaced blind stoppers.

- Check that the baskets and modules are correctly connected to the water supply.

Knocking noises during reprocessing

If load items are not attached to one or more module nozzles, you may hear knocking noises during reprocessing. The noises occur when water jets from the empty injector nozzles hit the top spray arm. This will not affect the cleaning performance in any way.

- Attach load items to all injector nozzles to prevent the knocking noises from occurring.

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