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
PG 8581 - 8583

Legend:

Circled, bold-type abbreviations:
Connection required

Abbreviations surrounded by broken circle:
Connection optional or required, depending on model

EL	Electrical connection (220 V)	Voltage (standard version) Connected load Fuse rating Connection cable, min. cross-sectional area Length of mains connection cable (H05(07)RN-F)	V/Hz kW A mm ² m	3 AC 220 – 230 V/60 6.9 – 7.5 3 x 20 4 x 2.5 1.6
EL	Electrical connection	Voltage (convertible) Connected load Fuse rating Connection cable, min. cross-sectional area Length of mains connection cable (H05(07)RN-F)	V/Hz kW A mm ² m	AC 220 – 230 V/60 3.0 -3.3 2 x 16 3 x 1.5 1.8
EL	Electrical connection (380 V)	Voltage (Standard version) Connected load Fuse rating Connection cable, min. cross-sectional area Length of mains connection cable (H05(07)RN-F)	V/Hz kW A mm ² m	3N AC 380 – 400 V/60 8.5 – 9.3 3 x 16 5 x 2.5 1.6
EL	Electrical connection	Voltage (convertible) Connected load Fuse rating Connection cable, min. cross-sectional area Length of mains connection cable (H05(07)RN-F)	V/Hz kW A mm ² m	AC 220 – 230 V/60 3.0 – 3.3 1 x 16 3 x 1.5 1.8
		We recommended connecting the washer-disinfector to the electrical supply via a plug and socket so that electrical safety checks can be carried out easily. The plug and socket must be easily accessible after installation. A hard-wired connection is possible. It must be made in accordance with IEE (IET) Wiring regulations BS 7671: 2008. Splash protection class IP 21		
PA	Equipotential bonding and protective earth conductor	There is a screw connection point at the back of the washer-disinfector, to which equipotential bonding can be connected if required. Connect equipotential bonding and protective earth conductor!		
NW	Network/ printer connection (optional)	Optional machine interface for transmission or process data output Communication module (not supplied) Ethernet module or RS232 module (install electrical sockets near machine if required) Connections/installation must comply with IEC 60950.		

(KW)	Cold water	Two cold water connections are required for models with a steam condenser (DK). The second cold water hose can be connected using the Y-piece provided.		
		Maximum temperature Max. water hardness Minimum flow pressure Minimum flow pressure with extended water intake Maximum pressure Flow rate On-site threaded union in accordance with DIN 44991 (flat seal) Length of cold water inlet hose Length of steam condenser hose Connection to a drinking water supply in the UK requires the installation of an additional WRAS approved non-return device.	°C mmol/l °dH kPa kPa kPa l/min Inch m m	20 12.6 70 200 100 1,000 7.5 3/4" male thread Approx. 1.7 Approx. 1.7
(WW)	Hot water (not available for PG 8581)	Maximum temperature	°C	65
		Max. water hardness Minimum flow pressure Minimum flow pressure with extended water intake Maximum pressure Flow rate On-site threaded union in accordance with DIN 44991 (flat seal) Length of hot water inlet hose Connection to a drinking water supply in the UK requires the installation of an additional WRAS approved non-return device.	mmol/l °dH kPa kPa kPa l/min Inch m	12.6 70 200 40 1,000 7.5 3/4" male thread Approx. 1.7
(VE)	Demineralised water  (PG 8581 note machine version) PG 8583 only:	Maximum temperature Minimum flow pressure (AD pressure resistant) Minimum flow pressure with extended water intake Maximum pressure (AD pressure resistant) Flow rate On-site threaded union in accordance with DIN 44991 (flat seal) Length of demin. water inlet hose Connection to a drinking water supply in the UK requires the installation of an additional WRAS approved non-return device.	°C kPa kPa kPa l/min Inch m	65 200 30 1,000 7.5 3/4" male thread Approx. 1.7
		Model with feed pump for unpressurised demineralised water (ADP) Minimum flow pressure (AD unpressurised) Maximum pressure (AD unpressurised) Machine connection nipple (outer dia. x l) Hose from supply container to machine to be fitted on site	kPa kPa mm	8.5 60 6 x 30
(AW)	Waste water	Drainage temperature Length of drain hose Max. drain hose extension Max. drain pump head height from base of machine Max. transient flow rate On-site sleeve for drain hose (dia. x length) It is advisable to lay the drain hose in a curve min. 600 mm above floor level to the drainage point.	°C m m m l/min mm	93 Approx. 1.4 Approx. 4.0 1.0 16 22 x 30

DS	External dispensing (optional) (note machine model)	Connection for up to two dispensing units for liquid agents on the back panel.		
	Machine feet	Height-adjustable Machine foot diameter Machine foot socket	mm mm M	-0/+60 35 8
	Heat dissipation rate to installation site	Heat dissipation through heat radiation Load heat dissipation on removal	kWh/h kWh/h	0.35 0.40
	Machine data	Undercounter height Height incl. lid Width Depth Door height Net weight Floor load in operation Min. access width, incl. transport pallet Min. access depth incl. transport pallet Min. access height, incl. transport pallet	mm mm mm mm mm kg N mm mm mm	820 835 598 598 622 74 1,200 670 740 920
	Installation conditions	Permissible ambient temperature Maximum relative humidity Max. installation above sea level	°C % m	5-40 80 1,500
<p>It is essential to observe the installation instructions when installing the washer-disinfector.</p> <p>In order to reduce the risk of steam damage, the area around the washer-disinfector should be limited to furniture and fittings that are designed for use in commercial environments.</p> <p>Installation should only be performed by qualified fitters in accordance with valid regulations, relevant standards and health and safety codes!</p> <p>All rights reserved! Measurements in mm</p>				

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