

## ProCare

### Dent 30 P


## Neutralizing agent

#### Use:

Liquid concentrate, acidic, based on phosphoric acid. For the automated reprocessing of dental instruments in Miele Professional washer-disinfectors.

*ProCare Dent 30 P* is not suitable for use on transmission instruments, non-ferrous metals, anodized aluminium or chrome- or nickel-plated items. The ProCare Dent 30 C neutralizing agent is recommended for use on transmission instruments.

#### Quality features:

*ProCare Dent 30 P* complies with European guidelines for medical devices. 

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

#### Product properties:


- Removes tarnishing and soils soluble in acids
- Excellent rinsing properties
- Does not contain surfactants

#### Use and dispensing:

Miele *ProCare Dent 30 P* is formulated for use in processes in Miele Professional washer-disinfectors.

- Precise and economical dispensing via dispenser pump.
- Suitable dosing dispensers should be used.
- Ensure that the load and drain pipework are sufficiently compatible with acids.
- Application concentration: 1 ml/l
- Program selection: Vario TD

#### Program sequence:



#### vario TD process:

- 1 Pre-wash with cold water to remove coarse soil (e.g. blood).
- 2 Cleaning Main wash
- 3 Neutralization with *ProCare Dent 30 P*
- 4 Intermediate rinse with cold water
- 5 Thermal disinfection at 200 °F / 93 °C with 5-minute holding time in the final rinse
- 6 Drying

The product in solution should be completely rinsed off using water. In Miele Professional washer-disinfectors, this is guaranteed by the Vario TD program when properly used (e.g. correct loading). Verification is provided by measuring the increase of the conductivity in the final rinse during commissioning and routine maintenance.

Miele recommends the use of fully demineralized water in the final rinse in order to protect the load and the washer-disinfectors against stains and deposits.

#### Please note!

- Do not mix with other products.
- Contact Miele Service when switching products.
- Reprocessing should comply with all ordinances pursuant to medical device regulations and should be performed with appropriate validated processes.
- Observe recommendations on reprocessing issued by medical device manufacturers with respect to the provisions of DIN EN ISO 17664.
- Operating instructions on Miele Professional washer-disinfectors must be observed.
- For commercial applications only.

## ProCare

### Dent 30 P

### Neutralizing agent


**Technical data:**

- Density (at 68 °F / 20 °C): approx. 1.5 g/cm<sup>3</sup>
- pH-value (determined using fully demineralized water at 68 °F / 20 °C) 1 ml/l: approx. 2.3

**Ingredients according to Regulation (EC) 648/2004 on detergents:**

- > 30 % phosphates

**Storage information:**

- Storage temperature: 5 °F to 86 °F / -15 °C to 30 °C
- Can be used for 3 years under proper storage conditions
- Use-by date: See label next to hourglass symbol 

**Warning and Safety notes:**

For safety and environmental information, cf. safety data sheets. These are available from Miele websites listed below.

Dispose of containers only when empty and sealed. Disposal of product residue: refer to safety data sheet.

All information contained herein is current and correct to our knowledge. This does not however release users from obligations to perform their own trials and experiments. This cannot be deemed as a legal warranty of product characteristics or properties.

**Packaging units and accessories:**

Designation	Packaging unit	Mat. no.
<i>ProCare Dent 30 P</i>	5 l (1.32 gal) canister	12 187 380
<i>ProCare Dent 30 P</i>	1 l (33.8 us fl oz) bottle Type 1	12 187 370
Canister key	For 5 l and 20 l canisters	11 851 600
DTD2 Drip tray	For 5 l and 10 l canisters	10 695 460

**Contact information:****Germany**

Miele & Cie. KG  
Carl-Miele-Straße 29,  
33332 Gütersloh  
Tel.: 49 (0) 8002244644  
vertrieb@miele-professional.de  
www.miele.de/professional

**USA**

Miele Inc.  
9 Independence Way  
Princeton, NJ 08540  
Tel.: 800-991-9380  
www.miele-professional.com

**Germany**

Chemische Fabrik  
Dr. Weigert GmbH & Co. KG  
Mühlenhagen 85  
20539 Hamburg  
www.drweigert.de