A 612

de	Gebrauchsanweisung Injektormodul
da	Brugsanvisning Injektormodul
el	Οδηγίες χρήσης μονάδα με εγχυτήρες
en	Operating instructions Injector module
es	Instrucciones de manejo Módulo inyector
fi	Käyttöohje Suorasuihkumoduuli
fr	Mode d'emploi Module à injection
it	Istruzioni d'uso Modulo iniettore
nl	Gebruiksaanwijzing injectormodule
no	Bruksanvisning injektormodul
pt	Instruções de utilização Módulo injetor
sv	Bruksanvisning injektormodul

de	4
da	20
el	34
en	50
es	65
fi	81
fr	95
it	110
nl	125
no	141
pt	155
sv	171

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Notes about these instructions	50
Appropriate use	51
Queries and technical problems	51
Items supplied	52
Load carriers.....	52
Loading dimensions.....	52
Disposal of the packing material	53
Optional accessories.....	53
Warnings and safety notes	54
Areas of application.....	55
To check before loading the machine and before starting a programme	55
Preparing the load	56
A 612/A 846 with flasks	56
A 612/A 846 with round flasks	56
A 612/A 847 with Erlenmeyer flasks.....	58
A 612/A 848 with measuring flasks.....	58
Installation.....	59
Tools required.....	59
Adjusting the height of the frame.....	59
Screwing in nozzles.....	60
A 846 for flasks and round flasks	60
A 847 for Erlenmeyer flasks	62
A 848 for measuring flasks	63

Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

This module can be used to reprocess machine-reprocessible laboratory glassware and utensils in a Miele washer-disinfector for laboratory glassware and utensils. Follow the operating instructions for your washer-disinfector as well as the instructions of the glassware and utensil manufacturer on how to reprocess their items by machine.

A 612 injector modules are suitable for reprocessing large-volume laboratory glassware.


The module can be used in an A 503 mobile unit.

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessible laboratory glassware and utensils are referred to as “items” if they are not more closely defined.

Queries and technical problems

If you have any queries or technical problems that you would like to discuss please contact Miele. Contact details can be found at the back of the Operating Instructions for your dishwasher or the washer-disinfector or go to www.miele-professional.com.


Load carriers

- A 612 module, height 361 mm, width 255 mm, depth 614 mm

Loading dimensions

The maximum diameter of the items depends on the number of nozzles used.

Ø max. 136 mm Ø max. 182 mm Ø max. 240 mm

The module is suitable for items with the following volumes:

- Laboratory flasks from 2 l to 10 l
- Round flasks from 2 l to 6 l
- Erlenmeyer flasks from 2 l to 5 l
- Measuring flasks from 2 l to 5 l

Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in land-fill sites.

Optional accessories

The following and other accessories are available to order from Miele:

- A 846, injector nozzle, length 127 mm, Ø 10 mm, support and 2 rods
This combination is particularly suitable for laboratory flasks and round flasks.
- A 847, injector nozzle, length 127 mm, Ø 10 mm and support Ø 142 mm
This combination is particularly suitable for Erlenmeyer flasks.
- A 848, injector nozzle, length 350 mm, Ø 10 mm and support Ø 142 mm
This combination is particularly suitable for measuring flasks.

en - Warnings and safety notes

Read the operating instructions carefully before using this load carrier. This will help protect users from personal injury, and help protect the load carrier from damage.
Keep these operating instructions in a safe place.

 Please also read the operating instructions for your dishwasher or the washer-disinfector and pay particular attention to the Warning and Safety instructions.

- ▶ The module is approved solely for the applications specified in the operating instructions. Components such as nozzles may only be replaced with Miele accessories or genuine Miele original spare parts. Alterations or conversion of the machine, or using it for purposes other than those for which it was intended, are not permitted and could be dangerous.
- ▶ New load carriers must be cleaned in the washer-disinfector without a load prior to first use.
- ▶ Inspect all load carriers as described in the “Maintenance” section in the operating instructions for the cleaning machine.
- ▶ Only items which have been declared by their manufacturer as suitable for machine reprocessing may be processed. The manufacturer's specific reprocessing instructions must be observed.
- ▶ Broken glass and ceramics can result in serious injury when loading or unloading. Damaged glass or ceramic load items must not be reprocessed in the machine.
- ▶ Only ever place empty, unloaded modules in the mobile unit. Check that they are correctly engaged before loading them. Modules must be completely emptied before removing them. Placing loaded modules in mobile units or removing them can cause damage to the wash load and cause injury should glassware be broken.
- ▶ The reprocessing result must be checked as necessary using a suitable test method and not just by a visual check.

Miele cannot be held liable for damage or injury caused by non-compliance with the warning and safety notes.

The operating instructions for the loading carriage describe how to insert and remove the modules.

To check before loading the machine and before starting a programme

- Check that fittings such as irrigation sleeves and injector nozzles are securely screwed into position.

⚠ Make sure that all screw connectors are fitted with injector nozzles, adapters, irrigation sleeves or blind stoppers to ensure that all fittings in use are supplied with sufficient standardised pressure. Damaged fittings such as injector nozzles, adapters and irrigation sleeves must not be used.

Fittings not equipped with wash items do not have to be replaced blind stoppers.

- Is the module correctly docked to the water supply in the load carrier?

Preparing the load

A 612/A 846 with flasks

The laboratory flasks rest on 4 points of the support; the neck of the flask should not touch the water supply pipe.

The height of the frame should be adjusted so that it comes to approximately halfway up the flasks.

The rods inserted in the frame prevent the flasks from touching each other if they are moved by the force of water during reprocessing.

2 l laboratory flasks are held securely on the support and do not need any additional supports at the sides. Rods are therefore not required.

A 612/A 846 with round flasks

The round flasks rest on 4 points of the support; their necks should not touch the water supply pipe.

The height of the frame should be adjusted so that it is approximately at the point of the largest diameter on the round flasks.

The rods inserted in the frame prevent the round flasks from touching each other if they are moved by the force of water during reprocessing.

en - Areas of application

A 612/A 847 with Erlenmeyer flasks

The Erlenmeyer flasks rest on the support in the module frame. The height of the frame should be adjusted so that the items do not touch the water supply pipe. The cone on the bottom of the injector nozzle should be inside the neck of the item. This protects the rim from becoming chipped.

A 612/A 848 with measuring flasks

The measuring flasks rest on the support in the module frame. The height of the frame should be adjusted so that the items do not touch the water supply pipe. The cone on the bottom of the injector nozzle should be inside the neck of the item. This protects the rim from becoming chipped.

Tools required


- 17 mm spanner (WAF 17)
- TORX T20 screwdriver

Adjusting the height of the frame

The height of the frame can be adjusted to 4 levels. The height of the frame is determined by the nozzles used and the items to be reprocessed:

The bottom hole in the struts of the frame cannot be used if the module is to be placed in the A 503 mobile unit. The frame cannot be placed under the mobile unit's water supply pipe at this height.


- If using laboratory flasks on the A 846 support/nozzle, the frame should be approximately halfway up the flask.
- If using round flasks on the A 846 support/nozzle, the frame should be approximately at the point of the largest diameter.
- If using Erlenmeyer flasks on the A 847 support/nozzle and measuring flasks on the A 848 support/nozzle, the frame should be adjusted so that the items do not touch the water supply pipe.


- Use a TORX T 20 screwdriver to remove the 4 screws on the struts of the module ①.
- Lift the frame up to the required height ② and secure back in place using the 4 Torx screws ①.

Screwing in nozzles

A 846 for flasks
and round flasks


- Remove the knurled nut from the screw attachment.
- Place the support on the screw attachment and turn the injector nozzle to insert it into the attachment ①.
- Use a 17 mm spanner (WAF 17) to tighten the injector nozzle ②.


- Place the rods on the frame above the injector nozzle.

Secure the rods on the right and left using the module's retaining bars.
To do this:


- Guide one end of the retaining bar through the relevant opening on the frame up to the corresponding opening on the opposite side ① and push it through.
- Turn the retaining bar so that the locking pin secures the position ②.

A 847 for Erlenmeyer flasks

- Remove the knurled nut from the screw attachment.
- Turn the injector nozzle to insert it into the screw attachment ①.
- Place the support on the frame above the injector nozzle.
- Use a 17 mm spanner (WAF 17) to tighten the injector nozzle ②.
- Secure the supports on the right and left using the module's retaining bars ③.

A 848 for measuring flasks

- Remove the knurled nut from the screw attachment.
- Turn the injector nozzle to insert it into the screw attachment ①.
- Place the support on the frame above the injector nozzle.
- Use a 17 mm spanner (WAF 17) to tighten the injector nozzle ②.
- Secure the supports on the right and left using the module's retaining bars ③.

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