A 201

de	Gebrauchsanweisung Injektorwagen
da	Brugsanvisning Injektorvogn
en	Operating instructions Injector unit
es	Instrucciones de manejo Carros inyectores
fi	Käyttöohje Suorasuihkuvaunu
fr	Mode d'emploi Chariot à injection
hr	Upute za uporabu Injektorska kolica
it	Istruzioni d'uso Carrello iniettore
nl	Gebruiksaanwijzing injectorwagen
no	Bruksanvisning injektorvogn
pt	Instruções de utilização Carro injetor
sl	Navodila za uporabo Injektorski voziček
sv	Bruksanvisning Injektorvagn

de	3
da	18
en	32
es	47
fi	62
fr	77
hr	92
it	107
nl	122
no	137
pt	151
sl	166
sv	181

Warnungen

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Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

This mobile unit can be used to reprocess machine-reprocessible medical devices in a Miele washer-disinfector. Follow the operating instructions for your washer-disinfector as well as the instructions issued by the medical device manufacturer on how to reprocess their items using a machine.

Internal cleaning and disinfection to an appropriate standard is only possible with instruments which can be sufficiently dismantled and with adequate use of injector nozzles, irrigation sleeves and adapters. The reprocessing result of the inside of instruments must be checked as necessary using a test other than a visual check.

The A 201 mobile injector unit is designed for reprocessing hollow anaesthetic and intensive care instruments and materials.

Special holders and injector nozzles are used for securing breathing tubes, breathing bags, Guedel tubes and endotracheal tubes, breathing masks and bellows.

The mobile injector unit is designed for the most commonly used silicone breathing tubes. The holder on the left hand side of the unit must be replaced before cleaning breathing tubes of a different diameter or material.

The mobile injector unit is supplied with one holder for paediatric breathing tubes and one for concertina breathing tubes.


The mobile unit can be used in the following washer-disinfectors:

- | | |
|---------------|--------------|
| - PWD 8682 | - PG 8582 |
| - PWD 8682 CD | - PG 8582 CD |
| - PWD 8692 | - PG 8592 |
| - DS 5010 | |
| - DS 5011 | |
| - DS 5020 | |

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessible medical devices are referred to as “load items” if they are not more closely defined.

Queries and technical problems

If you have any queries or technical problems that you would like to discuss please contact Miele. Contact details can be found at the back of the Operating Instructions for your dishwasher or the washer-disinfector or go to www.miele-professional.com.

- Mobile injector unit A 201, height 434 mm, width 531 mm, depth 546 mm.

Items supplied

- 5 x E 362 blind stoppers
- 1 x E 430/1, mesh tray, height 40 mm, width 150 mm, depth 445 mm
- 1 x E 431, injector nozzle for bellows
- 1 x E 432, holder for 3 - 4 concertina breathing tubes
- 2 x E 433, holder for 3 - 4 silicone breathing tubes (fitted)
- 1 x E 434, holder for 3 - 4 paediatric breathing tubes
- 6 x E 466, injector nozzle for breathing bag, length 330 mm, Ø 8 mm
- 15 x E 496, injector nozzle for intubation materials such as Guedel tubes or endotracheal tubes, length 120 mm, Ø 4 mm

Optional accessories

The following and other accessories are available to order from Miele:

- A 810/1 cover net, metal frame with replaceable silicone net A 811/1, 200 x 220 mm
- E 427, module insert for 6 laryngoscopes, height 92 mm, width 210 mm, depth 134 mm
- E 454, injector nozzle for trocar sleeves Ø 10 to 15 mm, with height adjustment clip, 150 mm long, Ø 8 mm
- E 464 holder for injector nozzle E 454 for hollow instruments Ø 20 mm, with height adjustment spring, length 65 mm, Ø 13 mm
- E 497, injector nozzle for single- and double-tubed larynx masks, length 195 mm, with side pipe, length 109 mm
- E 498, injector nozzle for single- and double-tubed larynx masks, length 238 mm, with side pipe, length 120 mm

Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in land-fill sites.

en - Warnings and safety notes

To avoid the risk of accidents and damage to this mobile unit please read these instructions carefully before using it for the first time. Keep these instructions in a safe place where they are accessible to users at all times.

 Please also read the operating instructions for your dishwasher or the washer-disinfector and pay particular attention to the Warning and Safety instructions.

► This mobile unit is only for use in the areas of application as defined in the Operating instructions. Components such as the nozzles must only be replaced with Miele accessories or original parts. Any other application, modification or change is unauthorised and could be dangerous.

Miele cannot be held liable for damage caused by unauthorised use or incorrect operation.

► New mobile units and baskets must be washed in the washer-disinfector without a wash load prior to first use.

► Inspect all mobile units, baskets, modules and inserts daily as described in the "Maintenance" section in the Operating instructions for the washer-disinfector.

► Only items which have been declared by their manufacturer as suitable for machine reprocessing may be processed. The manufacturer's specific reprocessing instructions must be observed.

Disposable instruments must not be put into the machine for processing.

► The reprocessing programme used should include thermal disinfection, and the process chemicals used must be suitable for the material quality of the load items.

► To clean and disinfect lumened instruments thoroughly from the inside, they must be attached to the relevant holders or nozzles on the injector manifolds. The connectors must not come loose during a programme sequence.

At the end of every programme sequence, check that the connectors have not come loose.

► The last rinse, as a minimum, must be carried out with DI water.

► Observe permissible disinfecting and drying temperatures to prevent the risk of premature material deterioration. This will be below 85 °C for some elastomer breathing bags and masks.

► It is essential that the drying result is checked at the end of a cleaning programme.

The inside of tubes, in particular must be completely dry to avoid the growth of water-borne bacteria.

The cleaning programme drying time must be adapted where necessary.

Miele cannot be held liable for damage or injury caused by non-compliance with the warning and safety notes.

en - Areas of application

For validation purposes please make sure that you follow the loading instructions given on the template.

Before inserting lumened instruments, please make sure that all hollow sections are free of blockages. Otherwise the suds solution will not be able to access all areas properly.

⚠ Risk of infection due to insufficient disinfection.
Load items that become detached from the irrigation connector during reprocessing will not be sufficiently disinfected on the inside. Any items that become detached from the irrigation connector during reprocessing must be reprocessed again.

To check before loading the machine and before starting a programme

- Check that fittings such as irrigation sleeves and injector nozzles are securely screwed into position.


⚠ Make sure that all screw connectors are fitted with injector nozzles, adapters, irrigation sleeves or blind stoppers to ensure that all fittings in use are supplied with sufficient standardised pressure. Damaged fittings such as injector nozzles, adapters and irrigation sleeves must not be used.

Fittings not equipped with wash items do not have to be replaced blind stoppers.

- Check that the mobile unit is correctly docked to the water path in the chamber.


Replacing holders

The injector unit can be adapted to take commonly used silicone tubes. If breathing tubes with different diameters or of different materials are used, the E 432 holder on the left hand side of the unit must be replaced with a different holder. The mobile unit is supplied with a variety of holders. Additional holders are available from Miele.

- Loosen the knurled nut on the holder.
- Remove the holder.
- Insert the holder you want to use and secure with the knurled nut.


Breathing tubes

- When attaching the breathing tubes, from the coil insert start with the last row of holders and work from the inside outwards.


To do this:

- Fit the breathing tubes onto the nozzles and press the cone 1 to 2 cm downwards to tension the spring.

- Fit each breathing tube between two retainer hoops. The breathing tubes have been attached correctly when the springs on the nozzles are under tension.


The breathing tubes must be secured by the hoops to prevent them being dislodged from the spring-loaded fittings by the force of the water pressure during the cleaning process.

- Arrange the breathing tubes parallel to each other along the spiral.


⚠ The breathing tubes must not hang down through the spiral as water would remain in the folds.

Breathing bags

- Fit the breathing bag onto the E 466 injector nozzle.
- Make sure the opening of the breathing bag fits securely over the plastic cone.
This will prevent the breathing bags from sticking together during re-processing.
The breathing bag is attached correctly if the plastic cone cannot dislodge itself from the opening of the breathing bag.

Breathing masks


- Remove the stopper from breathing masks which have a closable air cushion. Expel as much air as possible and then replace the stopper.

Concertina bellows


To reprocess a concertina bellows, one of the E 466 injector nozzles must be replaced with a special E 431 injector nozzle for concertina bellows which has a larger support ring.

To do this:

- Use a 7 mm spanner (WAF 7) to unscrew one of the E 466 injector nozzles.


- Insert the E 431 injector nozzle into the free screw attachment and screw tight with the 7 mm spanner (WAF 7) .


- Connect the concertina bellows to the E 431 injector nozzle.

Guedel tube

- Attach the Guedel tube to an E 496 injector nozzle and secure with the spring clip.

Endotracheal tube

- Attach the endotracheal tube to an E 496 injector nozzle and secure with the spring clip.
- Place the balloon lead with the valve behind the spring clip.

Larynx masks


A special combination of nozzles (optional extra) is required for reprocessing larynx masks. These nozzles can be purchased from Miele.

For reprocessing double-tubed larynx masks, a combination of the E 497 and E 498 injector nozzles is required. The nozzles are supplied with their own sets of operating instructions.

Single-tubed larynx masks can be reprocessed either with the E 497 and E 498 injector nozzles or with a combination of the E 454 and E 464 nozzles.


Installation of the E 454 and E 464

- Use a 7 mm spanner (WAF 7) to unscrew one of the E 466 injector nozzles.


- Remove the height adjustment clip from the E 454 injector nozzle.
- Screw the E 454 injector nozzle into the E 464 holder.
- Screw the two nozzles combined into the free screw holder and tighten with the 7 mm spanner.

Application


- Remove as much air as possible from the air cushion.
- Attach the larynx mask to the nozzle combination.
- Attach the balloon lead with the valve behind the spring clip.

Laryngoscopes


Up to 6 laryngoscopes of almost all sizes can be cleaned and disinfected using the E 427 insert. The insert is supplied with its own set of operating instructions.

Small parts

- Reprocess small instruments, e.g. Magill forceps and accessories such as Y-pieces, couplings or connectors in the E 430/1 mesh tray.
- Lightweight load items should be secured with a cover net such as the A 810/1 to prevent them blocking the spray arms.
- Place the mesh tray at the back of the mobile unit.

Test point for measuring water pressure

The left hand holder on the front injector manifold can be used as the test point for measuring water pressure for performance checks and for validation purposes in accordance with EN ISO 15883.

- To measure water pressure replace the existing fixture with a Luer Lock adapter, e.g. E 447.

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