A 204

de	Gebrauchsanweisung Injektorwagen
da	Brugsanvisning Injektorvogn
en	Operating instructions Injector unit
es	Instrucciones de manejo Carros inyectores
fi	Käyttöohje Suorasuihkuvaunu
fr	Mode d'emploi Chariot à injection
hr	Upute za uporabu Injektorska kolica
it	Istruzioni d'uso Carrello iniettore
nl	Gebruiksaanwijzing injectorwagen
no	Bruksanvisning injektorvogn
pt	Instruções de utilização Carro injetor
sl	Navodila za uporabo Injektorski voziček
sv	Bruksanvisning Injektorvagn

de	4
da	15
en	26
es	37
fi	48
fr	59
hr	70
it	81
nl	92
no	103
pt	114
sl	127
sv	138

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Notes about these instructions	26
Intended use	27
Queries and technical problems	27
Items supplied	28
Optional accessories	29
Disposal of the packing material	29
Warnings and safety notes	30
Areas of application.....	31
To check before loading the machine and before starting a programme	31
Loading ophthalmic instruments	32
Loading lumened instruments	32
Mesh tray E 142 for ophthalmic instruments	34
Mesh trays and cassettes from other manufacturers	35
Testing lumened instruments within the framework of quality assurance tests.....	35
Test point for measuring water pressure	35

Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

This mobile unit can be used to reprocess machine-reprocessible medical devices in a Miele washer-disinfector. Follow the operating instructions for your washer-disinfector as well as the instructions issued by the medical device manufacturer on how to reprocess their items using a machine.

Internal cleaning and disinfection to an appropriate standard is only possible with instruments which can be sufficiently dismantled and with adequate use of injector nozzles, irrigation sleeves and adapters. The reprocessing result of the inside of instruments must be checked as necessary using a test other than a visual check.

The A 204 mobile injector unit is intended for the reprocessing of ophthalmic instruments only.

This mobile injector unit has two levels.

The top level is equipped with different connections for reprocessing lumened instruments, e.g. for rinsing- and suction hand pieces and cannulae. Instruments are fitted in silicone holders and fittings in the rack in the mobile injector unit and secured to hose connectors.

The lower level can be loaded with inserts and mesh trays for reprocessing instruments without lumen. The mesh insert in the lower level must not be removed.

Narrow lumened instruments, e.g. Sauter cannulae or rinsing- and suction hand pieces, require the suds solution to be filtered. To do this the injector manifolds in the top level each have to be fitted with a Miele A 800 tubular filter. The tubular filter is supplied with its own fitting instructions.


The mobile unit can be used in the following washer-disinfectors:

- | | |
|---------------|--------------|
| - PWD 8682 | - PG 8582 |
| - PWD 8682 CD | - PG 8582 CD |
| - PWD 8692 | - PG 8592 |
| - DS 5010 | |
| - DS 5011 | |
| - DS 5020 | |

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessible medical devices are referred to as “load items” if they are not more closely defined.

Queries and technical problems

If you have any queries or technical problems that you would like to discuss please contact Miele. Contact details can be found at the back of the Operating Instructions for your dishwasher or the washer-disinfector or go to www.miele-professional.com.

- A 204 Mobile injector unit with two levels, height 365 mm, width 531 mm, depth 542 mm
- with spray arm magnets for spray arm monitoring

Items supplied

- 1 x E 476, 50 holders for micro-instruments, Ø 4 to 8 mm, for use in mesh trays with a 5 mm grid
- 1 x E 477, 20 holders for micro-instruments, for use in mesh trays with a 5 mm grid
- 11 x E 447, Luer Lock adapters ♀ for Luer Lock ♂
- 11 x E 517, Luer Lock adapters ♂ for Luer Lock ♀
- 6 x E 790, 4 Luer-Lock ♀ / Luer-Lock ♀ adapters
- 5 x E 791, 4 Luer-Lock ♀ / Luer ♂ adapters
- 11 silicone hoses with Luer Lock ♂ adapters for Luer Lock ♀, length 160 mm

Optional accessories

The following and other accessories are available to order from Miele:

- A 800, reusable tubular filter with a filtration grade of ≥ 0.07 mm, 405 mm long, $\varnothing$ 22.5 mm
- E 142, DIN mesh tray with a 5 mm grid and 2 hinged handles, height 45/55 mm, width, 255 mm, depth 480 mm
- E 479, 50 holders for $\varnothing$ 4 mm micro instruments, for use in mesh trays with a 5 mm grid

Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in land-fill sites.

en - Warnings and safety notes

To avoid the risk of accidents and damage to this mobile unit please read these instructions carefully before using it for the first time. Keep these instructions in a safe place where they are accessible to users at all times.

⚠ Please also read the operating instructions for your dishwasher or the washer-disinfector and pay particular attention to the Warning and Safety instructions.

► This mobile unit is only for use in the areas of application as defined in the Operating instructions. Components such as the nozzles must only be replaced with Miele accessories or original parts. Any other application, modification or change is unauthorised and could be dangerous.

Miele cannot be held liable for damage caused by unauthorised use or incorrect operation.

► New mobile units and baskets must be washed in the washer-disinfector without a wash load prior to first use.

► Inspect all mobile units, baskets, modules and inserts daily as described in the "Maintenance" section in the Operating instructions for the washer-disinfector.

► Only items which have been declared by their manufacturer as suitable for machine reprocessing may be processed. The manufacturer's specific reprocessing instructions must be observed.

Disposable instruments must not be put into the machine for processing.

► Ophthalmic instruments should only be reprocessed in a machine **specifically dedicated for this application**. This way there is no risk of residues from other disciplines being deposited in narrow lumened hollow instruments.

Synthetic cover nets, e.g. Miele cover nets A 2 or A 3, must not be used in this machine.

► When reprocessing ophthalmic instruments do not dispense rinsing agent.

► The reprocessing programme used should include thermal disinfection, and the process chemicals used must be suitable for the material quality of the load items.

► The last but one and the final rinse should be carried out with fully demineralised water wherever possible.

Miele cannot be held liable for damage or injury caused by non-compliance with the warning and safety notes.

⚠ Ophthalmic instruments should only be reprocessed in a washer-disinfector **specifically dedicated for this application**. This way there is no risk of residues from other disciplines being deposited in narrow lumened hollow instruments.
Synthetic cover nets, e.g. Miele A2 or A3 cover nets, must not be used in this washer-disinfector.

⚠ When reprocessing ophthalmic instruments do not dispense rinsing agent.

To check before loading the machine and before starting a programme

- Check that fittings such as irrigation sleeves and injector nozzles are securely screwed into position.

⚠ Make sure that all screw connectors are fitted with injector nozzles, adapters, silicone tubes with Luer-Lock connectors or blind stoppers to ensure that all rinsing fittings in use are supplied with sufficient standardised pressure.
Damaged fittings such as injector nozzles, adapters and silicone tubes with Luer-Lock connectors must not be used.

Rinsing positions not in use do not have to be fitted with blind stoppers.

- Check that spray arms are rotating freely.
- Check spray arm jets for blockages.
- Check that the mobile unit is correctly docked to the water path in the chamber.

Loading ophthalmic instruments

For validation purposes please make sure that you follow the loading instructions given on the template.

⚠ Risk of injury caused by load.

There is a risk of injury when loading and unloading due to possible sharp edges, rims or pointed ends.


To minimise the risk of injury, loading should take place from the back forwards and unloading should take place in the reverse order.

⚠ Before connecting hollow instruments, please make sure that all hollow sections and cannulae are free of blockages. Otherwise the cleaning solution will not be able to access all areas properly. The cannulae and other hollow instruments should be flushed through with a spray of water as soon as possible after use on patients, preferably whilst still in the theatre. This is to prevent residual blood or medication causing a blockage in the narrow lumen and endangering internal cleaning.

⚠ Narrow lumen instruments require the suds solution to be filtered prior to use. The injector manifolds should each be fitted with a A 800 tubular filter for this purpose. The tubular filter is supplied with its own instructions.


⚠ Instruments must not touch each other during reprocessing.

Loading lumened instruments


- Cannulae with Luer lock connectors should be screwed into the appropriate Luer lock adapters in the top level of the mobile injector unit.

Connect cannulae with Luer connectors to the Luer lock adapter in the injector unit using a connection piece e.g. Luer lock ♀ / Luer ♂.

- Screw long hollow instruments to flexible hose connectors using ♂ Luer-Lock adapters.
- When processing hollow instruments with several connectors, such as rinsing and suction hand pieces, each connector needs to be secured to a hose.
- Hollow instruments with a ♂ connector should be connected to the hose using an E 790 Luer lock ♀ / Luer lock ♀ connector.

- Fit hollow instruments securely into the silicone holders on the rack in the injector unit.
- Space the silicone holders in the rack according to the size of the instruments.
- Unscrew hollow instruments after successful reprocessing.

Risk of infection due to insufficient disinfection.

Load items that become detached from the irrigation connector during reprocessing will not be sufficiently disinfected on the inside. Any items that become detached from the irrigation connector during reprocessing must be reprocessed again.

en - Areas of application


Mesh tray E 142 for ophthalmic instruments

Instruments must not touch each other during reprocessing.


- Light, delicate instruments without lumen can be secured in the mesh tray using E 476 or E 479 silicone holders and E 477 silicone connectors.

If mesh trays are used in the lower level of the unit, E 479 silicone holders have to be used in order to hold the instruments securely.

To ensure that all surfaces are cleaned and rinsed thoroughly:

- Hinged instruments, e.g. scissors, clamps and forceps should be opened up and then placed in the tray.
- Diamond knives must be opened up and then placed in the tray.

- Straight instruments should also be secured using E 477 silicone connectors to prevent the risk of them slipping out during reprocessing.

Mesh trays and cassettes from other manufacturers

Non-Miele mesh trays, cassettes and other containers can be used for reprocessing instruments in Miele washer-disinfectors if they have been risk assessed and specified as suitable for such use and if reprocessing procedures have been successfully validated with their use.

Testing lumened instruments within the framework of quality assurance tests

⚠ Risk of chemical burns from alkaline residues.
Even minimal amounts of alkaline detergent residues can cause chemical burns to the retina.
After machine reprocessing, check the pH value of the residual water in the instruments.

- Residual water must be expelled from lumened instruments using medical compressed air. This procedure also verifies that lumens are clear.
- Check the pH value of the water droplets that come out of the instruments using indicator paper with increments of ≤ 0.5 . The pH value should be between 5 and 8.

If the test result is negative


- In the event of a blockage or reduced clearance, remove the offending particles or blockage manually.

Then, even if the pH value deviates from neutral:

- Reprocess the lumened instruments.

If the pH is still not neutral after reprocessing, contact Miele to arrange for further analysis of the problem.

Test point for measuring water pressure

The test point on the front right can be used for measuring water pressure for performance checks and for validation according to EN ISO 15883.

- To measure water pressure replace the existing fixture with a Luer Lock adapter, e.g. E 447.

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