



<b>de</b>	Dampfkondensator-Standard (SC)	<b>en</b>	Standard steam condenser (SC)	<b>cs</b>	Standardní kondenzátor par (SC)
<b>da</b>	Dampkondensator standard (SC)	<b>el</b>	Βασικός συμπυκνωτής ατμού (SC)	<b>en-US</b>	Standard steam condenser (SC)
<b>es</b>	Condensador de vapor estándar (SC)	<b>fr</b>	Condenseur vapeur standard (SC)	<b>hu</b>	Standard gőzkondenzátor (SC)
<b>it</b>	Condensatore di vapore standard (SC)	<b>nl</b>	Standaard voor dampcondensator (SC)	<b>no</b>	Standard dampkondensator (SC)
<b>pt</b>	Condensador de vapor standard (SC)	<b>ru</b>	Пароконденсатор стандартный (SC)	<b>sv</b>	Ångkondensator – standard (SC)


3

4

2

**5****6****de**

### Dampfkondensator-Standard (SC)

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## en

### Standard steam condenser (SC)

**Reason:** The conversion kit is used to install the standard steam condenser (SC) option.

Due to the weight and type of activities involved in the following work, a second person is required.

#### Parts required

No.	Mat. no.	Designation
1	10928990	Standard steam condenser PLW/PWD 86xx (SC)

An on-site extractor must be present for all connection types. The extractor must either be activated on site or can be activated via the AM-1 electronic module (connection module) of the appliance.

The condensate which forms in the wash chamber is directed into the steam condenser. There, the condensate is guided along the heat exchanger pipes and condenses as pure water. This water runs into the drain via the collector and the intermediate container.

Warm residual air is extracted via an external on-site extractor.

Cold water is channelled through the heat exchanger pipes. The cold water can be introduced via a normal cold-water connection and then directed to the drain via a free flow path, the collector and the intermediate container.

Various connection types for the water supply are available for the standard steam condenser (SC) option.

**If the steam condenser is connected to an on-site cooling circuit, the following points must be observed:**

- Suitable connection valves must be provided on site.
- The Miele solenoid valves cannot be used.
- The free flow path is not required if an on-site cooling circuit is used.

**In order to carry out this work, the appropriate Technical Service Documentation (TSD) is required.**

The warning and safety instructions given in the applicable model-specific TSD must be complied with.

Work must be carried out in accordance with the instructions given in the applicable model-specific TSD.

### Standard steam condenser (SC) – wiring diagram

#### Intake:

Intake connection DK1<sub>in</sub> is connected to DK2<sub>out</sub> via a hose bridge, see Fig. 1. This enables the cold water to flow through both heat exchanger circuits.

Intake connection DK2<sub>in</sub> is connected to cold-water intake valves (Y22), see Fig. 1. The double intake valve (Y22) increases the possible throughput.

An on-site cooling circuit can also be connected as an option. If you choose to do so, please observe the specific requirements for connecting to an on-site cooling circuit.

#### Drain:

Drain connection DK1<sub>out</sub> is connected to the intermediate container via a free flow path (DVGW) and the collector. The water then drains via the intermediate container, see Fig. 1.

Drain connection DK1<sub>out</sub> can also be connected to an on-site cooling circuit as an option. The free flow path is not required with this option.

Drain connection DK2<sub>out</sub> is connected to DK1<sub>in</sub> via a hose bridge. This enables the cold water to flow through both heat exchanger circuits, see Fig. 1.

### Standard steam condenser (SC) 2 circuits – wiring diagram

#### Intake:

Intake connection DK1<sub>in</sub> is connected to one of the outlets on the cold-water intake valve (Y22), see Fig. 2.

Intake connection DK2<sub>in</sub> is connected to the other outlet on the cold-water intake valve (Y22), see Fig. 2.

An on-site cooling circuit can also be connected as an option. If you choose to do so, please observe the specific requirements for connecting to an on-site cooling circuit.

**Drain:**

Drain connection DK1<sub>out</sub> is connected to the intermediate container via a free flow path (DVGW) and the collector. The water then drains via the intermediate container, see Fig. 2.

Drain connection DK2<sub>out</sub> is connected to the intermediate container via a free flow path (DVGW) and the collector. The water then drains via the intermediate connector, see Fig. 2.

Drain connections DK1<sub>out</sub> and DK2<sub>out</sub> can also be connected to an on-site cooling circuit as an option. The free flow path is not required with this option.

**Installing the steam condenser**

- ✚ Fit the frame for the top-box panelling.
- ✚ Fit the corrugated hose (1) on the water intake hood (3), see Fig. 3.
- ✚ Secure the corrugated hose (1) with a hose clip (2), see Fig. 3.
- ✚ Loosen the holding plate (13) fixing screws on the steam condenser, see Fig. 5.
- ✚ Pull the holding plates (13) of the steam condenser outwards as far as they will go, see Fig. 5.
- ✚ Tighten the holding plate (13) fixing screws on the steam condenser, see Fig. 5.
- ✚ Press down the corrugated hose in the direction of the water intake hood and hold. Slide the steam condenser into the frame.
- ✚ Insert the steam condenser into the support rail (3) on the frame, see Fig. 4.
- ✚ Insert the steam condenser into the frame until the connection for the corrugated hose is directly above the corrugated hose.
- ✚ Slide the corrugated hose over the steam condenser connection.
- ✚ Check that the holes and threads for the fixing screws are aligned.
- ✚ Fit the vent hose onto the top of the steam condenser and secure with a hose clip.
- ✚ Tighten the fixing screws.

**Fitting the water intake and drainage (standard, 1 circuit)****NOTE**

Incorrectly laid hoses and cables.

Lift-up sliding door may tear off hoses or cables.

- ✚ Lay all hoses and cables so that they cannot be disturbed by the movement of the lift-up sliding doors.
- ✚ Secure all hoses and cables with cable ties or Velcro tape at suitable points.
- ✚ Fit the cold-water intake hose (KW) on the connector (9) on the intake valve, see Fig. 5.
- ✚ Fit the Y connection hose between the intake valve (9) and connection DK2<sub>in</sub> (12) and secure with hose clips, see Fig. 5.
- ✚ Fit the hose bridge between DK2<sub>out</sub> (11) and DK1<sub>in</sub> (10) and secure with hose clips, see Fig. 5.
- ✚ Fit the hose between DK1<sub>out</sub> (7) and the free flow path intake (3) and secure with hose clips, see Fig. 5.
- ✚ Fit the hose (14) between the free flow path output and the intake of the intermediate container for the drain and secure with hose clips, see Fig. 5.
- ✚ Connect the condensate hose (15) to the condensate output of the steam condenser and the intermediate container for the drain and secure with hose clips, see Fig. 5.

**Connecting the steam condenser to an electrical supply****NOTE**

Incorrectly laid hoses and cables.

Lift-up sliding door may tear off hoses or cables.

- ✂ Lay all hoses and cables so that they cannot be disturbed by the movement of the lift-up sliding doors.
- ✂ Secure all hoses and cables with cable ties or Velcro tape at suitable points.

- ✂ Connect the minor wiring harness to the drying units between the intake valves and the connector (X5-4) and position so that the cables are not laid over hot surfaces or sharp edges.
- ✂ Secure the minor wiring harness to the frame of the top-box panelling so that the cables do not come into contact with the lift-up sliding door.
- ✂ Connect the steam condenser connection plug to the connector (X5-4), see Fig. 6.

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## en-US

### Standard steam condenser (SC)

**Reason:** This conversion kit is for installing the optional standard steam condenser (SC).

A second person is recommended for this procedure due to the weight and type of work.

#### Parts required

Quantity	Mat. no.	Designation
1	10928990	Standard steam condenser (SC) for PLW/PWD 86xx

All connection types will require an on-site ventilation fan. The fan can be internally controlled by the machine's AM-1 electronic (connection module), or externally activated.

Any condensate in the wash chamber is channeled out through the steam condenser. From there, the condensate is fed to the heat exchanger, where it condenses into pure water. This water goes through the sump and intermediate containers into the drain.

Any residual hot air is drawn out by an external on-site ventilation fan.

Cold water is fed through the heat exchanger. The cold water can be supplied by a standard cold-water connection and then fed into the drain via a free water path and the sump/intermediate containers.

There are different water connection types for the optional standard steam condenser (SC).

**If the steam condenser is connected to an on-site cooling circuit, the following must be taken into consideration:**

- Suitable water valves must be provided on site.
- The Miele solenoid valves cannot be used.

- The free water path eliminates the need for an on-site cooling circuit.

**The TSD will be needed for all service work**

The warning and safety instructions in the model-specific TSD must be observed.

Service work must be carried out following the instructions in the model-specific TSD.

**Standard steam condenser (SC) connection diagram****Intake:**

Intake connection DK1<sub>in</sub> is connected to DK2<sub>out</sub> via a hose bridge; see Fig. 1. This enables the cold water to flow through both heat exchanger circuits.

Intake connection DK2<sub>in</sub> is connected to cold-water valve Y22; see Fig. 1. Double-inlet valve Y22 boosts performance.

Alternatively, an on-site cooling circuit can be connected. If doing so, please observe the special requirements for connecting an on-site cooling circuit.

**Drain:**

Drain hose DK1<sub>out</sub> is connected to the intermediate containers via a free water path (DVGW) and the sump, where drainage takes place; see Fig. 1.

Alternatively, drain hose DK1<sub>out</sub> can be connected to an on-site cooling circuit. This option eliminates the free water path.

Drain hose DK2<sub>out</sub> is connected to DK1<sub>in</sub> via a hose bridge. This enables the cold water to flow through both heat exchanger circuits; see Fig. 1.

**Standard steam condenser (SC) 2-cooling-circuit connection diagram****Intake:**

Intake connection DK1<sub>in</sub> is connected to cold-water valve Y22; see Fig. .

Intake connection DK2<sub>in</sub> is connected to the other output of cold-water valve Y22; see Fig. 2.

Alternatively, an on-site cooling circuit can be connected. If doing so, please observe the special requirements for connecting an on-site cooling circuit.

**Drain:**

Drain hose DK1<sub>out</sub> is connected to the intermediate containers via a free water path (DVGW) and the sump, where drainage takes place; see Fig. 2.

Drain hose DK2<sub>out</sub> is connected to the intermediate containers via a free water path (DVGW) and the sump, where drainage takes place; see Fig. 2.

Alternatively, drain hoses DK1<sub>out</sub> and DK2<sub>out</sub> can be connected to an on-site cooling circuit. This option eliminates the free water path.

**Steam condenser installation**

- ✚ Install the framework for the outer casing.
- ✚ Install the corrugated hose (1) on top of the water intake (3); see Figure 3.
- ✚ Secure the corrugated hose (1) with the hose clamp (2); see Fig. 3.
- ✚ Remove the screws securing the steam condenser mounting bracket (13); see Fig. 5.
- ✚ Pull the steam condenser mounting bracket (13) out as far as it will go; see Fig. 5.
- ✚ Reinstall the steam condenser mounting bracket (13) screws; see Fig. 5.
- ✚ Press the corrugated hose down towards the top of the water intake, hold it in place, and push the steam condenser into the framework.
- ✚ Push the steam condenser into the support rail (3) on the framework; see Fig. 4.
- ✚ Push the steam condenser into the framework until the corrugated hose sits directly underneath its connection on the steam condenser.

- ↪ Connect the corrugated hose to the steam condenser.
- ↪ Check that the holes and threads for the retaining screws are on top of each other.
- ↪ Install the vent hose on top of the steam condenser and secure it with a hose clamp.
- ↪ Tighten the retaining screws.

#### Water intake and drain installation (standard, 1 circuit)

**NOTE**

Incorrectly laid hoses and wires.

Doors can tear hoses and wires.

- ↪ Lay hoses and wires so that they cannot be disturbed by the movement of the doors.
  - ↪ Secure all hoses and wires with cable ties or Velcro tape at suitable points.
- ↪ Install the cold-water hose (KW) on the inlet valve (9); see Fig. 5.
  - ↪ Install the Y-piece between the inlet valve (9) and the DK2<sub>in</sub> connection (12) and secure with hose clamps; see Fig. 5.
  - ↪ Install the hose bridge between DK2<sub>out</sub> (11) and DK1<sub>in</sub> (10) and secure with hose clamps; see Fig. 5.
  - ↪ Install a hose between DK1<sub>out</sub> (7) and the free water path intake (3) and secure with hose clamps; see Fig. 5.
  - ↪ Install a hose (14) between the free water path outlet and the intake from the intermediate container for drainage and secure with hose clamps; see Fig. 5.
  - ↪ Connect the condensate hose (15) to the condensate outlet on the steam condenser and the intermediate container and secure with hose clamps; see Fig. 5.

#### Steam condenser electrical connections

**NOTE**

Incorrectly laid hoses and wires.

Doors can tear hoses and wires.

- ↪ Lay hoses and wires so that they cannot be disturbed by the movement of the doors.
  - ↪ Secure all hoses and wires with cable ties or Velcro tape at suitable points.
- ↪ Connect the auxiliary wiring harness between the inlet hoses and the adapter (X5-4) on the drying unit and lay the wires so that they won't come into contact with hot surfaces or sharp edges.
  - ↪ Secure the auxiliary wiring harness to the framework so that its wires will not come into contact with the doors.
  - ↪ Connect the steam condenser to the adapter (X5-4); see Fig. 6.

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