Programmname **ROBOTVARIO X•tra** **Materialnummer** 10235070

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Programme name **ROBOTVARIO X•tra** **Material number** 10235070

Programme header	Conductivity/ Water intake	Last block but one	Maximum value [µS/cm]
		Last block	Maximum value [µS/cm]
	Conductivity/ Water drainage	Last block but one	Maximum value [µS/cm]
			Repeat
		Last block	Maximum value [µS/cm] 20
			Repeat
	Temperature compensation [%/°C]	2.164	
	Max. temp. 1st intake [°C]		
	Water intake time limit	Normal	
	Report process from [Block]	Main wash I	
	Maximum temperatures	Wash temperature [°C]	
		Drying temperature [°C]	
	Safety query	off	

Rinse through	Valve 1 cold [l]	Valve 2 hot [l]
	Valve 3 AD [l]	Valve 4 boiler [l]

Block	Pre wash I	Main wash I	Main wash II	Main wash III	Chemical disinfection	Rinse I	Neutralisation	Rinse II	Rinse III	Rinse IV	Final rinse	Thermal disinfection II
Parameter												
Block type	Main wash	Main wash	Main wash			Main wash		Main wash				Disin
Water intake	cold	AD	No			AD		AD				AD
+/- Water [l]	0	0	0			- 2		- 3				- 3
Dispensing phase 1	Dispensing temperature [°C]	45										
	Holding time before dispensing [s]											
	Dispenser pump	1										
	Dosage concentration [%]	0.5										
Dispensing phase 2	Dispensing temperature [°C]											
	Holding time before dispensing [s]											
	Dispenser pump											
	Dosage concentration [%]											
Dispensing phase 3	Dispensing temperature [°C]											
	Holding time before dispensing [s]											
	Dispenser pump											
	Dosage concentration [%]											
Temperature [°C]	10	50	55			10		55			93	
Holding time [s]	60	600	600			120		180			300	
A0 value												
Spray arm monitoring	On-stop	On-stop	On-stop			On-stop		On-stop				On-stop
Cooling down [l]	0	0	0			0		0				0
Water drainage	Standard	No	Standard			Standard		Standard				Standard
Water drain ext. contact												
Pause	Water intake											
	Water volume [l]											
	Time TA [s]											
	TA											
Print conductivity												On

Drying	Time 1 [min]	15
	Temperature 1 [°C]	110
	Time 2 [min]	
	Temperature 2 [°C]	
	Cooling down [min]	

Programme name	ROBOTVARIO X•tra	Material number	10235070
Machine type	PG8527/28	Date	2014-08-22