A 315/1

cs	Návod k obsluze Injektorový modul
el	Οδηγίες χρήσης μονάδα με εγχυτήρες
en	Operating instructions Injector module
hu	Használati utasítás Injektormodul
ja	取扱説明書 インジェクターモジュール
ko	사용 설명서 인젝터 모듈
pl	Instrukcja użytkowania Moduł iniekcyjny
ro	Instrucțiuni de utilizare Modul de injectie
ru	Инструкция по эксплуатации Инжекторный модуль
sk	Návod na použitie Injektorový modul
tr	Kullanım kılavuzu Enjektör modülü
zh	操作说明书 注射器模块

cs	4
el	18
en	33
hu	47
ja	61
ko	75
pl	89
ro	103
ru	117
sk	133
tr	147
zh	161

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en - Contents

Notes about these instructions	33
Queries and technical problems	33
Definition of terms	33
Intended use	34
IMPORTANT for the U.S.	34
Items supplied	35
Scope of delivery A 315/1.....	35
Disposal of the packing material	35
Optional components	36
Warnings and safety notes	37
Areas of application	38
Screwing in nozzles.....	38
Securing injector nozzles	38
Test point for measuring wash pressure	38
Programme selection.....	38
Injector modules.....	39
Inserting an injector module	39
Removing an injector module	39
Preparing the load	40
To check before loading the machine and before starting a programme	40
Reprocessing hollow instruments.....	40
A 803 holder for transmission instruments	41
Installation	41
Reprocessing transmission instruments.....	42
AUF 1 holder for transmission instruments	42
Optional components	42
Installation	43
Reprocessing transmission instruments.....	43
Replacing filter plates.....	43
A 814 holder for the tips of air scalers and piezo scalers	44
Fitting E 499 for suction and rinsing cannulae	45
Irrigation sleeves	45

Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

Queries and technical problems

In the event of queries or technical problems, please contact Miele. Contact details can be found at the end of the operating instructions for your cleaning machine or at www.miele.com/professional.

Definition of terms

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessable medical devices are referred to as “load items” if they are not more closely defined.

This module can be used to reprocess machine-reprocessable medical devices in special Miele washer-disinfectors. Follow the operating instructions for your washer-disinfector as well as the instructions provided by the medical device manufacturer on how to reprocess their items by machine.

The A 315/1 injector module is suitable for reprocessing dental and ENT instruments by machine.

The module has 2 injector rails with a total of 22 holders for nozzles and other irrigation connectors for reprocessing lumened instruments.

This module is suitable for use in the following:

- A 151 lower basket

Special reprocessing programmes are required for reprocessing instruments from the above areas of application by machine. These programmes are installed ex-works or can be installed, created or activated by the Miele Customer Service Department on the following machines:

- | | |
|---------------|------------------|
| - PWA 8672 | - PG 8581 |
| - PWD 8682 | - PG 8582 |
| - PWD 8682 CD | - PG 8582 CD |
| - PWD 8692 | - PG 8591 |
| - DS 5010 | - PG 8592 |
| - DS 5011 | - RID-100 teon |
| - DS 5020 | - RID-200 teon + |

The module can be used with or without an additional upper basket; e.g., as an addition to the A 105/1 upper basket.

If the module is being used with an upper basket with injector nozzles, a special programme with higher wash pressure must be selected.

- If the upper basket with injector manifolds is being used in combination with the A 102, the Vario TD Dental can continue to be used.
- The Vario TD Dental + programme is required for reprocessing dental instruments.
- The Vario TD ENT + is required for reprocessing ENT instruments.


Special fittings such as A 803 or AUF 1 are needed for reprocessing transmission instruments used in dental applications.

Disposable items must not be reprocessed.

IMPORTANT for the U.S.

Not all sections in this operating manual reflect the Intended Use for the USA. Please refer to the operating manual of the Miele Dental Washer-Disinfector or contact Miele Service Professional at +1 800-991-9380 with application questions.

Scope of delivery A 315/1

- A 315 injector module with 2 nozzle manifolds and a total of 22 connections for irrigation connectors
 - Height 145 mm, width 492 mm, depth 542 mm
 - Injector rail vertical clearance when using upper basket A 105/1 135 mm (± 30 mm)
- A 838, tool set, consisting of 1 x WAF 9 spanner and 1 x T20 offset screwdriver

Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in land-fill sites.

Optional components

- A 800, reusable tubular filter with a filtration grade of ≥ 0.07 mm, 405 mm long, $\varnothing$ 22.5 mm
- A 801, threaded adapter
- A 803, holder for dental transmission instruments, height 76.5 mm, width 40 mm, depth 32 mm
- A 813, two-piece adapter with 4 screw threads for various irrigation connectors, height approx. 48 mm, $\varnothing$ 45 mm, 4 nozzles with silicone hoses and 3 blind stoppers
- A 814, holder with screw cap for reprocessing the tips of air scalers and piezo scalers, height approx. 41 mm, $\varnothing$ 16 mm
- A 817, irrigation sleeve, length 65 mm, $\varnothing$ 11 mm, cap with opening $\varnothing$ 6 mm
- A 818, 12 caps for irrigation sleeves, opening $\varnothing$ 2 mm
- A 836, nozzle set for e.g. dental applications, consisting of 22 x A 833, injector nozzle with aperture on the side, length 110 mm, $\varnothing$ 4 mm
- A 837, irrigation sleeve set for e.g. ENT applications, consisting of 22 x A 817, irrigation sleeve for instruments $\varnothing$ 4–8 mm, length 67 mm, $\varnothing$ 11 mm, cap with aperture $\varnothing$ 6 mm
- AUF 1, two-piece adapter for transmission instruments (dental handpieces, for example), height approx. 48 mm, $\varnothing$ 45 mm; silicone adapter (ADS) required for reprocessing
- E 336, irrigation sleeve, length 121 mm, $\varnothing$ 11 mm
- E 442, irrigation sleeves for instruments $\varnothing$ 4 - 8 mm, 121 mm long, $\varnothing$ 11 mm, cap with aperture $\varnothing$ 6 mm
- E 443, irrigation sleeves for instruments $\varnothing$ 8 - 8.5 mm, 121 mm long, $\varnothing$ 11 mm, cap with aperture $\varnothing$ 10 mm
- E 445, 12 caps for irrigation sleeves, aperture $\varnothing$ 6 mm
- E 446, 12 caps for irrigation sleeves, aperture $\varnothing$ 10 mm
- E 447, Luer Lock Adapter ♀ for Luer Lock ♂
- E 499 for reprocessing suction and rinsing cannulae with a variety of connector types as well as Veress cannulae.
- FP 20, 20 filter plates, $\varnothing$ 30 mm

Further components are available from Miele as optional extras.


To avoid the risk of accidents and damage to this module please read these instructions carefully before using it for the first time. Keep these instructions in a safe place where they are accessible to users at all times.

 Please also read the operating instructions for your dishwasher or the machine and pay particular attention to the Warning and Safety instructions.

- ▶ The module is approved solely for the applications specified in the “Areas of application” chapter of these operating instructions. Components such as nozzles may only be replaced with Miele components or genuine original spare parts.
- ▶ New load carriers must be cleaned in the washer-disinfector without a load prior to first use.
- ▶ Only ever place empty, unloaded modules in the baskets. Check that they are correctly engaged before loading them. Modules must be completely emptied before they are removed. Placing loaded modules in baskets or removing them can cause damage to the load items.
- ▶ Inspect all mobile units, baskets, modules and inserts daily as described in the “Maintenance” section in the Operating instructions for the washer-disinfector.
- ▶ For your own safety, wear gloves and, if necessary, a surgical mask and protective goggles when sorting the instruments and unloading the machine. To avoid recontamination, wear clean gloves when unloading the machine.
- ▶ Load from back to front and unload from front to back.
- ▶ To clean and disinfect lumened instruments thoroughly from the inside, they must be attached to the relevant holders or nozzles on the injector manifolds. The connectors must not come loose during a programme sequence. At the end of every programme sequence, check that the connectors have not come loose.
- ▶ Any hollow items that have become disconnected from their injector nozzles, adapters or holders during processing must be re-processed.
- ▶ The last rinse, as a minimum, must be carried out with DI water.


Screwing in nozzles

Securing injector nozzles

- Screw the injector nozzles into the required positions.
- Use the 9 mm spanner (WAF 9) to tighten the injector nozzles.

Test point for measuring wash pressure

The holder at the top of the right-hand injector manifold is used as the test point for measuring wash pressure. This test point can be used for measuring wash pressure within the framework of performance and validation testing in accordance with EN ISO 15883.

- To measure the water pressure, replace the existing fixture with a Luer Lock adapter, e.g. E 447.

Programme selection

The module can be used with or without an additional upper basket; e.g., as an addition to the A 105/1 upper basket.

If the module is being used with an upper basket with injector nozzles, a special programme with higher wash pressure must be selected.


- The Vario TD Dental + programme is required for reprocessing dental instruments.
- The Vario TD ENT + is required for reprocessing ENT instruments.

Injector modules

⚠ Only ever place empty, unloaded modules in the baskets. Check that they are correctly engaged before loading them. Modules must be completely emptied before they are removed. Placing loaded modules in baskets or removing them can cause damage to the load items.

Inserting an injector module

Insert the injector module into the A 151 lower basket.

- The triangular bracket secures the module to prevent it from moving to the side. When inserting the module, make sure that the bracket does not protrude beyond the edge of the lower basket ①.
- Hook the module into the struts of the lower basket ②.
- Then engage the front bracket of the module by pushing it into the front cross strut of the basket ③.
- Carefully push the basket into the cleaning machine to check that the water connection is in the correct position.

Removing an injector module

⚠ Do not hold the basket by the injector module when removing it from the cleaning machine. The module could come loose or the cleaning machine's water connection could be damaged.

- Release the module from the locking mechanism and lift it out of the basket.

Preparing the load

For validation purposes please make sure that you follow the loading instructions given on the template.

⚠ Risk of injury caused by load.

There is a risk of injury when loading and unloading due to possible sharp edges, rims or pointed ends.

To minimise the risk of injury, loading should take place from the back forwards and unloading should take place in the reverse order.

⚠ Risk of infection due to insufficient disinfection.

Load items that become detached from the irrigation connector during reprocessing will not be sufficiently disinfected on the inside. Any items that become detached from the irrigation connector during reprocessing must be reprocessed again.

To check before loading the machine and before starting a programme

- Check that fittings such as irrigation sleeves and injector nozzles are securely screwed into position.

⚠ Make sure that all screw connectors are fitted with injector nozzles, adapters, irrigation sleeves or blind stoppers to ensure that all fittings in use are supplied with sufficient standardised pressure. Damaged fittings such as injector nozzles, adapters and irrigation sleeves must not be used.

Fittings not equipped with wash items do not have to be replaced blind stoppers.

- Is the module correctly docked to the water supply in the basket?

Reprocessing hollow instruments

Before inserting lumened instruments, please make sure that all hollow sections are free of blockages. Otherwise the suds solution will not be able to access all areas properly.


- Attach the lumened instruments to the nozzles.
- At the end of the programme, the insides of the instruments must be dried with sterile compressed air before being cleaned and sterilised as appropriate in accordance with the manufacturer's instructions. Please observe national regulations.

Where necessary, the nozzles supplied can be replaced with additional components, such as holders for dental transmission instruments, irrigation sleeves, etc. Further examples are given on the following pages.

Due to their length, surgical suction devices used in dental applications cannot be reprocessed using the A 315 module if an upper basket is being used in the cleaning machine at the same time. For this reason, surgical suction devices must be reprocessed in an upper basket with an injector manifold, such as the A 105/1.

⚠ The wash water used for reprocessing narrow-lumened instruments needs to be filtered. The injector manifolds can each be fitted with an A 800 tubular filter for this purpose. The tubular filter is supplied with its own operating instructions.

Cleaning the A 800 tubular filter

The tubular filter must be cleaned **at least once a week**.

We recommend cleaning the tubular filter after every 10 programme sequences.

- Follow the cleaning instructions in the operating instructions for the tubular filter.

A 803 holder for transmission instruments

A 803 holders must be retro-fitted for reprocessing dental transmission instruments. The holders are suitable for almost all standard transmission instruments.


- Remove a nozzle from the injector manifold and replace it with an A 803 holder.

When the A 803 holders are used, the wash water must be filtered. The injector manifold must be fitted with an A 800 tubular filter. The tubular filter is supplied with its own operating instructions. Follow in particular the instructions on how to perform regular cleaning on the tubular filter.

Installation

Reprocessing transmission instruments

- Check that internal sections of transmission instruments are free of blockages.

- Insert the transmission instrument into the holder. The instrument is then held in position by the spring clamp.
- After reprocessing pull the transmission instrument out of its holder. The spring clamp will automatically spring back to its starting position.
- After reprocessing, the insides of the transmission instruments must be dried with sterile compressed air before being cleaned and sterilised as appropriate in accordance with the manufacturer's instructions. Please observe national regulations.

Before using transmission instruments again following reprocessing, a function check must be carried out, e.g. by spraying into a basin, to ensure they are clear.

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
AUF 1 holders are designed for reprocessing dental transmission instruments.

Optional components

- A 801, threaded adapter
- ADS 1, silicone adapter, white, for holding transmission instruments with a Ø of approx. 20 mm
- ADS 2, silicone adapter, green, for holding transmission instruments with a Ø of approx. 16 mm
- ADS 3, silicone adapter, red, for holding transmission instruments with a Ø of approx. 22 mm
- FP 20, 20 filter plates, porosity 2, Ø 30 mm

Installation

An A 801 adapter is required.

- Unscrew the upper part ② from the AUF 1 holder.
- Fit a silicone adapter (ADS) ① in its place which has the same diameter as the transmission instrument being reprocessed.

⚠ Fitting a transmission instrument into a silicone adapter that is too big can cause damage to the filter plate.

- Remove an injector nozzle from the injector manifold and replace it with the A 801 adapter ⑤.
- Screw the lower part of the holder ④ into the adapter.
- Insert a filter plate ③ into the lower part.

⚠ Do not insert a filter plate into the holder if the injector manifold has been fitted with an A 800 tubular filter.

- Screw the upper part into the lower part.

Reprocessing transmission instruments

Before loading check that the silicone adapters are present and correctly fitted.

- Wipe away any deposits from the area around the silicone adapters that may be present depending on local water conditions using a damp cloth.
- Check that internal sections of transmission instruments are free of blockages.
- Fit the instruments into their holders using the respective silicone adapters.
- After reprocessing, the insides of the transmission instruments must be dried with sterile compressed air before being cleaned and sterilised as appropriate in accordance with the manufacturer's instructions. Please observe national regulations.


Before using transmission instruments again following reprocessing, a function check must be carried out, e.g. by spraying into a basin, to ensure they are clear.

en - Areas of application

Replacing filter plates

⚠ Filter plates must be replaced **approx. every 2 weeks** or after **20 reprocessing cycles** have been completed.

- Unscrew the holder completely.


- Dismantle the holder and dispose of the filter plate ③.
- Rinse out the lower part ④ thoroughly.
- Insert a new filter plate.


⚠ Damaged or used filter plates must not be inserted under any circumstances.

- Screw the upper and lower parts of the holder back together and fit the holder back onto the injector manifold.

A 814 holder for the tips of air scalers and piezo scalers

A 814 holder

A 814 holder with A 813 adapter

The A 814 holder enables the tips of air scalers and piezo scalers to be reprocessed. To increase capacity, the holder can be combined with the A 813 adapter. The holder and the adapter each come with their own operating instructions.

Fitting E 499 for suction and rinsing cannulae

The E 499 is for reprocessing rinsing and suction cannulae, including suction breakers. It comes with its own instructions.

Irrigation sleeves

Irrigation sleeves are suitable for reprocessing hollow instruments and instruments with the following external diameters:

- A 817, irrigation sleeve, length 65 mm, Ø 11 mm, cap with opening Ø 6 mm
- E 336, irrigation sleeve, length 121 mm, Ø 11 mm
- E 442, irrigation sleeves for instruments Ø 4 - 8 mm, 121 mm long, Ø 11 mm, cap with aperture Ø 6 mm
- E 443, irrigation sleeves for instruments Ø 8 - 8.5 mm, 121 mm long, Ø 11 mm, cap with aperture Ø 10 mm

The aperture can be varied using optional replacement caps:

- A 818, 12 caps for irrigation sleeves, opening Ø 2 mm
- E 445, 12 caps for irrigation sleeves, aperture Ø 6 mm
- E 446, 12 caps for irrigation sleeves, aperture Ø 10 mm

The diameter of the instruments and the irrigation sleeve apertures or caps must be compatible with each other, e.g. caps for irrigation sleeves with an aperture diameter of Ø 6 mm for dental surgical suction cannulae.

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