


Miele



Installation plan

PG 8582 CD

PG 8583 CD

To avoid the risk of accidents or damage to the machine, it is **essential** to read these instructions as well as the service documents before it is installed, commissioned and used for the first time.

en - GB, IE

M.-Nr. 10 399 151

Installation notes

| | |
|--------------------------------------|--|
| Installation requirements | <p>This machine must be installed by a suitably qualified person with the appropriate electrical and plumbing qualifications in accordance with the installation instructions supplied.</p> <p>This machine must be installed in accordance with all applicable standards and guidelines, including legal requirements and health and safety regulations.</p> <p>This machine must be commissioned and operatives trained in its use by Miele Service or by an approved Miele Service Partner only.</p> |
| Surrounding area | <p>Condensate can build up in the area surrounding the machine. Any cabinetry and fixtures in the room must be suitable. If the machine is installed under a countertop (built-under model), the protective foil supplied must be stuck to the underside of the worktop and a stainless steel panel must be installed above the door opening to provide protection from the steam. The stainless steel panel can be ordered from Miele.</p> |
| Preventing a build-up of heat | <p>During the active drying process (TA) there is a risk of heat building up. A build-up of heat at the back of the machine can cause the machine casing and electronic components to overheat. A build-up of heat can also cause condensate to settle on adjoining surfaces. This can lead to a reduced machine lifespan and cause damage to surrounding units and worktops.</p> <ul style="list-style-type: none">- To prevent a build-up of heat, ensure that there is sufficient space for air to circulate behind base units.- Leave a minimum safety gap of 10 mm for air flow between a built-under machine and the worktop above it.- Do not seal any gaps between base units or between base units and the machine.- If necessary, fit ventilation grilles in adjoining units and in the worktop. |
| Cooling the exhaust air | <p>Activate “Air cooling” if ambient room temperature exceeds 30 °C. This option is found in the System settings menu under “Further settings”.</p> |
| Plumbing | <p>This machine can be connected directly to a cold water supply or a hot water supply using the supplied double check valves, or to a demineralised water supply. Use the Y-piece supplied to connect the inlet hoses to the machine (KW) and to the steam condenser (DK). Alternatively, you may also install an additional cold water supply and connect the steam condenser to this supply.</p> <p>If access to a hot water supply is not available, a Y-piece must be used to connect both inlet hoses (KW cold water and WW hot water) to the cold water supply.</p> <p>The stopcocks must be easily accessible.</p> |
| Electrical connection | <p>Connection to the electrical supply must be carried out in accordance with local and national safety regulations. The power cord must be protected from the risk of thermal damage.</p> |

It is recommended to make electrical connection via a **plug and socket** so that electrical safety checks can be carried out easily. For **hard-wired** machines, connection must be via a main switch to be provided on site, which must completely isolate the machine from the power supply with a contact gap of at least 3 mm. The plug and socket as well as the main switch must be easily accessible after the machine has been installed.

Faulty components must only be replaced by genuine Miele original spare parts. Only when these parts are fitted can the safety standards of the machine be guaranteed. If the connection cable is faulty it must only be **replaced by a Miele approved service** technician to protect the user from danger.

Equipotential bonding and earthing

For added safety the machine should be protected with a **residual current device** with trip current of 30 mA.

Equipotential bonding should be carried out if required. A screw connection for equipotential bonding is provided at the back of the machine. Equipotential bonding and earthing must be completed before the machine is commissioned.



Liquid agents: position of external containers

The liquid agent container for external dispensing must be positioned either **next to or underneath** the machine only. The container may be placed on the floor or in an adjacent cabinet. Do not place the container on top of or above the machine. Make sure that the dispensing hose does not get kinked or trapped.


Communication modules

Ethernet and **RS232 modules** are available as optional extras. They are not supplied with the machine. The connection box for these should be installed near the machine for transferring and printing process data. These must be installed and connected in accordance with **IEC 60950**.



Figures



Dimensions, standard version



Dimensions, extended depth



EL Electrical connection

KW Cold water connection

DK Cold water connection, steam condenser

NW Network and printer connection (optional)


WW Hot water connection

AW Drain connection


VE Demineralised water connection

PA Equipotential bonding

Use the enclosed Y-piece to connect the steam condenser to the cold water supply using the supplied double check valve. In the UK the (DK) cold (KW) and hot (WW) water inlet hoses must be connected to the potable water using the supplied double check valves.




A - A



Connection points

A - A




on machines with increased depth


| | | | |
|----|---|----|--------------------------------|
| EL | Electrical connection | WW | Hot water connection |
| KW | Cold water connection | AW | Drain connection |
| DK | Cold water connection, steam condenser | VE | Demineralised water connection |
| NW | Network and printer connection (optional) | PA | Equipotential bonding |

Figures


Machine versions



Standard undercounter model



with lid
(optional)



with increased depth and lid
(optional)

Electrical connection

| | |
|--|-------------------------|
| Voltage (standard version) | 3N AC 400/50 |
| Total connected load | 9.3 kW |
| Fuse rating | 3 x 15-16 A |
| Connection cable, cross-section min. | 5 x 2.5 mm ² |
| Length of mains connection cable (H05(07)RN-F) | 1.8 m |
| | |
| Voltage (convertible) | AC 230/50 |
| Total connected load | 6.3 kW |
| Fuse rating | 1 x 30-32 A |
| Connection cable, cross-section min. | 3 x 4 mm ² |
| Length of mains connection cable (H05(07)RN-F) | 1.9 m |

Cold water connection

(UK: the supplied double check valve must be fitted)

| | |
|--|-------------|
| Maximum incoming temperature | 20 °C |
| Max. permitted water hardness | 12.6 mmol/l |
| | 70 °dH |
| Recommended flow pressure | 200 kPa |
| Minimum flow pressure with extended water intake | 100 kPa |
| Maximum pressure | 1,000 kPa |
| Flow rate | 7.5 l/min |
| On-site threaded union in accordance with DIN 44991 (flat sealing) | 3/4 Inch |
| Length of cold water inlet hose | 1.7 m |
| Length of steam condenser inlet hose | 1.7 m |

Hot water

(UK: the supplied double check valve must be fitted)

| | |
|--|-------------|
| Maximum incoming temperature | 60 °C |
| Max. permitted water hardness | 12.6 mmol/l |
| | 70 °dH |
| Recommended flow pressure | 200 kPa |
| Minimum flow pressure with extended water intake | 40 kPa |
| Maximum pressure | 1,000 kPa |
| Flow rate | 7.5 l/min |
| On-site threaded union in accordance with DIN 44991 (flat sealing) | 3/4 Inch |
| Length of hot water inlet hose | 1.7 m |

Technical data

Demin. water

| | |
|--|-----------|
| Maximum incoming temperature | 60 °C |
| Recommended flow pressure (AD pressure-resistant) | 200 kPa |
| Minimum flow pressure with extended water intake | 30 kPa |
| Maximum pressure (AD pressure resistant) | 1,000 kPa |
| Flow rate | 7.5 l/min |
| On-site threaded union in accordance with DIN 44991 (flat sealing) | 3/4 Inch |
| Length of demin. water inlet hose | 1.7 m |

PG 8583 CD: machines with a booster pump for demineralised water (ADP) are only suitable for connection/use with a non-pressurised system (PG 8583 CD only).

| | |
|--|-----------|
| Minimum flow pressure with extended water intake | 8.5 kPa |
| Maximum pressure (AD unpressurised) | 60 kPa |
| Machine connection (dia. x l) | 6 x 30 mm |

Waste water

| | |
|---|------------|
| Drainage water temperature | 93 °C |
| Drain hose length, standard | 1.4 m |
| Drain hose, max. drainage length | 4.0 m |
| Max. drain pump head height from bottom edge of machine | 1.0 m |
| Max. transient flow rate | 16 l/min |
| On-site sleeve for drain hose (dia. x length) | 22 x 30 mm |

Machine feet

| | |
|---|----------|
| Height adjustment at the front | 0 - 8 mm |
| Diameter of machine feet | 35 mm |
| Insert for machine foot, size of thread | M 8 |

Machine data

Technical data

| | |
|--|----------|
| Building-under height | 820 mm |
| Height including lid | 835 mm |
| Width | 898 mm |
| Depth | 598 mm |
| Depth incl. 10 cm extended building-in depth | 698 mm |
| Door height | 622 mm |
| Net weight | 98 kg |
| Floor load in operation | 2,000 N |
| Min. access width, incl. transport pallet | 1,170 mm |
| Min. access depth incl. transport pallet | 740 mm |
| Min. access height, incl. transport pallet | 1,000 mm |
| Noise level in dB (A), sound pressure LpA during cleaning and drying phases | < 70 dB |

Technical data

Heat dissipation rate to installation site

| | |
|--------------------------------------|----------|
| From heat radiation during operation | 0.35 kWh |
| From load whilst unloading | 0.40 kWh |

Installation requirements

| | |
|--|---------|
| Permissible ambient temperature | 40 °C |
| Max. relative humidity up to 31 °C | 80 % |
| Rel. humidity, declining proportionally to 40 °C | 50 % |
| Max. installation altitude above sea level | 2,000 m |

United Kingdom

Miele Co. Ltd., Fairacres, Marcham Road
Abingdon, Oxon, OX14 1TW
Professional Sales, Tel: 0845 365 6608
E-mail: professional@miele.co.uk
Internet: www.miele.co.uk/professional

Australia

Miele Australia Pty. Ltd.
ACN 005 635 398, ABN 96 005 635 398
1 Gilbert Park Drive, Knoxfield, VIC 3180
Tel: 1300 731 411
Internet: www.miele.com.au/professional
E-mail: info@miele-professional.com.au

China

Miele (Shanghai) Trading Ltd.
1-3 Floor, No. 82 Shi Men Yi Road
Jing' an District, 200040 Shanghai, PRC
Tel: +86 21 6157 3500, Fax: +86 21 6157 3511
E-mail: info@miele.cn, Internet: www.miele.cn

Hong Kong, China

Miele (Hong Kong) Ltd.
41/F - 4101, Manhattan Place
23 Wang Tai Road, Kowloon Bay, Hong Kong
Tel: (852) 2610 1025, Fax: (852) 3579 1404
Email: customerservices@miele.com.hk
Website: www.miele.hk

India

Miele India Pvt. Ltd.
Ground Floor, Copia Corporate Suites
Plot No. 9, Jasola, New Delhi - 110025
Tel: 011-46 900 000, Fax: 011-46 900 001
E-mail: customercare@miele.in, Internet: www.miele.in

Ireland

Miele Ireland Ltd.
2024 Bianconi Ave., Citywest Business Campus, Dublin 24
Tel: (01) 461 07 10, Fax: (01) 461 07 97
E-Mail: info@miele.ie, Internet: www.miele.ie

Manufacturer: Miele & Cie. KG, Carl-Miele-Straße 29, 33332 Gütersloh, Germany

Malaysia

Miele Sdn Bhd
Suite 12-2, Level 12
Menara Sapura Kencana Petroleum
Solaris Dutamas No. 1, Jalan Dutamas 1
50480 Kuala Lumpur, Malaysia
Phone: +603-6209-0288
Fax: +603-6205-3768

New Zealand

Miele New Zealand Limited
IRD 98 463 631
8 College Hill
Freemans Bay, Auckland 1011, NZ
Tel: 0800 464 353
Internet: www.miele.com.au/professional
E-mail: info@miele-professional.com.au

Singapore

Miele Pte. Ltd.
29 Media Circle, #11-04 ALICE@Mediapolis
Singapore 138565
Tel: +65 6735 1191, Fax: +65 6735 1161
E-Mail: info@miele.com.sg
Internet: www.miele.sg

South Africa

Miele (Pty) Ltd
63 Peter Place, Bryanston 2194
P.O. Box 69434, Bryanston 2021
Tel: (011) 875 9000, Fax: (011) 875 9035
E-mail: info@miele.co.za
Internet: www.miele.co.za

United Arab Emirates

Miele Appliances Ltd.
Showroom 1, Eiffel 1 Building
Sheikh Zayed Road, Umm Al Sheif
P.O. Box 114782 - Dubai
Tel. +971 4 3044 999, Fax. +971 4 3418 852
800-MIELE (64353)
E-Mail: info@miele.ae, Website: www.miele.ae