

ProCare Dent 11 A

Cleaning agent


Use:

Powder, alkaline. For the automated reprocessing of dental instruments in Miele Professional washer-disinfectors.

ProCare Dent 11 A supports the removal of organic soiling such as blood and saliva and prevents re-deposition. Coarse residues from composites, dental cement and dental impression compound should be removed before they have had time to set.

ProCare Dent 11 A is not suitable for cleaning transmission instruments. The use of the *ProCare Dent 10 MA* liquid cleaning agent is recommended for this purpose.

Quality features:

ProCare Dent 11 A complies with European guidelines for medical devices. 

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Product properties:

- Excellent material compatibility.
- Protection against corrosion for instruments and trays made from stainless steel.


Use and dispensing:

ProCare Dent 11 A is formulated for use in processes in Miele Professional washer-disinfectors.

- Dispensing via door dispenser for powder cleaning agents.
- Suitable dosing dispensers should be used.
- Application concentration: 0.375 oz/gal / 3 g/l
- Programme selection: Vario TD
- Minimum dispensing temperature: ≥ 104 °F / 40 °C

An identical quantity of cleaning agent should be used each time, preferably by using the enclosed dispensing spoon.

Programme sequence:



Vario TD process:

- 1 Pre-wash**
Cold water (softened), no heating, no process chemicals
- 2 Main wash**
Cold / hot water (softened), with heating, with *ProCare Dent 11 A*
- 3 Neutralisation**
Cold / hot water (softened), no heating, with *ProCare Dent 30 P*
- 4 Interediate rinse**
Cold / hot water (softened) water, no heating
- 5 Thermal disinfection**
at > 194 °F / 90 °C with 5-minute holding time in the final rinse
- 6 Drying**

The product in solution has to be completely rinsed off with water. No visible cleaning agent residue should be present at the end of the process. In Miele Professional washer-disinfectors, this is guaranteed by the Vario TD programme when properly used (e.g. correct loading). Verification is provided by measuring the increase of the conductivity in the final rinse during commissioning and routine maintenance. Miele recommends the use of fully demineralised water in the final rinse in order to protect the load and the thermal disinfectant against stains and deposits.

Please note!

- Do not mix with other products.
- Reprocessing should comply with all ordinances pursuant to medical device regulations and should be performed with appropriate validated processes.
- Observe recommendations on reprocessing issued by medical device manufacturers with respect to the provisions of DIN EN ISO 17664 (CAN/CSA-Z17664).
- Operating instructions for Miele Professional washer-disinfectors must be observed.
- For commercial applications only.

ProCare Dent 11 A

Cleaning agent


Technical data:

- Bulk density: approx. 1,015 g/l to 1,065 g/l
- pH value (determined using fully demineralised water, at 68 °F / 20 °C) 3 g/l: approx. 11.5

Ingredients according to Regulation (EC) No 648/2004 on detergents:

- < 5 % non-ionic tensides
- > 30 % phosphates

Storage information:

- Storage temperature: 32 °F / 0 °C to 86 °F / 30 °C
- Refit lid firmly after use. The product tends to get lumpy when exposed to dampness which can cause a loss of effectiveness.
- Can be used for 2 years under proper storage conditions.
- Use-by date: See label next to hourglass symbol 

Warning and Safety notes:

For safety and environmental information cf. safety data sheet. These are available from the Miele websites listed below.

Dispose of containers only when empty and sealed. Disposal of product residue: refer to safety data sheet.

All information contained herein is current and correct to our knowledge. This does not however release users from obligations to perform their own trials and experiments. This cannot be deemed as a legal warranty of product characteristics or properties.

Available packaging units:

Designation	Packaging unit
<i>ProCare Dent 11 A</i>	10 kg (22 lb) bucket

Contact information:**Germany**

Miele Vertriebsgesellschaft Deutschland KG
Carl-Miele-Straße 2
33332 Gütersloh
Tel.: +49 800 22 44 644
vertrieb@miele-professional.de
www.miele.de/professional

Canada

Miele Limited
161 Four Valley Drive
VAUGHAN
Ontario L4K 4V8
KANADA
Tel. +1 905 532-2270
customercare@miele.ca
www.miele.ca

USA

Miele Inc.
9 Independence Way
Princeton, NJ 08540
Tel: 800-991-9380
www.miele-professional.com

**Germany**

Chemische Fabrik
Dr. Weigert GmbH & Co. KG
Mühlenhagen 85
20539 Hamburg, Germany
www.drweigert.com