A 202

de	Gebrauchsanweisung Wagen
da	Brugsanvisning Vogn
en	Operating instructions: Mobile unit
es	Instrucciones de manejo Carros
fi	Käyttöohje Vaunu
fr	Mode d'emploi Chariot
hr	Upute za uporabu Kolica
it	Istruzioni d'uso Carrelli
nl	Gebruiksaanwijzing Wagen
no	Bruksanvisning vogn
pt	Instruções de utilização Carro
sl	Navodila za uporabo Voziček
sv	Bruksanvisning Vagn

de	4
da	17
en	29
es	41
fi	53
fr	65
hr	77
it	89
nl	101
no	114
pt	126
sl	139
sv	151

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en - Contents

Notes about these instructions	29
Intended use	30
Queries and technical problems	30
Items supplied	31
Optional accessories	31
Medical devices	31
Laboratory glassware and utensils	32
Disposal of the packing material	32
Warnings and safety notes	33
Medical devices.....	33
Laboratory glassware and utensils.....	33
Areas of application.....	34
To check before loading the machine and before starting a programme	34
Fitting the grille	34
Medical devices.....	35
Loading mesh trays	35
Processing hollow containers	36
Laboratory glassware and utensils.....	37
Microtiter plates and blood analysis plates	38
Test point for measuring water pressure	39

Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

This mobile unit can be used to reprocess machine-reprocessible medical devices, laboratory glassware and utensils in a Miele washer-disinfector. Follow the operating instructions for your washer-disinfector as well as the instructions of the medical device or glassware and utensil manufacturer on how to reprocess their items by machine.

The A 202 mobile unit can accommodate up to 4 DIN mesh trays or can be loaded with a variety of mesh trays and inserts for reprocessing medical devices, laboratory glassware and laboratory utensils.

The mobile unit has two levels. The lower level is equipped with a removable grille. This grille is needed for loading mesh trays in the lower level. It can be removed, if necessary, e.g. for reprocessing large items.

The Vario TD Inst 4trays programme is designed for processing 4 mesh trays with medical devices, e.g. surgical instruments.

The A 202 mobile unit, in combination with the A 103 upper basket, enables up to 6 DIN mesh trays to be processed at the same time. In this case medical devices must be processed using the Vario TD Inst 6trays programme.

The Injector plus programme is required for reprocessing laboratory glassware and utensils.

The mobile unit can be used in the following washer-disinfectors:

- | | |
|---------------|--------------|
| - PLW 8636 | - PG 8582 |
| - PLW 8683 | - PG 8582 CD |
| - PLW 8683 CD | - PG 8583 |
| - PLW 8693 | - PG 8583 CD |
| - PWD 8682 | - PG 8592 |
| - PWD 8682 CD | - PG 8593 |
| - PWD 8692 | |
| - DS 5010 | |
| - DS 5011 | |
| - DS 5020 | |

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessable medical devices, laboratory glassware and utensils are referred to as “items” if they are not more closely defined.

Queries and technical problems

If you have any queries or technical problems that you would like to discuss please contact Miele. Contact details can be found at the back of the Operating Instructions for your dishwasher or the washer-disinfector or go to www.miele-professional.com.

- A 202 Mobile unit with two levels, height 223 mm, width 531 mm, depth 542 mm
- with removable grille in the lower level.
- with spray arm magnets for spray arm monitoring

Items supplied

- 6 x E 362 blind stoppers

Optional accessories

The following and other accessories are available to order from Miele:

- E 142, DIN mesh tray with a 5 mm grid and 2 hinged handles, height 45/55 mm, width, 255 mm, depth 480 mm
- E 146, mesh tray with lid and hinged handles, height 55 mm, width 150 mm, depth 255 mm

Medical devices

- E 448, silicone tube with Luer Lock adapter ♂ for Luer-Lock ♀, 300 mm long
- E 473/2 mesh tray with lid for small items, height 85 mm, width 60 mm, depth 60 mm
- E 476, 50 holders for Ø 4 to 8 mm micro instruments, for use in mesh trays with a 5 mm grid
- E 477, 20 holders for micro-instruments, for use in mesh trays with a 5 mm grid
- E 479, 50 holders for Ø 4 mm micro instruments, for use in mesh trays with a 5 mm grid

en - Items supplied

Laboratory glass-ware and utensils

- E 492, insert with 9 holders for kidney dishes, spacing 49 mm, height 120 mm, width 256 mm, depth 474 mm
- A 11/1, stainless steel retaining frame for upper and lower baskets, width 429 mm, depth 429 mm.
- AK 12, insert for funnels, glass beakers, wide neck flasks, etc., height 67 (127) mm, width 225 mm, depth 442 mm.
- E 104/1, insert for approx. 200 test tubes up to 12 x 105 mm, including lid A 13, height 132 (152) mm, width 200 mm, depth 320 mm.
- E 106, insert with 10 spring hooks, height 175 mm and 16 spring hooks, height 105 mm, for wide neck flasks, measuring cylinders etc., spacing approx. 60 mm, height 186 mm, width 195 mm, depth 430 mm.
- E 109, insert for 21 glass beakers up to 250 ml, 21 x 3 rods, height 155 mm, width 230 mm, depth 460 mm.
- E 136, insert for 56 petri dishes Ø 100 mm, spacing approx. 26 mm, height 145 mm, width 485 mm, depth 445 mm.
- E 402, insert for 44 watch glasses Ø 80 - 125 mm, spacing approx. 15 mm, height 53 mm, width 200 mm, depth 445 mm.
- E 494, insert for 5 loose microtiter plates, height 35 mm, width 205 mm, depth 440 mm.

Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in land-fill sites.

To avoid the risk of accidents and damage to this mobile unit please read these instructions carefully before using it for the first time. Keep these instructions in a safe place where they are accessible to users at all times.

⚠ Please also read the operating instructions for your dishwasher or the washer-disinfector and pay particular attention to the Warning and Safety instructions.

► The mobile unit is approved solely for the applications specified in the “Areas of application” chapter of these operating instructions. Applications that go beyond those listed and use of the mobile unit with additional components from the Miele product range are described in the operating instructions for the components or are to be agreed with Miele.

► Inspect all mobile units, baskets, modules and inserts daily as described in the “Maintenance” section in the Operating instructions for the washer-disinfector.

► Only items which have been declared by their manufacturer as suitable for machine reprocessing may be processed. The manufacturer's specific reprocessing instructions must be observed. Disposable instruments must not be put into the machine for processing.

► Broken glass can result in serious injury when loading or unloading. Broken glass items must not be processed in the machine.

Medical devices

► The reprocessing programme used should include thermal disinfection, and the process chemicals used must be suitable for the material quality of the load items.

► The last rinse, as a minimum, must be carried out with DI water.

► Visual checks should be undertaken to check the cleaning result. Additional protein analysis checks should also be carried out e.g. weekly, on a random sampling basis.

Laboratory glassware and utensils

► The reprocessing result must be checked as necessary using a suitable test method and not just by a visual check.

Miele cannot be held liable for damage or injury caused by non-compliance with the warning and safety notes.

For validation purposes please make sure that you follow the loading instructions given on the template.

To check before loading the machine and before starting a programme

- Check that spray arms are rotating freely.
- Check spray arm jets for blockages.

⚠ Make sure that all screw connectors are fitted with injector nozzles, adapters, irrigation sleeves or blind stoppers to ensure that all fittings in use are supplied with sufficient standardised pressure. Damaged fittings such as injector nozzles, adapters and irrigation sleeves must not be used.

Fittings not equipped with wash items do not have to be replaced blind stoppers.

- Check that the mobile unit is correctly docked to the water path in the chamber.

Fitting the grille

The removable grille must be fitted when loading the lower level with DIN mesh trays. It can be removed, if necessary, to accommodate large items.

- Place the grille in the lower level. The flat side must face forwards.


The vertical clearance for the lower level is 95 mm with a top frame and 135 mm without a frame.

Instead of the top frame, a retaining frame can also be inserted into the lower level, e.g. A 11/1. The retaining frame can be loaded with items such as wide-necked glassware.


⚠ Do not place any retaining frames on the top frame. If you do, the laboratory glassware will not be secured by the struts of the mobile unit and could fall out.

Medical devices

Loading mesh trays

- Hinged instruments must be opened to minimise the risk of surfaces not being accessed during reprocessing.

- To prevent the risk of light instruments and other small items being moved about during reprocessing, place them in closable inserts such as the E 473/2 or secure them in silicone holders such as the E 476, E 477 and E 479. The E 473/2 insert is supplied with its own operating instructions.

Hollow instruments Up to 6 hollow instruments can be reprocessed in the A 202 mobile unit using a Luer lock connector.

- To do so replace the blind stopper on the water pipe with a silicone hose with Luer lock adapter, e.g. the E 448.
- Connect instruments with the silicone hoses and place the instruments in a mesh tray on the upper level of the mobile unit. Secure light instruments with silicone holders such as the E 476, E 477 and E 479.

en - Areas of application

Processing hollow containers

- Hollow containers such as beakers and kidney dishes should be loaded with the open side facing downwards. Lightweight containers should be secured with a cover net e.g. the A 810
- Special inserts should be used if necessary for processing hollow containers. These inserts can be purchased from Miele.

The containers must not hinder the processing of other instruments and must be secure so that they cannot be flipped over by the force of water during processing.

If hollow items are placed with the open side facing upwards, suds solution can collect in them. This can spill and impair cleaning results.

Laboratory glassware and utensils

Wide neck items, e.g. glass beakers, wide neck Erlenmeyer flasks and petri dishes, or cylindrical items, e.g. test tubes can be cleaned and rinsed inside and out using rotating spray arms.

- Arrange wide neck flasks with the open side facing downwards.

Laboratory glassware must not hinder the reprocessing of other items and must not be flipped over by the force of water during the programme. Secure light items with a cover net such as A 2 or re-process them in mesh trays with lids.

If laboratory glassware is placed with the open side facing upwards, suds solution can collect in it. This can spill and impair cleaning results in the washer-disinfector.

- Use special inserts where necessary to reprocess laboratory glassware and utensils. Various inserts for different types of items are available from Miele, for example:

Glass beakers,
wide neck Erlen-
meyer flasks

Glass beakers and wide neck Erlenmeyer flasks can be reprocessed in mesh trays, e.g. E 142, inserts, such as AK 12 or in a retaining frame, e.g. A 11/1.
Alternatively inserts with spring hooks or rods can be used, e.g. E 106 or E 109.

Test tubes

Special quarter segment baskets in a variety of sizes, e.g. E 104/1 can be used to reprocess test tubes.

Petri dishes, watch
glasses

Special inserts are available such as E 136 for processing up to 56 petri dishes or E 402 for up to 44 watch glasses.

Small items

Reprocess small items such as stoppers, lids, spatulas, magnetic stirring rods etc. in mesh trays with lids, e.g. E 146. Alternatively a mesh tray with a sufficiently large cover net could be used.

en - Areas of application

Microtiter plates and blood analysis plates

The E 494 insert is designed for the machine reprocessing of microtiter plates and blood analysis plates.

Plastic microtiter plates and blood analysis plates:

For microtiter plates / blood analysis plates made of polycarbonate, polyacrylic, polystyrene, etc., please note the following:


- Use cleaning agents with hypochlorite or hydrogen peroxide oxidising properties.

⚠ Do not use cleaning agents and defoamers containing surfactants and/or silicone oils.

- Use neutralising agents with a citric acid base.
- Select a cleaning programme with 1 - 2 cold pre-washes.
- The temperature must not exceed 55 °C in the washing and rinsing blocks.
- After the final rinse, manually tip any remaining water out of the cavities or wells before drying.
- The temperature when drying in the washer-disinfector must not exceed 65 °C.
- The temperature when drying in the drying cabinet must not exceed 55 °C.


The instructions of the manufacturer of the microtiter plates or the blood analysis plates should also be observed.

Arrange microtiter plates in the E 494 insert

- Push the microtiter plates into the insert holders with the open sides facing downwards.
- Up to two inserts can be used per level. Do not stack the inserts on top of each other.

Test point for measuring water pressure

Looking at the water connection at the back, the 4th connection on the water pipe can be used as the test point for measuring water pressure for performance checks and for validation according to EN ISO 15883.

- To measure water pressure replace the existing blind stopper with a Luer Lock adapter, e.g. E 447.

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