A 103

de	Gebrauchsanweisung Oberkorb
da	Brugsanvisning Overkurv
en	Operating instructions Upper basket
es	Instrucciones de manejo Cesto superior
fi	Käyttöohje Yläkori
fr	Mode d'emploi Panier supérieur
hr	Upute za uporabu Gornja košara
it	Istruzioni d'uso cesto superiore
nl	Gebruiksaanwijzing bovenrek
no	Bruksanvisning Overkurv
pt	Instruções de utilização Cesto superior
sl	Navodila za uporabo Zgornja košara
sv	Bruksanvisning Överkorg

de	4
da	16
en	28
es	40
fi	52
fr	64
hr	76
it	88
nl	100
no	112
pt	124
sl	136
sv	148

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Notes about these instructions	28
Intended use	29
Queries and technical problems	29
Items supplied	30
Optional accessories	30
Medical devices	30
Laboratory glassware and utensils	30
Disposal of the packing material	31
Warnings and safety notes	32
Medical devices.....	32
Laboratory glassware and utensils.....	32
Areas of application.....	33
To check before loading the machine and before starting a programme	33
Medical devices.....	34
Loading mesh trays	34
Processing hollow containers	35
Laboratory glassware and utensils.....	36
Microtiter plates and blood analysis plates	37
Test point for measuring water pressure	38

Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

Re-usable medical devices, laboratory glassware and laboratory utensils can be reprocessed using this unit in a Miele washer-disinfector. Follow the instructions for your washer-disinfector as well as the instructions of the medical device or glassware and utensil manufacturer on how to reprocess their items by machine.

The A 103 upper basket can accommodate up to 2 DIN mesh trays or can be loaded with a variety of mesh trays and inserts for processing medical devices, laboratory glassware and laboratory utensils.

The A 202 mobile unit, in combination with the A 103 upper basket, enables up to 6 DIN mesh trays to be processed at the same time. In this case medical devices must be processed using the Vario TD Inst 6trays programme.

The Injector plus programme is required for reprocessing laboratory glassware and utensils.


This upper basket can be used in the following washer-disinfectors:

- | | |
|---------------|--------------|
| - PLW 8636 | - PG 8582 |
| - PLW 8683 | - PG 8582 CD |
| - PLW 8683 CD | - PG 8583 |
| - PLW 8693 | - PG 8583 CD |
| - PWD 8682 | - PG 8592 |
| - PWD 8682 CD | - PG 8593 |
| - PWD 8692 | |
| - DS 5010 | |
| - DS 5011 | |
| - DS 5020 | |

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessable medical devices, laboratory glassware and utensils are referred to as “items” if they are not more closely defined.

Queries and technical problems

If you have any queries or technical problems that you would like to discuss please contact Miele. Contact details can be found at the back of the Operating Instructions for your dishwasher or the washer-disinfector or go to www.miele-professional.com.

- A 103 upper basket, height 133 mm, width 528 mm, depth 528 mm, Vertical clearance: 95 mm
- with spray arm magnets for spray arm monitoring

Optional accessories

The following and other accessories are available to order from Miele:

- E 142, DIN mesh tray with a 5 mm grid and 2 hinged handles, height 45/55 mm, width, 255 mm, depth 480 mm
- E 146, mesh tray with lid and hinged handles, height 55 mm, width 150 mm, depth 255 mm

Medical devices

- E 473/2 mesh tray with lid for small items, height 85 mm, width 60 mm, depth 60 mm
- E 476, 50 holders for Ø 4 to 8 mm micro instruments, for use in mesh trays with a 5 mm grid
- E 477, 20 holders for micro-instruments, for use in mesh trays with a 5 mm grid
- E 479, 50 holders for Ø 4 mm micro instruments, for use in mesh trays with a 5 mm grid

Laboratory glass-ware and utensils

- E 494, insert for 5 loose microtiter plates, height 35 mm, width 205 mm, depth 440 mm.

Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in land-fill sites.

en - Warnings and safety notes

Read the operating instructions carefully before using this upper basket.

Keep these operating instructions in a safe place.

 Please also read the operating instructions for your dishwasher or the washer-disinfector and pay particular attention to the Warning and Safety instructions.

- ▶ This upper basket is only for use in the areas of application as defined in the Operating instructions. Components such as injector nozzles must only be replaced with Miele accessories or original parts. Any other application, modification or change is unauthorised and could be dangerous.
Miele cannot be held liable for damage caused by unauthorised use or incorrect operation.
- ▶ Inspect all mobile units, baskets, modules and inserts daily as described in the “Maintenance” section in the Operating instructions for the washer-disinfector.
- ▶ Only items which have been declared by their manufacturer as suitable for machine reprocessing may be processed. The manufacturer's specific reprocessing instructions must be observed.
Disposable instruments must not be put into the machine for processing.
- ▶ Broken glass can result in serious injury when loading or unloading. Broken glass items must not be processed in the machine.

Medical devices

- ▶ The reprocessing programme used should include thermal disinfection, and the process chemicals used must be suitable for the material quality of the load items.
- ▶ The last rinse, as a minimum, must be carried out with DI water.
- ▶ Visual checks should be undertaken to check the cleaning result. Additional protein analysis checks should also be carried out e.g. weekly, on a random sampling basis.

Laboratory glassware and utensils

- ▶ The reprocessing result must be checked as necessary using a suitable test method and not just by a visual check.

Miele cannot be held liable for damage or injury caused by non-compliance with the warning and safety notes.


For validation purposes please make sure that you follow the loading instructions given on the template.

To check before loading the machine and before starting a programme


- Check that the basket is correctly connected to the water supply in the washer-disinfector.
- Check that spray arms are rotating freely.
- Check spray arm jets for blockages.

Medical devices

Loading mesh trays

- Hinged instruments must be opened to minimise the risk of surfaces not being accessed during reprocessing.

- To prevent the risk of light instruments and other small items being moved about during reprocessing, place them in closable inserts such as the E 473/2 or secure them in silicone holders such as the E 476, E 477 and E 479. The E 473/2 insert is supplied with its own operating instructions.

Hollow instruments Up to 6 hollow instruments can be reprocessed in the A 202 mobile unit using a Luer lock connector.

- To do so replace the blind stopper on the water pipe with a silicone hose with Luer lock adapter, e.g. the E 448.
- Connect instruments with the silicone hoses and place the instruments in a mesh tray on the upper level of the mobile unit. Secure light instruments with silicone holders such as the E 476, E 477 and E 479.

Processing hollow containers

- Hollow containers such as beakers and kidney dishes should be loaded with the open side facing downwards. Lightweight containers should be secured with a cover net e.g. the A 810
- Special inserts should be used if necessary for processing hollow containers. These inserts can be purchased from Miele.

The containers must not hinder the processing of other instruments and must be secure so that they cannot be flipped over by the force of water during processing.

If hollow items are placed with the open side facing upwards, suds solution can collect in them. This can spill and impair cleaning results.

Laboratory glassware and utensils

Items with wide necks, e.g. glass beakers or Erlenmeyer flasks can be cleaned and rinsed inside and out using rotating spray arms.

- Arrange wide neck flasks with the open side facing downwards.

Laboratory glassware must not hinder the reprocessing of other items and must not be flipped over by the force of water during the programme.

Secure light items with a cover net such as A 2 or reprocess them in mesh trays with lids.

If laboratory glassware is placed with the open side facing upwards, suds solution can collect in it. This can spill and impair cleaning results in the washer-disinfector.

- Use special inserts where necessary to reprocess laboratory glassware and utensils. Various inserts for different types of items are available from Miele, for example:

Glass beakers,
wide neck Erlen-
meyer flasks
Small items

Glass beakers and wide neck Erlenmeyer flasks are processed in mesh trays, e.g. E 142.

Reprocess small items such as stoppers, lids, spatulas, magnetic stirring rods etc. in mesh trays with lids, e.g. E 146. Alternatively a mesh tray with a sufficiently large cover net could be used.

Microtiter plates and blood analysis plates

The E 494 insert is designed for the machine reprocessing of microtiter plates and blood analysis plates.

Plastic microtiter plates and blood analysis plates:

For microtiter plates / blood analysis plates made of polycarbonate, polyacrylic, polystyrene, etc., please note the following:


- Use cleaning agents with hypochlorite or hydrogen peroxide oxidising properties.

⚠ Do not use cleaning agents and defoamers containing surfactants and/or silicone oils.

- Use neutralising agents with a citric acid base.
- Select a cleaning programme with 1 - 2 cold pre-washes.
- The temperature must not exceed 55 °C in the washing and rinsing blocks.
- After the final rinse, manually tip any remaining water out of the cavities or wells before drying.
- The temperature when drying in the washer-disinfector must not exceed 65 °C.
- The temperature when drying in the drying cabinet must not exceed 55 °C.


The instructions of the manufacturer of the microtiter plates or the blood analysis plates should also be observed.

Arrange microtiter plates in the E 494 insert

- Push the microtiter plates into the insert holders with the open sides facing downwards.
- Up to two inserts can be used per level. Do not stack the inserts on top of each other.

Test point for measuring water pressure

The test point for measuring water pressure for performance checks and for validation according to EN ISO 15883 is located at the side on the water pipe for the spray arm.

- To measure water pressure replace the blind stopper with a Luer Lock adapter, e.g. E 447.

Under no circumstances must wash load items, fixtures etc. be connected to the measuring point. After measuring the access point must be closed off again with the blind stopper.

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