en	Instructions for use / technical description Sterilizer baskets incl. accessories (silicone/plastic)
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusaifus.com . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung / technische Beschreibung Sterilisationskörbe inkl. Zubehör (Silikon/Kunststoff)
fr	Mode d'emploi / description technique Paniers d'autoclave avec accessoires (silicone/plastique)
es	Instrucciones de manejo / descripción técnica Cestas esterilizadoras incl. accesorios (silicona/plástico)
it	Istruzioni per l'uso / descrizione tecnica Cestelli per sterilizzatrice con accessori (silicone/plastica)
pt	Instruções de utilização / descrição técnica Cestos do esterilizador incl. acessórios (silicone/plástico)
nl	Gebruiksaanwijzing / technische beschrijving Sterilisatiemanden incl. accessoires (silicone/kunststof)
da	Brugsanvisning / teknisk beskrivelse Autoklavekurve inkl. tilbehør (silikone/plast)
no	Bruksanvisning / teknisk beskrivelse Steriliseringskurver inkl. tilbehør (silikon/plast)
sv	Bruksanvisning / teknisk beskrivning Steriliseringskorgar inkl. tillbehör (silikon/plast)
fi	Käyttöohjeet / tekninen kuvaus Sterilointikorit, sis. lisävarusteet (silikoni/muovi)
et	Kasutusjuhend / tehniline kirjeldus Steriliseerimiskorvid koos tarvikutega (silikoon/plast)
lv	Lietošanas instrukcija / tehniskais apraksts Sterilizatora grozi, ieskaitot piederumus (silikons/plastmasa)
lt	Naudojimo instrukcija / techninis aprašymas Sterilizatoriaus krepšeliai su priedais (silikoniniai / plastikiniai)
ru	Инструкция по применению / техническое описание Корзины для стерилизации с принадлежностями (силикон/пластик)
cs	Návod k použití / technický popis Sterilizační koše včetně příslušenství (silikon/plast)
pl	Instrukcja użytkowania / opis techniczny Kosze do sterylizacji wraz z akcesoriami (silikon / tworzywo sztuczne)
sk	Návod na používanie / technický opis Sterilizačné koše vrátane príslušenstva (silikón/plast)
hu	Használati utasítás / műszaki leírás Sterilizáló kosarak tartozékokkal (szilikon/műanyag)
sl	Navodila za uporabo / tehnični opis Košare sterilizatorjev z dodatno opremo (iz silikona/plastike)
hr	Upute za upotrebu / tehnički opis Košare za sterilizaciju uklj. dodatnu opremu (silikon/plastika)
ro	Manual de utilizare / descriere tehnică Coșuri pentru sterilizator, incl. accesorii (silicon/plastic)
bg	Инструкции за употреба / техническо описание Кошници за стерилизатор вкл. аксесоари (силикон/пластмаса)
tr	Kullanım kılavuzu / teknik açıklama Aksesuarlar dahil sterilizatör sepetleri (silikon/plastik)
el	Οδηγίες χρήσης / τεχνική περιγραφή Καλάθια αποστειρωτή συμπερ. εξαρτημάτων (σilikόνη/πλαστικό)
br	Instruções de uso / descrição técnica Cestos de esterilização, incluindo acessórios (silicone/plástico)

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.bbraun.com

AESCULAP® – a B. Braun brand

TA016545 2025-03 Change No. AE0064993

1 About this document

Note
General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply to all Aesculap sterilizer baskets incl. accessories made of silicone/plastic.
Note
Instructions for use and further information about B. Braun / AESCULAP products can be found on the B. Braun eIFU website at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

- ⚠ WARNING**
Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.
- ⚠ CAUTION**
Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2 Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended purpose
For intended purpose, see Intended use.

2.1.2 Intended use
Aesculap sterilizer baskets incl. accessories made of silicone/plastic or with silicone/plastic components are used by OR and CSSD personnel and are dimensionally stable containers intended for repeated use. They are used for holding goods to be sterilized, for transport and storage in a sterile barrier system (e.g. sterilization containers, etc.).

2.1.3 Indications
Note
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.
For indications, see Intended use.

2.1.4 Contraindications
No known contraindications.

2.2 Safety information

2.2.1 Clinical user
General safety information
To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:
► Use the product only according to these instructions for use.
► Follow the safety and maintenance instructions.
► Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
► Store any new or unused products in a dry, clean, and safe place.
► Prior to use, check that the product is in good working order.
► Keep the instructions for use accessible for the user.

Note
The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures
The medical professional will make decisions on concrete applicability based on the warranted properties and technical data.

2.2.2 Product-specific safety information
► Only combine Aesculap products with each other.

2.2.3 Sterility
The product is supplied non-sterile.
► Clean the new product after removing its transport packaging and prior to its initial sterilization.

2.3 Application

⚠ WARNING
Risk of injury and/or malfunction!
► Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
► Always carry out a function test prior to each use of the product.
Note
As required, elements made of silicone can be trimmed to have the desired shape.

3 Validated processing procedure

3.1 Safety information

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.
Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the processing of products.
Note
Mechanical processing should be favored over manual cleaning as it gives better and more reliable results.
Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.
Note
Up-to-date information about processing and material compatibility can be found on the B. Braun eIFU site at eifu.bbraun.com
The validated steam sterilization procedure was carried out in the AESCULAP sterile container system.

3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore, the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.
Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.
Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.
Additional drying, if necessary.
Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:
■ Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
■ Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
► Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
► Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to "AKI-Brochures", "Red brochure".

3.3 Service life

Materials for reusable surgical instruments are generally chosen to be suitable for repeated processing. However, it should be noted that each mechanical, chemical and thermal treatment can lead to stress and thus to aging of the material.
The service life of the product is limited by damage, normal wear, type and duration of use, as well as by handling, storage and transport of the product.
Influences of processing using the validated procedure that lead to damage to the product are not known. Careful visual and functional inspection before each use is the best way to detect a product that is no longer functional, see Visual inspection and see Functional test.

3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.5 Cleaning/Disinfection

3.5.1 Product-specific safety information on the processing method
The reprocessing procedure described below applies only to unloaded trays.
► For loaded trays, read and follow the instructions for reprocessing the products in their respective instructions for use.
Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!
► Use cleaning and disinfecting agents according to the manufacturer's instructions.
► Observe specifications regarding concentration, temperature and exposure time.
► Do not exceed the maximum allowable disinfection temperature of 95 °C.
► Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

3.5.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with ultrasound and immersion disinfection	—	see Manual cleaning/disinfecting and subsection: ■ see Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection	—	see Mechanical cleaning/disinfection and subsection: ■ see Mechanical alkaline cleaning and thermal disinfecting

3.6 Manual cleaning/disinfecting

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfectant.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfecting process if necessary.

3.6.1 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/ °F]	t [min]	Conc. [%]	Water quality	Chemical/Note
I	Ultrasonic cleaning	RT (cold)	> 15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ≈ 9*
II	Intermediate rinse	RT (cold)	1	–	D-W	–
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ≈ 9*
IV	Final rinse	RT (cold)	1	–	FD-W	–
V	Drying	RT	–	–	–	–

D-W	Drinking water
FD-W	Fully demineralized water (low-germ, max. 10 CFU / 100 ml, as well as low endotoxin contamination, max. 0,25 endotoxin units/ml)
RT	Room temperature

- Phase I
- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
 - Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
 - Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - Drain any remaining water fully.

- Phase III
- Fully immerse the product in the disinfectant solution.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

- Phase IV
- Rinse/flush the product thoroughly (all accessible surfaces).
 - Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
 - Rinse lumens with an appropriate disposable syringe at least five times.
 - Drain any remaining water fully.

- Phase V
- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.7 Mechanical cleaning/disinfection

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. compliance with EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.7.1 Mechanical alkaline cleaning and thermal disinfecting

Phase	Step	T [°C/ °F]	t [min]	Water quality	Chemical/Note
I	Prerinse	< 25/ 77	3	D-W	–
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: – pH ≈ 13 – less than 5 % anionic surfactant ■ 0,5 % working solution – pH = 11*
III	Intermediate rinse	> 10/ 50	> 1	FD-W	–
IV	Thermal disinfection	90/194	5	FD-W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D-W	Drinking water
FD-W	Fully demineralized water (low-germ, max. 10 CFU / 100 ml, as well as low endotoxin contamination, max. 0,25 endotoxin units/ml)

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning/disinfecting process if necessary.

3.8 Inspection

- Allow the product to cool down to room temperature.
- Dry the product if it is wet or damp.

3.8.1 Visual inspection

- Ensure that all soiling has been removed. In particular, pay attention to hinges and hinge handles.
- If the product is dirty: repeat the cleaning and disinfection process.
- Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- Check the product for missing or faded labels.
- Check the surfaces for rough spots.
- Check the product for burrs that could damage tissue or surgical gloves.
- Check the product for loose or missing parts.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.8.2 Functional test

- Check that the product functions correctly.
- Check for compatibility with associated products.
- Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.9 Steam sterilization

- Check to ensure that the sterilizing agent comes into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- Use validated sterilization process.
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to EN 285 and validated according to ISO 17665
 - Allowed sterilization parameters, see table below
- If several devices are sterilized simultaneously in the same steam sterilizer: Make sure that the maximum allowed load according to the specifications of the manufacturer is not exceeded.

Note
The material resistance of the products is not impaired by a sterilization temperature of up to 134 °C and/or a holding time of 18 minutes.

Allowed sterilization parameters

Sterilization process	T [°C]	Holding time [min]	Drying time (at least recommended) [min]
Steam sterilization (fractionated vacuum process)	134	3 – 18	20

The sterilization of products approved for 134 °C is permissible in the temperature range from 134 °C to 137 °C.

3.10 Storage

- Shelf life depends on the quality of the packaging system or material, the tightness of the seals, and the storage conditions.
- Store sterile products at room temperature in a dust-free, clean, dry, and pest-free environment.
 - Follow the storage instructions provided by the sterile barrier system manufacturer.

3.11 Transport

- Transport and storage must not adversely affect the characteristics of the processed medical device.
- Use appropriate transport systems and aids to prevent damage or recontamination.

4 Technical service

⚠ CAUTION
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.
► Do not modify the product.
► For service and repairs, please contact your national B. Braun/AESCLAP agency.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1601
Fax: +49 7461 16-2887
E-Mail: ats@aesculap.de
Other service addresses can be obtained from the address indicated above.

5 Disposal

⚠ WARNING
Risk of infection due to contaminated products!
► Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

Note
The user institution is obliged to process the product before its disposal, see Validated processing procedure.

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