

ProCare

Dent 11 A

Detergent

Use:

Powder, alkaline. For the automated reprocessing of dental instruments in Miele Professional washer-disinfectors.

ProCare Dent 11 A supports the removal of organic soiling such as blood and saliva and prevents redeposition. Coarse residues from composites, dental cement and dental impression compound should be removed before they have had time to set.

ProCare Dent 11 A is not suitable for cleaning transmission instruments. The use of the *ProCare Dent 10 MA* liquid detergent is recommended for this purpose.

Quality features:

ProCare Dent 11 A complies with European guidelines for medical devices. 

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Product properties:

- Excellent material compatibility
- Protection against corrosion for instruments and trays made from stainless steel.


Use and dispensing:

Miele *ProCare Dent 11 A* is formulated for use in processes in Miele Professional washer-disinfectors.

- Application concentration: 3 g/l via door dispenser for powder detergents
- Program selection: Vario TD

An identical quantity of detergent should be used each time, preferably by using the enclosed dispensing spoon.

Program sequence:



vario TD process:

- 1 Pre-wash** with cold water to remove coarse soiling (e.g. blood).
- 2 Main wash** with *ProCare Dent 11 A* at 130 °F / 55 °C (dispensing temperature 105 °F / 40 °C)
- 3 Neutralization** with *ProCare Dent 30 P*
- 4 Intermediate rinse** with cold water
- 5 Thermal disinfection** at 200 °F / 93 °C with 5-minute holding time in the final rinse
- 6 Drying**

The product in solution is to be completely rinsed off using water. No visible detergent residue should be present at the end of the process. In Miele Professional washer-disinfectors, this is guaranteed by the Vario TD program when properly used (e.g. correct loading). Verification is provided by measuring the increase of the conductivity in the final rinse during commissioning and routine maintenance. Miele recommends the use of fully demineralized water in the final rinse in order to protect the load and the washer-disinfector against stains and deposits.

Please note!

- Do not mix with other products.
- Reprocessing should comply with all ordinances pursuant to medical device regulations and should be performed with appropriate validated processes.
- Observe recommendations on reprocessing issued by medical device manufacturers with respect to the provisions of DIN EN ISO 17664.
- Operating instructions on Miele Professional washer-disinfector must be observed.
- For commercial applications only.

ProCare

Dent 11 A

Detergent


Technical data:

- Bulk weight: approx. 1015 to 1065 g
- pH-value (determined using fully demineralized water at 68 °F / 20 °C) 3 g/l: approx. 11.5

Ingredients according to Regulation (EC) 648/2004 on detergents:

- < 5 % non-ionic tensides
- > 30 % phosphates

Storage information:

- Storage temperature: 32 °F to 86 °F / 0 °C to 30 °C
- Refit lid firmly after use. The product tends to get lumpy when exposed to dampness which can cause a loss of effectiveness.
- Can be used for 2 years under proper storage conditions
- Use-by date: See label next to hourglass symbol 

Warning and Safety notes:

For safety and environmental information cf. safety data sheet. These are available from Miele websites listed below.

Dispose of containers only when empty and sealed. Disposal of product residue: refer to safety data sheet.

All information contained herein is current and correct to our knowledge. This does not however release users from obligations to perform their own trials and experiments. This cannot be deemed as a legal warranty of product characteristics or properties.

Available packaging unit:

Designation	Packaging unit	Mat. no.
ProCare Dent 11 A	10 kg (22.0 lbs) bucket	12 186 990

Contact information:

Germany

Miele & Cie. KG
Carl-Miele-Straße 29,
33332 Gütersloh
Tel: +49 (0) 8002244644
vertrieb@miele-professional.de
www.miele.de/professional

USA

Miele Inc.
9 Independence Way
Princeton, NJ 08540
Tel: 800-991-9380
www.miele-professional.com



Germany

Chemische Fabrik
Dr. Weigert GmbH & Co. KG
Mühlenhagen 85
20539 Hamburg, Germany
www.drweigert.de