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**en / fr - CA**

## Installation instructions

**Please read and observe the instructions in the operating instruction manual and the service documents regarding installation and setting up of this machine.**

### Installation requirements

This machine must be installed by a suitably qualified person with **the appropriate electrical and plumbing qualifications** in accordance with the installation instructions supplied.

This machine must be installed in accordance with all applicable standards and guidelines, including legal requirements and health and safety regulations.

The machine must be **commissioned** and operatives trained in its use by Miele Service or by an approved Miele Service Partner only.

### Environmental requirements

Condensate can build up in the area surrounding the machine. Any furniture and fittings in the room must therefore be suitable for purpose. If the machine is fitted under a suitable worktop (built-under model), a stainless steel panel must be fitted above the door aperture to protect the worktop from moisture damage. A suitable stainless steel cover plate is available to order from Miele.

### Electrical connection

Connection to the electrical supply must be carried out in accordance with local and national safety regulations. The connection cable must be protected from the risk of thermal damage.

The machine is supplied with electrical cord and plug as this will make it easier to carry out the electrical safety tests during commissioning and after any service work. For hard-wired machines, connection should be made via a suitable mains switch with all-pole isolation, which in the off position ensures a 3 mm gap between all open contacts.

The plug and socket as well as the mains switch must be easily accessible after the machine has been installed.

### Equipotential bonding and earthing

For added safety the machine should be protected with a residual current device with a trip current of 30 mA.

Equipotential bonding should be carried out if required. The screw connection point for equipotential bonding is located at the back of the machine. Equipotential bonding and earthing must be carried out before the machine is commissioned.

### Plumbing



The machine can be plumbed in to the cold water, hot water, demineralised water and drainage points without a backflow protection device.

Use the **Y-connector supplied with the machine to connect the water intake to the machine (KW) and the steam condenser (DK)**.



Alternatively you can provide an additional cold water supply to connect the steam condenser to.

### Communication modules


Ethernet and RS232 communication modules are available as optional extras. They are not supplied with the machine. The junction box for these should be installed near the machine for transferring and printing process data. These must be installed and connected in accordance with IEC 60950.


Dimensions

EL	Electrical supply cable	WW	Hot water connection
KW	Cold water connection	AW	Drain water
DK	Cold water connection, steam condenser	VE	Demineralised water connection
NW	Network and printer connection (optional)	PA	Equipotential bonding

### A - A

Connection points

EL	Electrical supply cable	WW	Hot water connection
KW	Cold water connection	AW	Drain water
DK	Cold water connection, steam condenser	VE	Demineralised water connection
NW	Network and printer connection (optional)	PA	Equipotential bonding

## Electrical connection (3 x 20 A)

Voltage <b>(default)</b>	3 AC 208V 60 Hz	3 AC 208V 60 Hz
Rated load	6.2 kW	6.2 kW
Fuse rating	3 x 20 A	3 x 20 A
Minimum cross-section	4 x AWG 12	4 x AWG 12
Length of electrical cord (UL 921) with plug	approx. 5'9"	approx. 1.8 m
Electrical connection	NEMA L15-20	NEMA L15-20

Voltage <b>(convertible)</b>	2 AC 208V 60 Hz	2 AC 208V 60 Hz
Rated load	6.2 kW	6.2 kW
Fuse rating	2 x 30 A	2 x 30 A
Minimum cross-section	3 x AWG 10	3 x AWG 10
Length of electrical cord (UL 921) with plug	approx. 5'9"	approx. 1.8 m
Electrical connection	NEMA L6-30	NEMA L6-30

## Cold water connections\* (KW+DK)

Max. temperature	68 °F	20 °C
Max. water hardness	73 gpg	70 °d
Recommended flow pressure	29 psi	200 kPa
Min. flow pressure	14.5 psi	100 kPa
Max. pressure	145 psi	1,000 kPa
Throughput	2 gal/min	7.5 l/min
Connection thread - on site (flat sealing)	3/4" Male Garden Hose Thread (USA: 11.5 NH)	3/4" Male Garden Hose Thread
Length of connection hose (supplied)	approx. 5' 7"	approx. 1.7 m
Length of the steam condensor hose	approx. 5' 7"	approx. 1.7 m

\* Two cold water connections required: use supplied Y-adapter if only one cold water connection is available.

## Hot water connection (WW)

Max. temperature	150 °F	65 °C
Max. water hardness	73 gpg	70 °dH
Recommended flow pressure	29 psi	200 kPa
Min. flow pressure	5.8 psi	40 kPa
Max. pressure	145 psi	1,000 kPa
Throughput	2 gal/min	7.5 l/min
Connection thread - on site (flat sealing)	3/4" Male Garden Hose Thread (USA: 11.5 NH)	3/4" Male Garden Hose Thread
Length of connection hose (supplied)	approx. 5' 7"	approx. 1.7 m

## Demineralised water (VE)

Max. temperature	150 °F	65 °C
Max. water hardness	29 psi	200 kPa
Recommended flow pressure	4.4 psi	30 kPa
Min. flow pressure	145 psi	1,000 kPa
Max. pressure	2 gal/min	7.5 l/min
Connection thread - on site (flat sealing)	11.5 NH	3/4" Male Garden Hose Thread
Length of connection hose (supplied)	approx. 5' 7"	approx. 1.7 m

**Machines with a booster pump for demineralised water (ADP) are only suitable for connection/use with a non-pressurised system.**

Minimum flow pressure (AD unpressurised)	1,3 psi	8.5 kPa
Maximum pressure (AD unpressurised)	8.7 psi	60 kPa
Machine connection nipple (outer dia. x l)	1/4" x 3/16"	6 x 30 mm

## Drain water (AW)

Drainage temperature	200 °F	93 °C
Length of Commercial Dishwasher drain hose	approx. 4' 11"	approx. 1.5 m
Drain hose, max. drainage length	approx. 13' 1"	approx. 4.0 m
Drain pump head height from floor level, max.	approx. 3' 3"	approx. 1 m
Drain hose delivery, max.	4.2 gal/min	16 l/min
On-site hose connector (external diameter x length)	7/8" x 1 3/16"	22 x 30 mm

## Machine feet

Height adjustable	0 - 5/16"	0-8 mm
Foot diameter	1 3/8"	35 mm
Entry point into housing	8 M	8 M

## Machine data

Unit height	32 1/4"	820 mm
Height including cover	32 7/8"	835 mm
Width	23 9/16"	598 mm
Depth	27 1/2"	698 mm
Door height	24 1/2"	622 mm
Net weight	216 lbs	98 kg
Dynamic floor load	2.000 N	2.000 N
Width including transport pallet	46 1/16"	1.170 mm
Depth including transport pallet	29 1/8"	740 mm
Height including transport pallet	39 3/8"	1.000 mm
Sound emissions in dB (A), Sound pressure LpA, washing and drying	< 70 dB	< 70 dB


## Heat dissipation to room

Heat transfer by emission	1194 btu/h	0.35 kWh/h
Load on removal	1365 btu/h	0.40 kWh/h

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**Miele Limited / Miele Limitée**