# A 620, A 621, A 622

<b>de</b>	Gebrauchsanweisung Injektormodul
<b>da</b>	Brugsanvisning Injektormodul
<b>en</b>	Operating instructions Injector module
<b>es</b>	Instrucciones de manejo Módulo inyector
<b>fi</b>	Käyttöohje Suorasuihkumoduuli
<b>fr</b>	Mode d'emploi Module à injection
<b>it</b>	Istruzioni d'uso Modulo a iniezione
<b>nl</b>	Gebruiksaanwijzing injectormodule
<b>no</b>	Bruksanvisning Injektormodul
<b>pt</b>	Instruções de utilização Módulo injetor
<b>sv</b>	Bruksanvisning Injektormodul

<b>de</b>	.....	4
<b>da</b>	.....	19
<b>en</b>	.....	34
<b>es</b>	.....	49
<b>fi</b>	.....	64
<b>fr</b>	.....	79
<b>it</b>	.....	94
<b>nl</b>	.....	109
<b>no</b>	.....	124
<b>pt</b>	.....	139
<b>sv</b>	.....	154

---

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---

<b>Notes about these instructions .....</b>	34
<b>Appropriate use .....</b>	35
Queries and technical problems.....	35
<b>Items supplied .....</b>	36
Load carriers .....	36
A 620 .....	36
A 620 + A 850.....	36
A 621 .....	36
A 621 + A 851.....	37
A 622 .....	37
A 622 + A 852.....	38
Disposal of the packing material .....	39
Optional accessories.....	39
<b>Warning and Safety instructions .....</b>	40
<b>Areas of application .....</b>	41
To check before loading the machine and before starting a programme .....	41
Sample loads.....	41
Shapes of items.....	41
A 620 .....	42
A 621 .....	42
A 622 .....	42
Preparing the load .....	43
<b>Installation.....</b>	46
Tools required:.....	46
Installation .....	46
Screwing the injector nozzle onto the module .....	46
Securing the frame .....	47
Removing the frame .....	47

## Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property.

Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

## Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

## Additional information and comments

Additional information and comments are contained in a simple frame.

## Operating steps

Operating steps are indicated by a black square bullet point.

### Example:

- Select an option using the arrow buttons and save your choice with OK.

## Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

### Example:

Settings  menu.

This module can be used to reprocess machine-reprocessable laboratory glassware and utensils in a Miele washer-disinfector for laboratory glassware and utensils. Follow the operating instructions for your washer-disinfector as well as the instructions of the glassware and utensil manufacturer on how to reprocess their items by machine.

The A 620, A 621 and A 622 injector modules are suitable for reprocessing narrow-necked laboratory glassware.

The modules are to be used in an A 503 mobile unit.

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessable laboratory glassware and utensils are referred to as “items” if they are not more closely defined.

### **Queries and technical problems**

If you have any queries or technical problems that you would like to discuss please contact Miele. Contact details can be found at the back of the Operating Instructions for your washer-disinfector or go to [www.miele-professional.com](http://www.miele-professional.com).

## en - Items supplied

### Load carriers

A 620


- A 620 injector module, height 113 mm, width 142 mm, depth 614 mm,  
sample load on injector nozzles A 840 and A 841

A 620 + A 850


- A 620 injector module with A 850 frame and centring aids, sample load on injector nozzles A 840 and A 841

A 621

- A 621 injector module, height 113 mm, width 184 mm, depth 614 mm,  
sample load on injector nozzles A 842 and A 843

A 621 + A 851


- A 621 injector module with A 851 frame and centring aids, sample load on injector nozzles A 842 and A 841

## en - Items supplied


---

**A 622**

- A 622 injector module, height 113 mm, width 208 mm, depth 614 mm,  
sample load on injector nozzles on A 844 and A 845

**A 622 + A 852**

- A 622 injector module with A 852 frame, sample load on injector nozzles A 844 and A 845

<b>Disposal of the packing material</b>	The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.  Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in landfill sites.
<b>Optional accessories</b>	The following and other accessories are available to order from Miele: <ul style="list-style-type: none"><li>- A 843, injector nozzle, length 185 mm, Ø 4 mm The injector nozzle is suitable for items with an opening width of 10 to 70 mm.</li><li>- A 840, injector nozzle, length 130 mm, Ø 6 mm The injector nozzle is suitable for items with an opening width of 12 to 85 mm.</li><li>- A 841, injector nozzle, length 210 mm, Ø 6 mm The injector nozzle is suitable for items with an opening width of 12 to 85 mm.</li><li>- A 842, injector nozzle, length 90 mm, Ø 4 mm The injector nozzle is suitable for items with an opening width of 10 to 70 mm.</li><li>- A 844, injector nozzle, length 80 mm, Ø 2.5 mm The injector nozzle is suitable for items with an opening width of 6 to 55 mm.</li><li>- A 845, injector nozzle, length 125 mm, Ø 2.5 mm The injector nozzle is suitable for items with an opening width of 6 to 55 mm.</li></ul>
A 850	- A 850 frame for A 620, with centring aids
A 851	- A 851 frame for A 621, with centring aids
A 852	- A 852 frame for A 622

## en - Warning and Safety instructions

Read the operating instructions carefully before using this load carrier. This will help protect users from personal injury, and help protect the load carrier from damage.

Keep these operating instructions in a safe place.

 Please also read the operating instructions for your washer-disinfector and pay particular attention to the Warning and Safety instructions.

- The module is approved solely for the applications specified in the operating instructions. Components such as nozzles may only be replaced with Miele accessories or genuine Miele original spare parts. Alterations or conversion of the machine, or using it for purposes other than those for which it was intended, are not permitted and could be dangerous.
- New load carriers must be cleaned in the washer-disinfector without a load prior to first use.
- Inspect all load carriers as described in the “Maintenance” section in the operating instructions for the cleaning machine.
- Only items which have been declared by their manufacturer as suitable for machine reprocessing may be processed. The manufacturer's specific reprocessing instructions must be observed.
- Broken glass can result in serious injury when loading or unloading. Damaged glass items must not be reprocessed in the machine.
- Only ever place empty, unloaded modules in the mobile unit. Check that they are correctly engaged before loading them. Modules must be completely emptied before removing them. Placing loaded modules in mobile units or removing them can cause damage to the wash load and cause injury should glassware be broken.
- The reprocessing result must be checked as necessary using a suitable test method and not just by a visual check.

Miele cannot be held liable for damage caused by non-compliance with these Warning and Safety Instructions.

The operating instructions for the loading carriage describe how to insert and remove the modules.

### To check before loading the machine and before starting a programme

- Check that fittings such as irrigation sleeves and injector nozzles are securely screwed into position.

**⚠** Make sure that all screw connectors are fitted with injector nozzles, adapters, irrigation sleeves or blind stoppers to ensure that all fittings in use are supplied with sufficient standardised pressure.

Damaged fittings such as injector nozzles, adapters and irrigation sleeves must not be used.

Fittings not equipped with wash items do not have to be replaced blind stoppers.

- Is the module correctly docked to the water supply in the injector unit?

### Sample loads

#### Shapes of items

Shapes of items	Laboratory flasks	Round flasks	Erlenmeyer flasks	Measuring flasks
				

## en - Areas of application

---

### A 620

The A 620 module has 10 positions for reprocessing items. It is designed for items with a volume from 200 ml to 1000 ml.

#### Capacity

Volume [ml]	Laboratory flasks	Round flasks	Erlenmeyer flasks	Measuring flasks
200–500	Max. 10	Max. 10	Max. 10	Max. 10
1000	Max. 10	Max. 5*	Max. 5*	Max. 5*

\* Plus 5 x items with a smaller volume

### A 621

The A 621 module has 20 positions for reprocessing items. It is designed for items with a volume from 50 ml to 250 ml.

#### Capacity

Volume [ml]	Laboratory flasks	Round flasks	Erlenmeyer flasks	Measuring flasks
50	Max. 20	Max. 20	Max. 20	–
100–200	Max. 20	Max. 20	Max. 20	Max. 20
250	Max. 20	Max. 10*	Max. 10*	Max. 10*

\* Plus 10 x items with a smaller volume

### A 622


The A 622 module has 36 positions for reprocessing items. It is designed for items with a volume from 20 ml to 100 ml.

#### Capacity

Volume [ml]	Laboratory flasks	Round flasks	Erlenmeyer flasks	Measuring flasks
20–50	Max. 36	Max. 36	Max. 36	Max. 36
100	Max. 36	Max. 18*	Max. 18*	Max. 18*


\* Plus 18 x items with a smaller volume

## Preparing the load

Ideally, a gap should be left between neighbouring items when they are reprocessed.

An additional frame can be mounted when reprocessing particularly fragile items. This prevents them from touching each other if they are moved by the force of water during reprocessing.

If bulbous items are loaded in neighbouring positions in the module, the items may touch each other.

If the cleaning and rinsing requirements are particularly high, a different position must be selected on the module for these items.


 Damage to the items.

If items touch each other during reprocessing, damage can occur at the points of contact, for example, the items may become scratched or the glass may break.

If the items are fragile:

- Select a different position on the module
- Use a frame
- Select a module with a greater distance between the nozzles

## en - Areas of application


### Damage to the items.

The rims of the items should rest fully on the support of the injector nozzle. If items dislodge neighbouring nozzles, damage can occur during reprocessing, for example, the items may become scratched or the glass may break.

Use:

- A different position on the module
- A module with a greater distance between the nozzles


The tip of the nozzle is designed so that it can touch the bottom of a piece of glassware during reprocessing. The shape of the protective cap ensures that water can flow out of the nozzle during reprocessing.


### Damage to the items.

If items are reprocessed frequently, their surfaces may be damaged at the point where the nozzle touches them, for example, the items may become scratched.

If the items are fragile, select a shorter nozzle to prevent the nozzle from coming into contact with the item.

If the cleaning and rinsing requirements are particularly high, a shorter nozzle must be selected.


Glassware must not extend beyond the sides of the load carrier.

## Tools required:

- 9 mm Spanner (WAF 9)


## Installation

### Screwing the injector nozzle onto the module


- Lift up the support on the injector nozzle and screw the injector nozzles into the required position ①.
  - Use the spanner to tighten the injector nozzle ②.
  - Slide the support on the injector nozzle downwards until it rests on the water supply pipe ③.
- Screw the other injector nozzles into the required positions.

### Securing the frame

- Position the rods for the frame on the module ①.
  - Turn the supports on the rods of the frame ②.
- The locking retainers on the supports lock into place ③.
- Repeat the process for the other 3 rods.

### Removing the frame

- Pull the locking retainers downwards slightly ① and turn the supports ②.
- Repeat the process for the other 3 rods.
- Lift the frame off the module ③.

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