# O 801/2 upper injector basket

## Items supplied:

- O 801/2 upper injector basket, H 287 mm, W 381 mm, D 475 mm,
- 10 x Ø 4 mm injector jets, 30 mm long


## Intended use:

The O 801/2 upper basket is designed for machine reprocessing dental, gynaecological and ENT instruments.

O 801/2 upper basket is for use in the G 7831 washer-disinfector.

The middle axle of the upper basket has an injector bar with 10 holders for processing hollow instruments. There is space at the front and back for inserts and mesh trays for processing instruments without lumens. Loading examples and further tips can be found in "Areas of application" in the Operating instructions for the washer-disinfector.

Use only fixtures specifically designed for the preparation of transmission instruments, such as the AUF 2 holders.

⚠ Please observe the Warning and Safety instructions in the Operating instructions for your washer-disinfector carefully. You should only process instruments/ medical products in this basket if they have been declared as suitable for machine processing by the manufacturer. Please also observe the manufacturer's care instructions. The programme selected must be suitable for the materials from which the items to be processed are made. Where possible the programme should include thermal disinfection. The final rinse should, if possible, always be carried out with fully de-mineralised water.

## Application:

For validation purposes please make sure that you follow the loading instructions given on the template.


⚠ Instruments which have become disconnected from the jets/fixtures during processing must be re-processed.

Before connecting hollow instruments, please make sure that all hollow sections are free of blockages. Otherwise the cleaning solution will not be able to access all areas properly.

- Place hollow instruments onto the jets. The height-adjustable frame stabilises the instruments.
- At the end of the programme, instruments must be dried with sterile compressed air before further processing or sterilisation in accordance with the manufacturer's instructions. Please observe local and national health and safety regulations.

## To adjust the height of the retaining frame:

- Lift the retaining frame upwards
- fix it again on both sides on level ① or ②.

## Before loading the machine and before starting a programme, please check that:

- the silicone holders are all present and undamaged
- injector jets and fixtures are secure in their silicone holders.

Missing or water-damaged silicone holders must be replaced with new ones. Only then may the insert be used again. If these conditions are present silicone holders standardise the pressure themselves.

Damaged jets, adapters or sleeves etc. must not be used.

- that the spray arms (basket and machine upper spray arm) can rotate freely.

## Optional accessories:

If several accessories are being used to process hollow instruments at the same time, check that the number is permitted with regard to validation. The maximum number specified for validation must not be exceeded.

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For processing transmission instruments up to 6 additional AUF 2 adapters are required. The holders come with their own instructions.

Before using hand pieces and contra angles again following reprocessing, a function check must be carried out, e.g. by spraying into a basin, to ensure they are clear.

### E 499 fixture:

The E 499 fixture is for processing rinsing and suction cannulae, including suction breakers. The fixture comes with its own instructions.

### E 442 and E 443 sleeves for instruments with lumens:

Injector sleeves E 442 and E 443 are for processing hollow instruments and instruments less than Ø 10 mm.

- E 442 for instruments of 4 - 8 mm Ø,
- E 443 for instruments of 8 - 8.5 mm Ø,

The apertures can be varied using optional replacement caps:

- E 445 cap für sleeves; aperture Ø 6 mm,
- E 446 cap for sleeves; aperture Ø 10 mm,

The diameters of the instruments and the sleeve apertures or caps must be compatible with each other, e.g. caps for sleeves with an aperture diameter of Ø 6 mm for dental surgical suction cannulae.

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