O 177/1 upper injector basket

Items supplied:

- O 177/1 upper injector basket, H 281 mm, W 531 mm, D 475 mm.
- 26 x Ø 4 mm injector jets, 30 mm long,
- 7 funnels.

Intended use:

The O 177/1 upper injector basket is designed for the machine reprocessing of medical instruments, e.g. dental instruments.

The O 177/1 upper injector basket can be used in G 7881, G 7882, G 7882 CD, G 7891, G 7892, PG 8535 and PG 8536 washer-disinfectors.


On the left hand side there are two injector bars with a total of 26 holders for reprocessing hollow instruments. The right hand side is free for inserts and mesh trays for reprocessing instruments without lumens.

Loading examples and further tips can be found in "Areas of application" in the Operating instructions for the washer-disinfector.

Only use fixtures specifically designed for the reprocessing of transmission instruments, such as the AUF 1 holders.

Spray arm monitoring

A magnetic spray arm has to be retro-fitted to PG 8535 and PG 8536 if spray arm monitoring is required.

⚠ Please observe the Warning and Safety instructions in the Operating instructions for your washer-disinfector carefully. You should only process instruments / medical products in this basket if they have been declared as suitable for machine processing by the manufacturer. Please also observe the manufacturer's care instructions. The programme selected must be suitable for the materials from which the items to be processed are made. Where possible the programme should include thermal disinfection. The final rinse should, if possible, always be carried out with fully de-mineralised water.

Application:

For validation purposes please make sure that you follow the loading instructions given on the template.


⚠ Instruments which have become disconnected from the jets/holders during processing must be re-processed.

Before connecting hollow instruments, please make sure that all hollow sections are free of blockages. Otherwise the cleaning solution will not be able to access all areas properly.

- Place hollow instruments onto the jets. The height-adjustable frame stabilises the instruments.
- At the end of the programme, instruments must be dried with sterile compressed air before further processing or sterilisation in accordance with the manufacturer's instructions. Please observe local and national health and safety regulations.

To adjust the height of the retaining frame:

- Lift the retaining frame upwards,
- then fit it on both sides on level ① or ②.

Before loading the machine and before starting a programme, please check:

- that the silicone holders are all present and undamaged.
- that the injector jets and fixtures are secure in their silicone holders.

Missing or water-damaged silicone holders must be replaced with new ones. Only then may the basket be used again. If these conditions are present silicone holders standardise the pressure themselves.

Damaged jets, adapters or sleeves etc. must not be used.

- that the spray arms (basket and machine upper spray arm) rotate freely.

Optional accessories:

If several accessories are being used to reprocess hollow instruments at the same time, check that the number is permitted with regard to validation. The maximum number specified for validation must not be exceeded.

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Up to 11 additional AUF 1 adapters are required for reprocessing transmission instruments. The holders are supplied with their own Operating instructions.

Before using hand pieces and contra angles again following reprocessing, a function check must be carried out, e.g. by spraying into a basin, to ensure they are clear.

E 499 fixture:

The E 499 fixture is for reprocessing rinsing and suction cannulae, including suction breakers. The fixture comes with its own instructions.

E 442 and E 443 sleeves for hollow instruments:

Injector sleeves E 442 and E 443 are for reprocessing hollow instruments and instruments with an external diameter less than 10 mm.

- E 442 for 4 - 8 mm Ø instruments,
- E 443 for 8 - 8.5 mm Ø instruments.

The aperture can be varied using optional replacement caps:

- E 445 cap for sleeves; aperture Ø 6 mm,
- E 446 cap for sleeves; aperture Ø 10 mm.

The diameter of the instruments and the sleeve apertures or caps must be compatible with each other, e.g. caps for sleeves with an aperture diameter of Ø 6 mm for dental surgical suction cannulae.

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