Mobile unit E 327/1 for surgical instruments

Miele
PROFESSIONAL

⚠ Please read the Warning and Safety instructions in the Operating Instructions for your washer-disinfector carefully. You should only process instruments/medical products in this mobile unit if they have been declared as suitable for machine processing by the manufacturer. Please also observe the manufacturer's care instructions. The programme selected must be suitable for the materials from which the medical products to be processed are made. Where possible the programme should include thermal disinfection. The final rinse should, if possible, always be carried out with fully de-mineralised water.
If using automatic mobile unit recognition please make sure that the coding is correctly set.

Items supplied:

Mobile unit E 327/1 with 2 levels, H 502, W 535, D 515 mm, with a holder for the magnetic strip used in automatic mobile unit recognition (G 7835, G 7836, PG 8535 and PG 8536) and spray arm magnets used for spray arm monitoring (PG 8535 and PG 8536 only).

Accessories available to order:

- Magnet strip ML/2 for automatic mobile unit recognition
- E 142 DIN mesh tray
- E 476 holders for micro instruments
- E 473 insert for very small items

Areas of application:

The E 327/1 mobile unit is suitable for use in the following washer-disinfectors: G 7882, G 7882 CD, G 7892, G 7835, PG 8535, G 7836 and PG 8536.

It enables surgical instruments to be cleaned and disinfected in 4 DIN mesh trays. 2 DIN mesh trays can be accommodated on each level.



Loading height (levels numbered from the bottom upwards) [H x W x D in mm]

Level 1	112 x 520 x 510
Level 2	105 x 512 x 480

Before loading the machine and before starting a programme please check:

- that the spray arms rotate freely,
- that the spray arm jets are clean and free of any blockages,
- that the mobile unit is correctly fitted into the cabinet and that it has connected correctly with the water intake.

For validation purposes please make sure that you follow the loading instructions given on the template.


■ Hinged instruments must be opened up before being placed in the mesh trays to minimise the risk of surfaces not being accessed during processing.

■ To prevent the risk of light instruments such as bulldog clips and other small items being moved about during processing, silicone holders e.g. E 473 and closable inserts e.g. E 476 can be used to keep them secure.

■ Containers e.g. beakers and kidney dishes should be placed in the mobile unit, preferably on the top level, with the open side facing downwards, ensuring that all surfaces can be accessed easily and that cleaning solution can flow freely off them. They must also be secure so that they cannot be flipped over by the force of water during processing.

If hollow items are placed such that cleaning solution can collect in recesses etc. they will not get properly clean and rinsing results will be impaired.

Access port for wash pressure measurements, e.g. within context of performance testing and validation in accordance with DIN ISO 15883.

■ To carry out wash pressure measurements, exchange the screw stopper for a Luer lock adapter, e.g. E 447.

⚠ Under no circumstances may the access port be used to connect instruments, spray rails, etc.

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