XKM 3100 WCDS

de	Gebrauchsanweisung Kommunikationsmodul
da	Brugsanvisning Kommunikationsmodul
en	Operating instructions Communication module
es	Instrucciones de manejo Módulo de comunicación
fr	Mode d'emploi Module de communication
it	Istruzioni d'uso Modulo di comunicazione
nl	Gebruiksaanwijzing communicatiemodule
pl	Instrukcja użytkowania Moduł komunikacyjny
pt	Instruções de utilização Módulo de comunicação
sv	Bruksanvisning Kommunikationsmodul

de	4
da	18
en	32
es	46
fr	60
it	74
nl	88
pl	102
pt	116
sv	130

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Notes about these instructions	32
Intended use	33
Queries and technical problems	33
Items supplied	34
Items supplied	34
Technical data	34
Guide to the appliance	35
Machine overview	35
Switching on and off	35
Control LEDs	36
Button	36
WiFi connection (WPS)	36
Restore factory settings (Reset)	36
Firmware update	36
Warning and Safety instructions	37
Fitting and installation	38
Range and availability of WiFi signal	38
Module slot	38
Installation	39
Configuring the interface	40
Activating the interface on the washer-disinfector	40
Changing the report language	40
Downloading the app	41
Establishing the WiFi connection	41
Operation	42
Retrospective output of reports	42
Firmware Update	42
Firmware update via the Miele app	42
Firmware update via the module	43
Caring for the environment	44
Disposal of the packing material	44

Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

Additional information and comments

Additional information and comments are contained in a simple frame.

Operating steps

Operating steps are indicated by a black square bullet point.

Example:

- Select an option using the arrow buttons and save your choice with *OK*.

Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

Example:

Settings  menu.

The communication module supplied XKM 3100 WCDS enables a WiFi connection to be established between one of the Miele washer-disinfectors listed below and a tablet. The Miele app DataDiary must be installed on the tablet for this.

The Miele app DataDiary provides the process documentation. The app is available to download for free from the respective app stores. Further information on process documentation can be obtained from Miele.

Only use tablets that are approved according to EN/IEC 60950. The interface must only be configured by suitably qualified and competent persons.

The communication module can be used in the following washer-disinfectors:


- PG 8562
- PG 8581
- PG 8582
- PG 8582 CD
- PG 8583
- PG 8583 CD
- PG 8591
- PG 8592
- PG 8593

The communication module is generally referred to as the module in these Operating instructions.

Queries and technical problems

If you have any queries or technical problems that you would like to discuss please contact Miele. Contact details can be found at the back of the Operating Instructions for your washer-disinfector or go to www.miele-professional.com.

Items supplied


- Communication module XKM 3100 WCDS

Technical data

Voltage, connected load, fuse rating	See data plate
Baud rate	9600 (Setting on washer-disinfectors)
Bit	8
Parity	None (Setting on washer-disinfectors)
Stop bits	1
WiFi standard	802.11 b/g/n
Temperature during operation	Up to max. 80 °C
Test certification	UL Recognised Component cULus , ⓓ
CE -Labelling	2004/108/EG EMC Directive, 2006/95/EG Low Voltage Directive , 2011/65/EU RoHS Directive 2014/53/EU Radio Installation Directive
Manufacturer's address	Miele & Cie. KG, Carl-Miele-Straße 29, 33332 Gütersloh, Germany

Machine overview

- ① Reset button
- ② Indicator LED
indicates operational readiness after installation and active communication when in operation
- ③ Data plate (located on the side)

Switching on and off

The module does not have a separate On/Off switch. It is switched on and off together with the washer-disinfector.

The module requires approx. 30 s for initialisation after switch-on. Only then can files be transmitted or settings be changed on the interface.

Control LEDs

The control LEDs (light emitting diodes) give status information about the communication module.

LED	Status
LIT UP GREEN	The module is ready for use and a WiFi connection has been established.
LIT UP YELLOW	The module is ready for use and a WiFi connection can be established.
FLASHES GREEN	A WiFi connection is being established via WPS procedure.
FLASHES YELLOW	Resetting in progress.
FLASHES alternately GREEN and RED	Firmware update in progress.
FLASHES RED	Communication fault has occurred. Check the configuration of the RS232 interface on the washer-disinfector.

Button

The button on the module has the following functions:

WiFi connection (WPS)

WPS log in is required to connect the module directly to the internet, e.g. for updating the module Firmware.
Instructions for logging in to WPS can be found in “Firmware update”.

Restore factory settings (Reset)

The module can be reset to factory default settings by continuously pressing the button.

- Press and hold the reset button for approx. 10-15 seconds until the LED begins to flash yellow.

Restoring the settings takes approx. 30 seconds. The module will then carry out an automatic restart and will restart in AP mode. The restart is completed when the LED on the module lights up steady green or yellow.

Firmware update

An update of the module Firmware can be initialised from the Miele app or by pressing the button.
Instructions for this can be found in “Firmware update”.

To avoid the risk of accidents and damage to this communication module please read these instructions carefully before using it for the first time.
Keep these instructions in a safe place.

 Please also read the operating instructions for your washer-disinfector and pay particular attention to the Warning and Safety instructions.

- ▶ This communication module must only be used for applications specified in the operating instructions. Changes or modifications to it and incorrect use of it are not permitted as this can be dangerous.
- ▶ The communication module must only be installed, commissioned and serviced by Miele service or by a suitably qualified person. Repairs are only to be undertaken by Miele Service. Repair work by unqualified persons could be dangerous.
- ▶ Check the communication module carefully before installation for any visible damage. A damaged module could be dangerous and must not be installed.
- ▶ Before installing the communication module compare the connection data (voltage and frequency) on the data plate with that on the data plate of your washer-disinfector. If the data does not match, the module could sustain damage. If in any doubt, consult an electrician.
- ▶ The communication module housing must not be opened on any account.
Touching live electrical components or altering the electrical or mechanical components is dangerous and could cause the module to malfunction.
- ▶ The WiFi technology supported by the module can be temporarily or completely disrupted e.g. by structural conditions or the set-up location. For this reason availability of featured functions cannot always be guaranteed.

Range and availability of WiFi signal


The WiFi range or signal strength is dependent on the structural environment. Reinforced concrete walls and ceilings, metal surfaces, metal-vaporised films, glass screens and fire protection constructions can muffle or divert the signal, therefore reducing the range.

For example, if the washer-disinfector is located under a metal worktop, or is surrounded by metal cabinets, the signal will be almost completely blocked.

Module slot

The module is inserted into the module slot. The module slot is located on the back of the washer-disinfector.

- If the washer-disinfector is installed in a niche:
 - Pull the washer-disinfector forwards out of the niche.
- If the washer-disinfector is installed under a worktop:
 - Open the door.

- Unscrew the fixing screws.
- Pull the washer-disinfector out from under the worktop.

Installation

⚠ Before fitting the communication module compare the connection data (voltage and frequency) on the data plate with that on the data plate of your washer-disinfector. If the data does not match, the module could sustain damage. If in any doubt, consult an electrician.

- Take a note of the module serial number. This is required as the password for establishing the WiFi connection. The serial number can be found on the data plate located on the side of the module.


The serial number of the module, without the initial zeroes, acts as the password. Example:

12/0001234567890

Password: **1234567890**

- Disconnect the washer-disinfector from the mains electricity supply.
- Remove any stickers or protective caps as necessary from the slot aperture.

- Push the module into the module slot. Gently push the module until it clicks into place.
- Connect the washer-disinfector to the electrical supply.
The module is supplied with mains power by the washer-disinfector. An additional power connection is not required.

Activating the interface on the washer-disinfector

The module XKM 3100 WCDS uses the configuration settings of the RS232 interface of the washer-disinfector.

- Open the menu for configuring the RS232 interface via Further settings/Interface/RS232.

- In the Mode menu, select the Terminal option and confirm your choice with *OK*.
- Change the Baud rate to 9600 (factory default setting).

The RS232 interface is now activated for the WiFi module. After approx. 30 seconds the module is ready for operation.

Changing the report language

The RS232 interface only supports a limited number of languages. Therefore the output language for cycle reports is set separately from the language setting of the washer-disinfector and the Miele app DataDiary.

- Open the menu for configuring the RS232 interface via Further settings/Interface/RS232.
- Select a language in the Language  menu and confirm your choice with *OK*.

Downloading the app

Download the Miele app DataDiary for free from the relevant App Store for process documentation.

The operating system of your tablet needs to be at least iOS 10[®] or Android 8[™].

- Search in App Store[®] or Google Play[™] for the Miele app **DataDiary**.

- Download the app to your tablet.

Establishing the WiFi connection

To connect via the app you will need:

- A tablet with the Miele app DataDiary installed on it.
- The serial number of your Miele module XKM 3100 WCDS. This can be found on the data plate.
- Switch the washer-disinfector on. After approx. 30 seconds the module is ready for operation.
- Connect the tablet with the module WiFi.
 - The WiFi SSID is DataDiary-AP-XXXXXX, where XXXXXX corresponds to the last 6 letters of the module MAC address.
 - The WiFi password is the serial number of the module without the initial zeroes.

Retrospective output of reports

If there are WiFi problems during data transmission, the cycle report can be transmitted retrospectively.

- Open the menu as follows:
Further settings/Interface/RS232/Print reports

Further instructions and additional information can be found in “Process documentation” in the Operating instructions for your washer-disinfector.

Firmware Update

An update of the module Firmware may be necessary at irregular intervals. Two different procedures are available:

1. Update via the Miele app DataDiary.
2. Update directly via the module XKM 3100 WCDS in your washer-disinfector. The module must be connected to the internet via a WPS-enabled router for this.

 Under no circumstances must the internet be accessed via the washer-disinfector, e.g. via port forwarding!
Make sure that your internet connection or your router is configured as appropriate.

The current module Firmware version can be selected via the app.

Firmware update via the Miele app

The Miele app DataDiary has an Update wizard which you can use to update the module Firmware.

For the update you will need:

- Administrator rights for the Miele app.
- The WiFi access code to establish a connection with the internet.
- Switch the washer-disinfector on. After approx. 30 seconds the module is ready for operation.
- Log into the app as administrator and follow the wizard instructions.

Firmware update via the module

For the update you will need:

- A WPS-enabled router that is connected to the internet.
- Switch the washer-disinfector on. After approx. 30 seconds the module is ready for operation.
- Press the button on the module for approx. 3-8 seconds until the LED begins to flash green.

 Do not press it for any longer than this as this would reset the module to its factory default settings.

When the button is released the LED lights up steady yellow.

- Activate the WPS function on your WiFi router within 120 seconds, e.g. by pressing the WPS button (if available).

If there is a connection, the LED lights up steady green.

- Then press the button again for approx. 18-23 seconds until the LED begins to flash green.

During the update the LED flashes alternately green and red. After the update the module carries out a restart. The module starts in AP mode. The restart is completed when the LED on the module lights up steady green or yellow.

- Green: The update has completed successfully and the connection with the tablet has been established.
- Yellow: The update has completed successfully but the connection with the tablet has not been established. Check whether the tablet is switched on and the WiFi connection is activated.


Disposal of the packing material

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in landfill sites.

Disposing of your old appliance

Electrical and electronic appliances contain many valuable materials. They also contain certain materials, compounds and components which were essential for their correct functioning and safety. These could be hazardous to human health and to the environment if disposed of with household waste or if handled incorrectly. Please do not, therefore, dispose of your old appliance with household waste.

Instead, please make use of officially designated collection and disposal points to dispose of and recycle electrical and electronic appliances in your local community, with your dealer or with Miele, free of charge. By law, you are solely responsible for deleting any personal data from the old appliance prior to disposal. You are legally obliged to remove any old batteries which are not securely enclosed by the appliance and to remove any lamps without destroying them, where this is possible. These must be taken to a suitable collection point where they can be handed in free of charge. Please ensure that your old appliance poses no risk to children while being stored for disposal.

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