

Printing date 01.10.2015 Revision: 01.10.2015

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- · 1.1 Product identifier
- · Trade name: Miele ProCare Universal 61
- · CAS Number:

7647-14-5

· EC number:

231-598-3

 1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

- · Application of the substance / the mixture Auxiliary
- · 1.3 Details of the supplier of the safety data sheet
- · Manufacturer/Supplier:

Miele & Cie. KG

Carl-Miele-Straße 29

33332 Gütersloh

Tel.: +49 (0)5241/89-0

Fax.: +49 (0)5241/89-2090

www.miele-professional.com

claro products GmbH

Vogelsangstraße 1

A-5310 Mondsee

Tel.: +43 (0)6232/2626

Fax.: +43 (0)6232/2626-575

office@claro.at www.claro.at

· Further information obtainable from:

Labor: Tel.: +43 (0) 6232 / 2626 - 568

Contact: Mr. Dr. Lukasser; josef.lukasser@claro.at

· 1.4 Emergency telephone number:

Poison control centre Vienna: Tel.: +43 (0) 1/4064-343

SECTION 2: Hazards identification

- · 2.1 Classification of the substance or mixture
- · Classification according to Regulation (EC) No 1272/2008

 The substance is not classified according to the CLP regulation.
- · 2.2 Label elements
- · Labelling according to Regulation (EC) No 1272/2008 Void
- · Hazard pictograms Void
- · Signal word Void
- · Hazard statements Void
- · 2.3 Other hazards
- · Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.

SECTION 3: Composition/information on ingredients

- · 3.1 Chemical characterisation: Substances
- · CAS No. Description

7647-14-5 Sodium Chloride

- Identification number(s)
- **EC** number: 231-598-3

EU -



Printing date 01.10.2015 Revision: 01.10.2015

Trade name: Miele ProCare Universal 61

(Contd. of page 1)

SECTION 4: First aid measures

- · 4.1 Description of first aid measures
- · General information: No special measures required.
- · After inhalation:

Supply fresh air; consult doctor in case of complaints.

- · After skin contact: Generally the product does not irritate the skin.
- · After eye contact:

Rinse opened eye for several minutes under running water.

- · After swallowing: If symptoms persist consult doctor.
- · 4.2 Most important symptoms and effects, both acute and delayed No further relevant information available.
- · 4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

SECTION 5: Firefighting measures

- · 5.1 Extinguishing media
- · Suitable extinguishing agents:

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

- · 5.2 Special hazards arising from the substance or mixture No further relevant information available.
- · 5.3 Advice for firefighters
- · Protective equipment: No special measures required.

SECTION 6: Accidental release measures

- · 6.1 Personal precautions, protective equipment and emergency procedures
 Not required.
- · 6.2 Environmental precautions:

Do not allow to enter sewers/ surface or ground water.

- 6.3 Methods and material for containment and cleaning up: Pick up mechanically.
- · 6.4 Reference to other sections

No dangerous substances are released.

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

SECTION 7: Handling and storage

- · 7.1 Precautions for safe handling No special measures required.
- Information about fire and explosion protection:
 No special measures required.
- · 7.2 Conditions for safe storage, including any incompatibilities
- · Storage:
- · Requirements to be met by storerooms and receptacles:
 No special requirements.
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: None.
- · 7.3 Specific end use(s) No further relevant information available.

EU



Printing date 01.10.2015 Revision: 01.10.2015

Trade name: Miele ProCare Universal 61

(Contd. of page 2)

SECTION 8: Exposure controls/personal protection

- · Additional information about design of technical facilities:
 No further data; see item 7.
- · 8.1 Control parameters
- · Ingredients with limit values that require monitoring at the workplace:
 Not required.
- · Additional information:

The lists valid during the making were used as basis.

- · 8.2 Exposure controls
- · Personal protective equipment:
- · General protective and hygienic measures:

The usual precautionary measures are to be adhered to when handling chemicals.

- · Respiratory protection: Not required.
- · Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation ${\sf consideration}$

· Material of gloves

Butyl rubber, BR Nitrile rubber, NBR Natural rubber, NR

· Penetration time of glove material

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection: Not required.

SECTION 9: Physical and chemical properties

· 9.1 Information on basic physical · General Information	and chemical properties
· Appearance:	
Form:	Solid
Colour:	White
· Odour:	Odourless
· Odour threshold:	Not determined.
· pH-value (10 g/l) at 20 °C:	4,5 - 7
· Change in condition Melting point/Melting range: Boiling point/Boiling range:	
· Flash point:	Not applicable.
· Flammability (solid, gaseous):	Product is not flammable.
· Ignition temperature:	
Decomposition temperature:	Not determined.
· Self-igniting:	Not determined.
· Danger of explosion:	Product does not present an explosion hazard.
· Explosion limits:	
Lower:	Not determined.
Upper:	Not determined.

(Contd. on page 4)



Printing date 01.10.2015 Revision: 01.10.2015

Trade name: Miele ProCare Universal 61

	(Contd. of page 3)
· Vapour pressure:	Not applicable.
· Density at 20 °C:	2,16 g/cm³
 Bulk density at 20 °C: Relative density Vapour density Evaporation rate 	~ 1140 kg/m³ Not determined. Not applicable. Not applicable.
· Solubility in / Miscibility with water at 20 °C:	358 g/l
 Partition coefficient (n-octanol, water): 	/ Not determined.
· Viscosity: Dynamic: Kinematic:	Not applicable. Not applicable.
· Solvent content: Organic solvents: VOC (EC)	0,0 % 0,00 %
Solids content: · 9.2 Other information	100,0 % No further relevant information available.

SECTION 10: Stability and reactivity

- · 10.1 Reactivity No further relevant information available.
- · 10.2 Chemical stability
- · Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

- · 10.3 Possibility of hazardous reactions
 - No dangerous reactions known if used according to specifications.
- · 10.4 Conditions to avoid No further relevant information available.
- · 10.5 Incompatible materials: No further relevant information available.
- · 10.6 Hazardous decomposition products:
 - No dangerous decomposition products known.

SECTION 11: Toxicological information

- · 11.1 Information on toxicological effects
- · Acute toxicity

Based on available data, the classification criteria are not met.

- · LD/LC50 values relevant for classification:
 - 7647-14-5 Sodium chloride

Oral LD50 3000 mg/kg (rat)

- · Primary irritant effect:
- · Skin corrosion/irritation

Based on available data, the classification criteria are not met.

- · Serious eye damage/irritation
 - Based on available data, the classification criteria are not met.
- · Respiratory or skin sensitisation

Based on available data, the classification criteria are not met.

- · CMR effects (carcinogenity, mutagenicity and toxicity for reproduction)
- · Germ cell mutagenicity

Based on available data, the classification criteria are not met.

(Contd. on page 5)



Printing date 01.10.2015 Revision: 01.10.2015

Trade name: Miele ProCare Universal 61

(Contd. of page 4)

· Carcinogenicity

Based on available data, the classification criteria are not met.

· Reproductive toxicity

Based on available data, the classification criteria are not met.

· STOT-single exposure

Based on available data, the classification criteria are not met.

· STOT-repeated exposure

Based on available data, the classification criteria are not met.

· Aspiration hazard

Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

- · 12.1 Toxicity
- · Aquatic toxicity: No further relevant information available.
- · 12.2 Persistence and degradability

No further relevant information available.

· 12.3 Bioaccumulative potential

No further relevant information available.

- · 12.4 Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes:

Water hazard class 1 (German Regulation) (Assessment by list): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · 12.5 Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.
- · 12.6 Other adverse effects No further relevant information available.

SECTION 13: Disposal considerations

- · 13.1 Waste treatment methods
- · Recommendation

Smaller quantities can be disposed of with household waste.

- · Uncleaned packaging:
- · Recommendation:

Disposal must be made according to official regulations.

SECTION I	L4:	Transport	information

•	14.1	<i>UN-Number</i>

· ADR, ADN, IMDG, IATA Void

· 14.2 UN proper shipping name

· ADR, ADN, IMDG, IATA Void

· 14.3 Transport hazard class(es)

· ADR, ADN, IMDG, IATA

· Class Void

· 14.4 Packing group

· ADR, IMDG, IATA Void

· 14.5 Environmental hazards:

· Marine pollutant: No

(Contd. on page 6)



Printing date 01.10.2015 Revision: 01.10.2015

Trade name: Miele ProCare Universal 61

(Contd. of page 5)

· 14.6 Special precautions for user Not applicable.

· 14.7 Transport in bulk according to Annex II of Marpol and the IBC

Code Not applicable.

· UN "Model Regulation":

SECTION 15: Regulatory information

· 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Void

- · National regulations:
- · Additional classification according to Decree on Hazardous Materials, Annex II:
- · Waterhazard class:

Water hazard class 1 (Assessment by list): slightly hazardous for

· 15.2 Chemical safety assessment:

A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- · Department issuing MSDS: Labor: Tel.: +43 (0)6232/2626-568
- · Contact: Mr. Josef Lukasser
- Abbreviations and acronyms:

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by

ICAO: International Civil Aviation Organisation

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

 ${\tt vPvB:} \ \, {\tt very} \ \, {\tt Persistent} \ \, {\tt and} \ \, {\tt very} \ \, {\tt Bioaccumulative}$

· * Data compared to the previous version altered.

EU