

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Miele ProCare Lab 30 P

1.2. Relevant identified uses of the substance or mixture and uses advised against

Identified Uses

PC35

Washing and cleaning products (including solvent based products)

1.3. Details of the supplier of the safety data sheet

Address/Manufacturer

Miele & Cie. KG

Carl-Miele-Straße 29

33241 Gütersloh

Telephone no.

+49 5241 89 0

Fax no.

+49 5241 89 2090

E-mail address of
person responsible

sida@drweigert.de

for this SDS

Address/Supplier

Miele Inc.

9 Independence Way

Princeton, NJ 08540

Telephone no.

+1-800-991-9380

1.4. Emergency telephone number

GBK/ Infotrac: (USA domestic) +1 800 535 5053 or international +1 352 323 3500

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Classification according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Skin Corr. 1B

H314

Eye Dam. 1

H318

Met. Corr. 1

H290

2.2. Label elements

Labelling according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Hazard pictograms



Signal word

Danger

Safety Data Sheet according to Appendix D, OSHA Hazard Communication Standard 29 CFR 1910:1200

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

Hazard statements

H314 Causes severe skin burns and eye damage.
H290 May be corrosive to metals.

Precautionary statements**Prevention**

P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

Storage/Disposal

Dispose only when container is empty and closed. For disposal of product residues, refer to safety data sheet.

Hazardous component(s) to be indicated on label

contains Phosphoric acid

2.3. Other hazards**Supplemental information**

No special hazards have to be mentioned. The product contains no PBT or vPvB substances.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Hazardous ingredients according to OSHA Hazard Communication Standard 29 CFR 1910:1200****phosphoric acid**

CAS No. 7664-38-2
Concentration \geq 50 %

citric acid

CAS No. 77-92-9
Concentration \geq 1 < 10 %

Other information

Complete text of hazard statements in chapter 16

SECTION 4: First aid measures**4.1. Description of first aid measures****General information**

Remove contaminated, soaked clothing immediately and dispose of safely. Clean body thoroughly (bath, shower). In any case show the physician the Safety Data Sheet.

After inhalation

Ensure supply of fresh air. When spray fog inhaled, seek medical aid.

After skin contact

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

After contact with skin, wash immediately with plenty of water. Take medical treatment.

After eye contact

In case of contact with the eyes, rinse immediately for at least 15 minutes with plenty of water. Summon a doctor immediately.

After ingestion

If swallowed, seek medical advice immediately and show this container or label. Rinse mouth thoroughly with water. Let plenty of water be drunk in small gulps. Do not induce vomiting.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed**Hints for the physician / hazards**

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Product itself is non-combustible; adapt fire extinguishing measures to surrounding areas.

Non suitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixture

In case of combustion evolution of dangerous gases possible.

5.3. Advice for firefighters**Special protective equipment for fire-fighting**

Do not inhale explosion and/or combustion gases. In case of combustion use a suitable breathing apparatus.

Other information

Collect contaminated fire-fighting water separately, must not be discharged into the drains. Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures**

Avoid contact with skin, eyes and clothing. Refer to protective measures listed in Sections 7 and 8.

6.2. Environmental precautions

Do not discharge into the drains/surface waters/groundwater.

6.3. Methods and material for containment and cleaning up

Pick up with absorbent material. Dispose of absorbed material in accordance with the regulations.

6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Avoid formation of aerosols. Observe the usual precautions for handling chemicals. Keep container tightly closed.

Advice on protection against fire and explosion

The product is not combustible.

7.2. Conditions for safe storage, including any incompatibilities

Recommended storage temperature

Value > -15 < 30 °C

Requirements for storage rooms and vessels

Keep in original packaging, tightly closed. Storage rooms must be properly ventilated. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Storage class according to TRGS 510

Storage class according to TRGS 510 8B Non-combustible corrosive hazardous substances

7.3. Specific end use(s)

no data

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Other information

There are not known any further control parameters.

8.2. Exposure controls

General protective and hygiene measures

Hold eye wash fountain available. Hold emergency shower available. Do not inhale gases/vapours/aerosols. Avoid contact with skin and eyes. Do not eat, drink or smoke during work time. Wash hands before breaks and after work. Clean skin thoroughly after work; apply skin cream.

Respiratory protection

If workplace limits are exceeded, a respiratory protection approved for this particular job must be worn. Particle filter P2

Hand protection

Chemical resistant gloves	
Use	Permanent hand contact
Appropriate Material	neoprene
Material thickness	>= 0,65
Breakthrough time	> 480
Appropriate Material	nitrile
Material thickness	>= 0,4
Breakthrough time	> 480
Appropriate Material	butyl
Material thickness	>= 0,7
Breakthrough time	> 480

Safety Data Sheet according to Appendix D, OSHA Hazard Communication Standard 29 CFR 1910:1200

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

Use Short-term hand contact
Appropriate Material nitrile
Material thickness \geq 0,11
Hand protection must comply with EN 374.

Eye protection

Safety glasses with side protection shield; Eye protection must comply with EN 166.

Body protection

Clothing as usual in the chemical industry. Protective shoes

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Form liquid
Colour colourless
Odour characteristic
Odour threshold
Remarks not determined
pH value
Value appr. 0,7
Temperature 20 °C

Melting point
Remarks not determined

Freezing point
Remarks not determined

Initial boiling point and boiling range
Remarks not determined

Flash point
Remarks Not applicable

Evaporation rate (ether = 1) :
Remarks not determined

Flammability (solid, gas)
Not applicable

Upper/lower flammability or explosive limits
Remarks Not applicable

Vapour pressure
Remarks not determined

Vapour density
Remarks not determined

Density
Value 1,47 g/cm³
Temperature 20 °C

Solubility in water
Remarks miscible in all proportions

Solubility(ies)

Safety Data Sheet according to Appendix D, OSHA Hazard Communication Standard 29 CFR 1910:1200

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

Remarks not determined

Partition coefficient: n-octanol/water

Remarks not determined

Ignition temperature

Remarks Not applicable

Decomposition temperature

Remarks not determined

Viscosity

dynamic

Value < 50 mPa.s
Temperature 20 °C

Explosive properties

evaluation not determined

Oxidising properties

evaluation None known

9.2. Other information

Other information

None known

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

No hazardous reactions known.

Decomposition temperature

Remarks not determined

10.5. Incompatible materials

Reactions with metals, with evolution of hydrogen. Reactions with alkalies.

10.6. Hazardous decomposition products

No hazardous decomposition products known.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute oral toxicity

ATE > 2000 mg/kg
Method calculated value (Regulation (EC) No. 1272/2008)
Remarks Based on available data, the classification criteria are not met.

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

Acute oral toxicity (Components)

citric acid

Species	rat		
LD50		11700	mg/kg

citric acid

Species	mouse		
LD50		5040	mg/kg

Phosphoric acid

Species	rat		
LD50		2600	mg/kg

Acute dermal toxicity

Remarks Based on available data, the classification criteria are not met.

Acute dermal toxicity (Components)

Phosphoric acid

Species	rabbit		
LD50		2740	mg/kg

Acute inhalational toxicity

Remarks Based on available data, the classification criteria are not met.

Skin corrosion/irritation

evaluation	corrosive
Remarks	The classification criteria are met.

Serious eye damage/irritation

evaluation	corrosive
Remarks	The classification criteria are met.

Sensitization

Remarks Based on available data, the classification criteria are not met.

Subacute, subchronic, chronic toxicity

Remarks Based on available data, the classification criteria are not met.

Mutagenicity

Remarks Based on available data, the classification criteria are not met.

Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

National Toxicology Program (NTP)

No component of this product present at levels greater than or equal to 0,1% is identified as a known or reasonably anticipated to be human carcinogen by NTP

International Agency for research on Cancer (IARC)

No component of this product present at levels greater than or equal to 0,1% is identified as probable, possible or confirmed human carcinogen by IARC

Specific Target Organ Toxicity (STOT)

Single exposure

Remarks Based on available data, the classification criteria are not met.

Repeated exposure

Remarks Based on available data, the classification criteria are not met.

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

Aspiration hazard

Based on available data, the classification criteria are not met.

Experience in practice

Inhalation may lead to irritation of the respiratory tract.

Other information

There is no data available on the product apart from the information given in this subsection.

SECTION 12: Ecological information

12.1. Toxicity

General information

not determined

Fish toxicity (Components)

citric acid

Species	golden orfe (<i>Leuciscus idus</i>)	
LC50	440 to 706	mg/l
Duration of exposure	96 h	

Phosphoric acid

Species	mosquito fish	
LC50	138	mg/l
Duration of exposure	96 h	

Daphnia toxicity (Components)

citric acid

Species	Daphnia magna	
EC50	120	mg/l
Duration of exposure	72 h	

Phosphoric acid

Species	Daphnia magna	
EC50	> 100	mg/l
Duration of exposure	48 h	
Method	OECD 202	

Algae toxicity (Components)

Phosphoric acid

Species	Scenedesmus subspicatus	
EC50	> 100	mg/l
Duration of exposure	72 h	
Method	OECD 201	

12.2. Persistence and degradability

General information

not determined

Ready degradability (Components)

citric acid

12.3. Bioaccumulative potential

General information

not determined

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

Partition coefficient: n-octanol/water

Remarks

not determined

12.4. Mobility in soil**General information**

not determined

12.5. Results of PBT and vPvB assessment**General information**

not determined

Evaluation of persistence and bioaccumulation potential

The product contains no PBT or vPvB substances.

12.6. Other adverse effects**General information**

not determined

General information / ecology

Do not allow to enter soil, waterways or waste water canal. Avoid release into the atmosphere.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations for the product**

The listed waste code numbers, according to the European Waste Catalogue (EWC), are to be understood as a recommendation. A final decision must be made in agreement with the regional waste disposal company.

Disposal recommendations for packaging

Completely emptied packagings can be given for recycling.

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

SECTION 14: Transport information




Safety Data Sheet according to Appendix D, OSHA Hazard Communication Standard 29 CFR 1910:1200

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

	Ground transport DOT	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
Tunnel restriction code	E		
IMDG-Code segregation group		1 Acids	
14.1. UN number	1805	1805	1805
14.2. UN proper shipping name	PHOSPHORIC ACID, SOLUTION	PHOSPHORIC ACID, SOLUTION	PHOSPHORIC ACID, SOLUTION
14.3. Transport hazard class(es)	8	8	8
Label			
14.4. Packing group	III	III	III
Limited Quantity	5 l		
Transport category	3		
14.5. Environmental hazards		no	

Information for all modes of transport

14.6. Special precautions for user

See Sections 6 to 8

Other information

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Other regulations, restrictions and prohibition regulations

Observe employment restrictions for young people.

Other information

The product does not contain substances of very high concern (SVHC).

US. EPA Emergency Planning and Community Right-To-Know Act (EPCRA) SARA Title III Section 302Extremely Hazardous Substance (40 CFR 355)

The product does not contain any listed components.

US. EPA Emergency Planning and Community Right-To-Know Act (EPCRA) SARA Title III Section 313Toxic Chemicals (40 CFR 372.65) - Supplier Notification Required

Safety Data Sheet according to Appendix D, OSHA Hazard Communication Standard 29 CFR 1910:1200

Trade name: Miele ProCare Lab 30 P

Version: 2 / US

Replaces Version: - / US Date revised: 16.09.2022

Print date: 28.09.22

The product does not contain any listed components.

California Safe Drinking Water and Toxic Enforcement Act (Proposition 65)

Warning! This product may contain trace quantities of substance(s) known to the state of California to cause cancer and/or reproductive toxicity – not added as a part of the formulation but remaining as residuals from the manufacturing process of our raw material suppliers.

SECTION 16: Other information

Abbreviations

ADR: Accord européen relatif au transport international des marchandises Dangereuses par Route

RID: Règlement concernant le transport international ferroviaire de marchandises dangereuses

IMDG: International Maritime Code for Dangerous Goods

ICAO: International Civil Aviation Organization

IATA: International Air Transport Association

VOC: Volatile Organic Compound

LD: Lethal dose

LC: Lethal concentration

PBT: Persistent, Bioaccumulative and Toxic

vPvB: Very persistent and very bioaccumulative

SVHC: Substances of very high concern

MARPOL 73/78: International Convention for the Prevention of Pollution From Ships, 1973 as modified by the Protocol of 1978 (MARPOL: Marine Pollution)

IBC: Intermediate Bulk Container

CAS: Chemical Abstracts Service

ISO: International Organization for Standardization

OEL: Occupational exposure limit

OECD: Organisation for Economic Co-operation and Development

UN: United Nations

IMO: International Maritime Organization

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: ***

This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.