

**Safety data sheet
according to Regulation (EC) No 1907/2006, Article 31**

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

**SECTION 1: Identification of the substance/mixture and of the company/
undertaking**

- **1.1 Product identifier**
- **Trade name: Miele ProCare Universal 77A**
- **UFI:** H880-Y08G-X00T-A1A6
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the mixture** Cleaning material/ Detergent
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Miele & Cie. KG
Carl-Miele-Straße 29
33332 Gütersloh
Tel.: +49 (0)5241/89-0
Fax.: +49 (0)5241/89-2090
www.miele-professional.com
- claro products GmbH
Sonystraße 20
A-5081 Anif
Tel.: +43 (0)6232/2626
Fax.: +43 (0)6232/2626-575
office@claro.at
www.claro.at
- **Further information obtainable from:**
Labor: Tel.: +43 (0) 6232 / 2626 - 568
Contact: Mr. Dr. Lukasser; josef.lukasser@claro.at
- **1.4 Emergency telephone number:**
National Poisons Information Centre: +353 (1) 809 2166 (8.00 a.m. to 10.00 p.m. 7 days a week)
Healthcare Professionals: +353 (1) 809 2566 (24 hour service)
Emergency CONTACT (24-Hour-Number): GBK GmbH +49 (0)6132-84463

SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**



Eye Irrit. 2 H319 Causes serious eye irritation.

- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008**
The product is classified and labelled according to the CLP regulation.
- **Hazard pictograms**



- **Signal word** Warning

(Contd. on page 2)

Safety data sheet according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 1)

· **Hazard statements**

H319 Causes serious eye irritation.

· **Precautionary statements**

P101 If medical advice is needed, have product container or label at hand.

P102 Keep out of reach of children.

P280 Wear eye protection / face protection.

P301+P312 IF SWALLOWED: Call a POISON CENTER doctor if you feel unwell.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

· **2.3 Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

· **3.2 Mixtures**

· **Description:** Mixture of substances listed below with nonhazardous additions.

· **Dangerous components:**

CAS: 497-19-8 EINECS: 207-838-8	sodium carbonate ☠ Eye Irrit. 2, H319	>30%
CAS: 15630-89-4 EINECS: 239-707-6	SODIUM CARBONATE PEROXIDE ☠ Ox. Sol. 2, H272; ☠ Eye Dam. 1, H318; ☠ Acute Tox. 4, H302 Specific concentration limits: Eye Dam. 1; H318: C ≥ 25 % Eye Irrit. 2; H319: 10 % ≤ C < 25 %	≥5-<10%
CAS: 9038-95-3	Oxirane, methyl-, polymer with oxirane, monobutyl ether ☠ Acute Tox. 4, H302	0.2-5%
CAS: 1344-09-8 EINECS: 215-687-4	Silicic acid, sodium salt ☠ Acute Tox. 4, H332; Skin Irrit. 2, H315; Eye Irrit. 2, H319; STOT SE 3, H335	0.2-5%
CAS: 29329-71-3 EINECS: 249-559-4	(1-Hydroxyethyliden) - bisphosphonsäure, Natriumsalz ☠ Acute Tox. 4, H302; Eye Irrit. 2, H319	0.2-5%
CAS: 10543-57-4 EINECS: 234-123-8	N,N'-Ethylenebis[N-acetylacetamide] ☠ Acute Tox. 4, H332	0.2-5%

· **Regulation (EC) No 648/2004 on detergents / Labelling for contents**

phosphates	≥15 - <30%
oxygen-based bleaching agents	≥5 - <15%
non-ionic surfactants, polycarboxylates, phosphonates, perfumes	<5%

· **Additional information:** For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

· **4.1 Description of first aid measures**

· **After inhalation:** Supply fresh air; consult doctor in case of complaints.

· **After skin contact:**

Wash off with plenty of water. Remove clothing contaminated with product. In case of complaints, consult a doctor.

(Contd. on page 3)

Safety data sheet
according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 2)

- **After eye contact:**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** If symptoms persist consult doctor.
- **4.2 Most important symptoms and effects, both acute and delayed**
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, powder or water spray. Fight larger fire with alcohol resistant foam.
Use fire extinguishing methods suitable to surrounding conditions.
- **5.2 Special hazards arising from the substance or mixture**
No further relevant information available.
- **5.3 Advice for firefighters**
- **Protective equipment:** The usual measures for fire protection are to be taken.
- **Additional information**
Dispose of fire debris and contaminated fire fighting water in accordance with official regulations.

SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures**
Wear protective equipment. Keep unprotected persons away.
- **6.2 Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **6.3 Methods and material for containment and cleaning up:**
Dispose contaminated material as waste according to section 13.
Pick up mechanically (wipe up, sweep up) and collect in suitable containers for disposal.
- **6.4 Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling**
Ensure good ventilation of the workplace. Wear personal protective equipment. Avoid contact with eyes and skin.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep container tightly sealed.
- **Storage class:** 13
- **7.3 Specific end use(s)** No further relevant information available.

(Contd. on page 4)

Safety data sheet
according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 3)

SECTION 8: Exposure controls/personal protection

· **8.1 Control parameters**

· **Ingredients with limit values that require monitoring at the workplace:**

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

· **DNELs**

497-19-8 sodium carbonate

Inhalative	DNEL Langzeit - lokale Wirkung	10 mg/m ³ (worker)
	DNEL Akut - lokale Wirkung	10 mg/m ³ (Consumer)

15630-89-4 SODIUM CARBONATE PEROXIDE

Dermal	DNEL Akut - lokale Wirkung	12.8 mg/cm ² (worker)
		6.4 mg/cm ² (Consumer)
Inhalative	DNEL Langzeit - lokale Wirkung	5 mg/m ³ (worker)

1344-09-8 Silicic acid, sodium salt

Oral	DNEL Langzeit - lokale Wirkung	1.59 mg / kg KG /Tag (worker)
Inhalative	DNEL Langzeit - lokale Wirkung	5.61 mg/m ³ (worker)

10543-57-4 N,N'-Ethylenebis[N-acetylacetamide]

Oral	DNEL Langzeit - systemische Wirkung	0.45 mg/kg KG/Tag (human)
Dermal	DNEL Langzeit - systemische Wirkung	10 mg / kg KG / Tag (Consumer)
		20 mg / kg KG / Tag (worker)
Inhalative	DNEL Langzeit - lokale Wirkung	75 mg/m ³ (Consumer)
		6.4 mg/m ³ (worker)

· **PNECs**

15630-89-4 SODIUM CARBONATE PEROXIDE

PNEC STP	16.24 mg/l (watertreatment plant)
PNEC aquatic	35 µg/l (Algae)

1344-09-8 Silicic acid, sodium salt

PNEC STP	348 mg/l (watertreatment plant)
PNEC aquatic	7,500 µg/l (algae)

10543-57-4 N,N'-Ethylenebis[N-acetylacetamide]

PNEC STP	10 mg/l (watertreatment plant)
PNEC aquatic	10 µg/l (algae)

· **Additional information:** The lists valid during the making were used as basis.

· **8.2 Exposure controls**

· **Appropriate engineering controls** Ensure good ventilation of the workplace.

· **Individual protection measures, such as personal protective equipment**

· **General protective and hygienic measures:**

- Keep away from foodstuffs, beverages and feed.
- Immediately remove all soiled and contaminated clothing
- Wash hands before breaks and at the end of work.
- Avoid contact with the eyes.

· **Respiratory protection:**

- Not necessary if room is well-ventilated.
- Wear a dust mask if dust is generated: dust mask with filter type P2 (EN143) or fine dust mask FFP2 / FFP3 (EN149).

(Contd. on page 5)

Safety data sheet according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 4)

· **Hand protection**



Protective gloves according to EN 374

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

Recommended material thickness: $\geq 0.4\text{mm}$

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Butyl rubber, BR
Nitrile rubber, NBR
Natural rubber, NR

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

For the mixture of chemicals mentioned below the penetration time has to be at least 480 minutes (Permeation according to EN 16523-1:2015: Level 6).

· **Eye/face protection**



Safety goggles according to EN166

· **Body protection:**

The above advice on protective equipment refers to commercial handling of larger quantities.

SECTION 9: Physical and chemical properties

· **9.1 Information on basic physical and chemical properties**

· **General Information**

· **Physical state**

Solid

· **Colour:**

According to product specification

· **Odour:**

Characteristic

· **Odour threshold:**

Not determined.

· **Melting point/freezing point:**

Undetermined.

· **Boiling point or initial boiling point and boiling range**

Undetermined.

· **Flammability**

Not determined.

· **Lower and upper explosion limit**

· **Lower:**

Not determined.

· **Upper:**

Not determined.

· **Flash point:**

Not applicable.

· **Decomposition temperature:**

Not determined.

· **pH at 20 °C**

10.9 (1%)

(Contd. on page 6)

Safety data sheet according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 5)

<ul style="list-style-type: none"> · Viscosity: · Kinematic viscosity · Dynamic: · Solubility · water: · Partition coefficient n-octanol/water (log value) · Vapour pressure: · Density and/or relative density · Density: · Relative density · Vapour density · Particle characteristics 	<p>Not applicable.</p> <p>Not applicable.</p> <p>Soluble.</p> <p>Not determined.</p> <p>Not applicable.</p> <p>Not determined.</p> <p>Not determined.</p> <p>Not applicable.</p> <p>See section 3.</p>
<ul style="list-style-type: none"> · 9.2 Other information · Appearance: · Form: · Important information on protection of health and environment, and on safety. · Ignition temperature: · Explosive properties: · Solvent content: · Organic solvents: · Solids content: · Change in condition · Evaporation rate 	<p>Solid</p> <p>Product is not selfigniting.</p> <p>Product does not present an explosion hazard.</p> <p>0.0 %</p> <p>100.0 %</p> <p>Not applicable.</p>
<ul style="list-style-type: none"> · Information with regard to physical hazard classes · Explosives · Flammable gases · Aerosols · Oxidising gases · Gases under pressure · Flammable liquids · Flammable solids · Self-reactive substances and mixtures · Pyrophoric liquids · Pyrophoric solids · Self-heating substances and mixtures · Substances and mixtures, which emit flammable gases in contact with water · Oxidising liquids · Oxidising solids · Organic peroxides · Corrosive to metals · Desensitised explosives 	<p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p> <p>Void</p>

SECTION 10: Stability and reactivity

- **10.1 Reactivity** The product is not reactive under normal conditions of use, storage and transport.

(Contd. on page 7)

Safety data sheet
according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 6)

- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions**
No dangerous reactions known if used according to specifications.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.

- **LD/LC50 values relevant for classification:**

497-19-8 sodium carbonate

Oral	LD50	4,090 mg/kg (rat)
------	------	-------------------

15630-89-4 SODIUM CARBONATE PEROXIDE

Oral	LD50	1,034 mg/kg (rat)
Dermal	LD50	>2,000 mg/kg (rabbit)

9038-95-3 Oxirane, methyl-, polymer with oxirane, monobutyl ether

Oral	LD50	>300-2,000 mg/kg (rat)
------	------	------------------------

1344-09-8 Silicic acid, sodium salt

Oral	LD50	3,400 mg/kg (rat)
Dermal	LD50	>5,000 mg/kg (rat)
Inhalative	LC50/4 h	2.06 mg/l (rat)

10543-57-4 N,N'-Ethylenebis[N-acetylacetamide]

Oral	LD50	>2,000 mg/kg (rat)
Dermal	LD50	>2,000 mg/kg (rat)
Inhalative	LC50/4 h	>2.08 mg/l (rat)

- **Primary irritant effect:**
- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Causes serious eye irritation.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **11.2 Information on other hazards**

- **Endocrine disrupting properties**

None of the ingredients is listed.

(Contd. on page 8)

Safety data sheet
according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 7)

SECTION 12: Ecological information

· **12.1 Toxicity**

· **Aquatic toxicity:**

497-19-8 sodium carbonate

LC50/96 h (static)	300 mg/l (Lepomis macrochirus)
EC50/48 h	200-227 mg/l (Ceriodaphnia dubia)

15630-89-4 SODIUM CARBONATE PEROXIDE

LC50/96 h	70.7 mg/l (Pimephales promelas)
EC50/48 h	4.9 mg/l (Daphnia pulex)
NOEC/48h	2 mg/l (Daphnia pulex)
NOEC/96h	7.4 mg/l (Pimephales promelas)

9038-95-3 Oxirane, methyl-, polymer with oxirane, monobutyl ether

LC50/96 h	>100 mg/l (Brachydanio rerio)
EC50/48 h	>100 mg/l (Daphnia magna)
EC50/72h	>100 mg/l (Scenedesmus subspicatus)

1344-09-8 Silicic acid, sodium salt

LC50/96 h	310 mg/l (Oncorhynchus mykiss)
ErC50/72h	>345.4 mg/l (algae)
EC50/48 h	1,700 mg/l (Daphnia magna)

· **12.2 Persistence and degradability**

The surfactants contained in this preparation fulfil the conditions of biodegradability as laid down in Regulation (EC) No 648/2004 on detergents. Documentation confirming this will be held at the disposal of the competent authorities of the Member States and will only be made available to them either at their direct request or at the request of a detergent manufacturer.

· **12.3 Bioaccumulative potential** No further relevant information available.

· **12.4 Mobility in soil** No further relevant information available.

· **12.5 Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

· **12.6 Endocrine disrupting properties**

The product does not contain substances with endocrine disrupting properties.

· **12.7 Other adverse effects**

· **Additional ecological information:**

· **General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
 Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

SECTION 13: Disposal considerations

· **13.1 Waste treatment methods**

· **Recommendation**

Must not be disposed together with household garbage. Do not allow product to reach sewage system.

Take the contents and container to a recognised hazardous waste disposal company or to a hazardous waste collection point.

(Contd. on page 9)

Safety data sheet
according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 8)

Small amounts may be diluted with plenty of water and washed away. Dispose of bigger amounts in accordance with Local Authority requirements.

· **European waste catalogue**

07 06 01*	aqueous washing liquids and mother liquors
-----------	--

· **Uncleaned packaging:**

· **Recommendation:**

Disposal must be made according to official regulations.

Empty contaminated packagings thoroughly. They may be recycled after thorough and proper cleaning.

Recycle or dispose of in accordance with the applicable legal requirements.

· **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

SECTION 14: Transport information

· **14.1 UN number or ID number**

· **ADR, ADN, IMDG, IATA** Void

· **14.2 UN proper shipping name**

· **ADR, ADN, IMDG, IATA** Void

· **14.3 Transport hazard class(es)**

· **ADR, ADN, IMDG, IATA**

· **Class** Void

· **14.4 Packing group**

· **ADR, IMDG, IATA** Void

· **14.5 Environmental hazards:**

· **Marine pollutant:** No

· **14.6 Special precautions for user** Not applicable.

· **14.7 Maritime transport in bulk according to IMO instruments** Not applicable.

· **UN "Model Regulation":** Void

SECTION 15: Regulatory information

· **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** None of the ingredients is listed.

· **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II**

None of the ingredients is listed.

· **REGULATION (EU) 2019/1148**

· **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3))**

None of the ingredients is listed.

(Contd. on page 10)

Safety data sheet
according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 9)

· **Annex II - REPORTABLE EXPLOSIVES PRECURSORS**

None of the ingredients is listed.

· **Regulation (EC) No 273/2004 on drug precursors**

None of the ingredients is listed.

· **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors**

None of the ingredients is listed.

· **National regulations:**

· **Additional classification according to Decree on Hazardous Materials, Annex II: -**

· **Waterhazard class:** Water hazard class 1 (Self-assessment): slightly hazardous for water.

· **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

This Safety Data Sheet is in compliance with Regulation (EC) No 1907/2006, Article 31 as amended by Regulation (EU) 2020/878.

· **Relevant phrases**

H272 May intensify fire; oxidiser.

H302 Harmful if swallowed.

H315 Causes skin irritation.

H318 Causes serious eye damage.

H319 Causes serious eye irritation.

H332 Harmful if inhaled.

H335 May cause respiratory irritation.

· **Classification according to Regulation (EC) No 1272/2008**

Serious eye damage/irritation

The classification of the mixture is generally based on the calculation method using substance data according to Regulation (EC) No 1272/2008.

· **Department issuing SDS:** Labor: Tel.: +43 (0)6232/2626-568

· **Contact:** Mr. Josef Lukasser

· **Date of previous version:** 18.01.2023

· **Version number of previous version:** 2.1

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

DNEL: Derived No-Effect Level (REACH)

PNEC: Predicted No-Effect Concentration (REACH)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Ox. Sol. 2: Oxidizing solids – Category 2

Acute Tox. 4: Acute toxicity – Category 4

Skin Irrit. 2: Skin corrosion/irritation – Category 2

(Contd. on page 11)

Safety data sheet
according to Regulation (EC) No 1907/2006, Article 31

Printing date 18.12.2025

Version number 2.2 (replaces version 2.1)

Revision: 18.12.2025

Trade name: Miele ProCare Universal 77A

(Contd. of page 10)

Eye Dam. 1: Serious eye damage/eye irritation – Category 1
Eye Irrit. 2: Serious eye damage/eye irritation – Category 2
STOT SE 3: Specific target organ toxicity (single exposure) – Category 3
· * **Data compared to the previous version altered.**

IE