

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Miele ProCare Lab 30 C

1.2. Relevant identified uses of the substance or mixture and uses advised against **Identified Uses**

PC35 Washing and cleaning products (including solvent based products)

1.3. Details of the supplier of the safety data sheet

Address/Manufacturer

Miele & Cie. KG
Carl-Miele-Straße 29
33241 Gütersloh
Telephone no. +49 5241 89 0
Fax no. +49 5241 89 2090
E-mail address of
person responsible sida@drweigert.de
for this SDS

Address/Supplier

Miele Inc.
9 Independence Way
Princeton, NJ 08540
Telephone no. +1-800-991-9380

1.4. Emergency telephone number

GBK/ Infotrac: (USA domestic) +1 800 535 5053 or international +1 352 323 3500

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Classification according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Eye Irrit. 2	H319
Skin Sens. 1A	H317

2.2. Label elements

Labelling according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Hazard pictograms



Signal word

Warning

Hazard statements

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

H319 Causes serious eye irritation.
H317 May cause an allergic skin reaction.

Precautionary statements**Prevention**

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P264.1 Wash hands thoroughly after handling.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P337 If eye irritation persists:
P362+P364 Take off contaminated clothing and wash it before reuse.

Storage/Disposal

P501.1 Dispose of contents/container to industrial incineration plant.

2.3. Other hazards**Supplemental information**

No special hazards have to be mentioned. The product contains no PBT or vPvB substances.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Hazardous ingredients according to OSHA Hazard Communication Standard 29 CFR 1910:1200****citric acid**

CAS No. 77-92-9
Concentration \geq 25 < 50 %

2-octyl-2H-isothiazol-3-one

CAS No. 26530-20-1
Concentration \geq 0,0025 < 0,025 %

Other information

Complete text of hazard statements in chapter 16

SECTION 4: First aid measures**4.1. Description of first aid measures****General information**

Remove contaminated, soaked clothing immediately and dispose of safely.

After inhalation

Ensure supply of fresh air. In the event of symptoms take medical treatment.

After skin contact

After contact with skin, wash immediately with plenty of water. Consult a doctor if skin irritation persists.

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

After eye contact

In case of contact with the eyes, rinse immediately for at least 15 minutes with plenty of water. In case of irritation consult an oculist.

After ingestion

Rinse mouth thoroughly with water.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed**Hints for the physician / hazards**

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Product itself is non-combustible; adapt fire extinguishing measures to surrounding areas.

Non suitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixture

In case of combustion evolution of dangerous gases possible.

5.3. Advice for firefighters**Special protective equipment for fire-fighting**

Do not inhale explosion and/or combustion gases. In case of combustion use a suitable breathing apparatus.

Other information

Collect contaminated fire-fighting water separately, must not be discharged into the drains. Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures**

Avoid contact with skin, eyes and clothing. Refer to protective measures listed in Sections 7 and 8.

6.2. Environmental precautions

Do not discharge into the drains/surface waters/groundwater.

6.3. Methods and material for containment and cleaning up

Pick up with absorbent material. Dispose of absorbed material in accordance with the regulations.

6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

SECTION 7: Handling and storage

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

7.1. Precautions for safe handling

Advice on safe handling

Avoid formation of aerosols. Observe the usual precautions for handling chemicals. Keep container tightly closed.

Advice on protection against fire and explosion

The product is not combustible.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep in original packaging, tightly closed. Storage rooms must be properly ventilated. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Storage class according to TRGS 510

Storage class according to TRGS 510 10-13 Other combustible and non-combustible substances

7.3. Specific end use(s)

no data

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Other information

There are not known any further control parameters.

8.2. Exposure controls

General protective and hygiene measures

Hold eye wash fountain available. Do not inhale gases/vapours/aerosols. Avoid contact with skin and eyes. Do not eat, drink or smoke during work time. Wash hands before breaks and after work.

Respiratory protection

Not necessary, but do not inhale vapours. If workplace limits are exceeded, a respiratory protection approved for this particular job must be worn.

Hand protection

Chemical resistant gloves

Use Permanent hand contact
Appropriate Material neoprene

Material thickness >= 0,65

Breakthrough time > 480

Appropriate Material nitrile

Material thickness >= 0,4

Breakthrough time > 480

Appropriate Material butyl

Material thickness >= 0,7

Breakthrough time > 480

Use Short-term hand contact

Appropriate Material nitrile

Material thickness >= 0,11

Hand protection must comply with EN 374.

Eye protection

Safety glasses with side protection shield; Eye protection must comply with EN 166.

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

Body protection

Clothing as usual in the chemical industry.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Form	liquid
Colour	colourless
Odour	odourless
Odour threshold	
Remarks	not determined
pH value	
Value	appr. 1,0
Temperature	20 °C
Melting point	
Remarks	not determined
Freezing point	
Remarks	not determined
Initial boiling point and boiling range	
Value	> 100 °C
Flash point	
Remarks	Not applicable
Evaporation rate (ether = 1) :	
Remarks	not determined
Flammability (solid, gas)	
Not applicable	
Upper/lower flammability or explosive limits	
Remarks	Not applicable
Vapour pressure	
Remarks	not determined
Vapour density	
Remarks	not determined
Density	
Value	1,17 g/cm ³
Temperature	20 °C
Solubility in water	
Remarks	miscible in all proportions
Solubility(ies)	
Remarks	not determined
Partition coefficient: n-octanol/water	
Remarks	not determined
Ignition temperature	
Remarks	Not applicable

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

Decomposition temperature

Remarks not determined

Viscosity**dynamic**Value < 10 mPa.s
Temperature 20 °C**Explosive properties**

evaluation no

Oxidising properties

evaluation None known

9.2. Other information**Other information**

None known

SECTION 10: Stability and reactivity**10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

No hazardous reactions known.

Decomposition temperature

Remarks not determined

10.5. Incompatible materials

Reactions with alkalis.

10.6. Hazardous decomposition products

No hazardous decomposition products known.

SECTION 11: Toxicological information**11.1. Information on toxicological effects****Acute oral toxicity**

Remarks Based on available data, the classification criteria are not met.

Acute oral toxicity (Components)**citric acid**Species rat
LD50 11700 mg/kg**citric acid**Species mouse
LD50 5040 mg/kg

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

Acute dermal toxicity

Remarks Based on available data, the classification criteria are not met.

Acute inhalational toxicity

Remarks Based on available data, the classification criteria are not met.

Skin corrosion/irritation

Remarks Based on available data, the classification criteria are not met.

Serious eye damage/irritation

evaluation irritant

Remarks The classification criteria are met.

Sensitization

evaluation May cause sensitization by skin contact.

Remarks The classification criteria are met.

Subacute, subchronic, chronic toxicity

Remarks Based on available data, the classification criteria are not met.

Mutagenicity

Remarks Based on available data, the classification criteria are not met.

Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

National Toxicology Program (NTP)

No component of this product present at levels greater than or equal to 0,1% is identified as a known or reasonably anticipated to be human carcinogen by NTP

International Agency for research on Cancer (IARC)

No component of this product present at levels greater than or equal to 0,1% is identified as probable, possible or confirmed human carcinogen by IARC

Specific Target Organ Toxicity (STOT)**Single exposure**

Remarks Based on available data, the classification criteria are not met.

Repeated exposure

Remarks Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Experience in practice

Inhalation may lead to irritation of the respiratory tract.

Other information

There is no data available on the product apart from the information given in this subsection.

SECTION 12: Ecological information**12.1. Toxicity****General information**

not determined

Fish toxicity (Components)

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

citric acid

Species	golden orfe (Leuciscus idus)			
LC50	440	to	706	mg/l
Duration of exposure	96	h		

Daphnia toxicity (Components)**citric acid**

Species	Daphnia magna			
EC50	120			mg/l
Duration of exposure	72	h		

12.2. Persistence and degradability**General information**

not determined

Ready degradability (Components)

citric acid

12.3. Bioaccumulative potential**General information**

not determined

Partition coefficient: n-octanol/water

Remarks not determined

12.4. Mobility in soil**General information**

not determined

12.5. Results of PBT and vPvB assessment**Evaluation of persistence and bioaccumulation potential**

The product contains no PBT or vPvB substances.

12.6. Other adverse effects**General information**

not determined

General information / ecology

Do not allow to enter soil, waterways or waste water canal. Avoid release into the atmosphere.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations for the product**

Allocation of a waste code number, according to the European Waste Catalogue (EWC), should be carried out in agreement with the regional waste disposal company.

Disposal recommendations for packaging

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

SECTION 14: Transport information

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
14.1. UN number	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.

Information for all modes of transport

14.6. Special precautions for user

See Sections 6 to 8

Other information

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

US. EPA Emergency Planning and Community Right-To-Know Act (EPCRA) SARA Title III Section 302 Extremely Hazardous Substance (40 CFR 355)

The product does not contain any listed components.

US. EPA Emergency Planning and Community Right-To-Know Act (EPCRA) SARA Title III Section 313 Toxic Chemicals (40 CFR 372.65) - Supplier Notification Required

Components: Triethylene glycol

California Safe Drinking Water and Toxic Enforcement Act (Proposition 65)

Warning! This product may contain trace quantities of substance(s) known to the state of California to cause cancer and/or reproductive toxicity – not added as a part of the formulation but remaining as residuals from the manufacturing process of our raw material suppliers.

SECTION 16: Other information

Abbreviations

ADR: Accord européen relatif au transport international des marchandises Dangereuses par Route

RID: Règlement concernant le transport international ferroviaire de marchandises dangereuses

IMDG: International Maritime Code for Dangerous Goods

ICAO: International Civil Aviation Organization

IATA: International Air Transport Association

IBC: Intermediate Bulk Container

CAS: Chemical Abstracts Service

VOC: Volatile Organic Compound

LD: Lethal dose

LC: Lethal concentration

PBT: Persistent, Bioaccumulative and Toxic

vPvB: Very persistent and very bioaccumulative

SVHC: Substances of very high concern

MARPOL 73/78: International Convention for the Prevention of Pollution From Ships, 1973 as modified by the Protocol of 1978 (MARPOL: Marine Pollution)

ISO: International Organization for Standardization

OECD: Organisation for Economic Co-operation and Development

Safety Data Sheet according to Appendix D, OSHA Hazard Communication Standard 29 CFR 1910:1200

Trade name: Miele ProCare Lab 30 C

Version: 4 / US

Replaces Version: 3 / US Date revised: 02.03.2022

Print date: 02.03.22

IMO: International Maritime Organization

UN: United Nations

EU: European Union

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: ***

This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.