

Trade name: Miele ProCare Dent 11 A

Version: 3 / US

Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Miele ProCare Dent 11 A

1.2. Relevant identified uses of the substance or mixture and uses advised against **Identified Uses**

PC35 Washing and cleaning products (including solvent based products)

1.3. Details of the supplier of the safety data sheet

Address/Manufacturer

Miele & Cie. KG
Carl-Miele-Straße 29
33241 Gütersloh
Telephone no. +49 5241 89 0
Fax no. +49 5241 89 2090
E-mail address of
person responsible sida@drweigert.de
for this SDS

Address/Supplier

Miele Inc.
9 Independence Way
Princeton, NJ 08540
Telephone no. +1-800-991-9380

1.4. Emergency telephone number

GBK/ Infotrac: (USA domestic) +1 800 535 5053 or international +1 352 323 3500

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Classification according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Skin Corr. 1B	H314
Eye Dam. 1	H318
STOT SE 3	H335

2.2. Label elements

Labelling according to OSHA Hazard Communication Standard 29 CFR 1910:1200

Hazard pictograms



Signal word

Danger

Trade name: Miele ProCare Dent 11 A

Version: 3 / US

Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

Hazard statements

H314 Causes severe skin burns and eye damage.
H335 May cause respiratory irritation.

Precautionary statements**Prevention**

P260 Do not breathe dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER or doctor.

Storage/Disposal

Dispose only when container is empty and closed. For disposal of product residues, refer to safety data sheet.

Hazardous component(s) to be indicated on label

contains disodium metasilicate

2.3. Other hazards**Supplemental information**

No special hazards have to be mentioned. The product contains no PBT or vPvB substances.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Hazardous ingredients according to OSHA Hazard Communication Standard 29 CFR 1910:1200****disodium metasilicate pentahydrate**

CAS No. 10213-79-3
Concentration >= 25 < 50 %

fatty alcohol, ethoxylated

CAS No. 146340-16-1
Concentration >= 1 < 10 %

Other information

Complete text of hazard statements in chapter 16

SECTION 4: First aid measures**4.1. Description of first aid measures****General information**

Remove contaminated clothing immediately and dispose of safely. In any case show the physician the Safety Data Sheet.

After inhalation

Ensure supply of fresh air. When dust is intensively inhaled, seek medical help immediately.

Trade name: Miele ProCare Dent 11 A

Version: 3 / US Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

After skin contact

Wash off immediately with soap and water. Take medical treatment.

After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). Summon a doctor immediately.

After ingestion

If swallowed, seek medical advice immediately and show this container or label. Rinse mouth thoroughly with water. Let plenty of water be drunk in small gulps. Do not induce vomiting.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed**Hints for the physician / hazards**

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Product itself is non-combustible; adapt fire extinguishing measures to surrounding areas.

Non suitable extinguishing media

Compatible with all usual extinguishing media.

5.2. Special hazards arising from the substance or mixture

In case of combustion evolution of dangerous gases possible.

5.3. Advice for firefighters**Special protective equipment for fire-fighting**

Do not inhale explosion and/or combustion gases. In case of combustion use a suitable breathing apparatus.

Other information

Collect contaminated fire-fighting water separately, must not be discharged into the drains. Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures**

Avoid contact with skin, eyes and clothing.

6.2. Environmental precautions

Do not discharge into the drains/surface waters/groundwater. Knock down dust with water spray jet.

6.3. Methods and material for containment and cleaning up

Pick up mechanically. Dispose of absorbed material in accordance with the regulations.

6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

Trade name: Miele ProCare Dent 11 A

Version: 3 / US Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Avoid the formation and deposition of dust. Keep container tightly closed.

7.2. Conditions for safe storage, including any incompatibilities

Recommended storage temperature

Value > 0 < 30 °C

Requirements for storage rooms and vessels

Keep in original packaging, tightly closed. Provide solvent-resistant and impermeable floor.

Storage class according to TRGS 510

Storage class according to TRGS 510 8B Non-combustible corrosive hazardous substances

7.3. Specific end use(s)

no data

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Other information

There are not known any further control parameters.

8.2. Exposure controls

General protective and hygiene measures

Do not inhale dust/fumes/aerosols. Avoid contact with skin and eyes. Do not eat, drink or smoke during work time. Wash hands before breaks and after work. Clean skin thoroughly after work; apply skin cream.

Respiratory protection

Use breathing apparatus in dust-laden atmosphere. Particle filter P2

Hand protection

Chemical resistant gloves

Use Permanent hand contact

Appropriate Material neoprene

Material thickness >= 0,65

Breakthrough time > 480

Appropriate Material nitrile

Material thickness >= 0,4

Breakthrough time > 480

Appropriate Material butyl

Material thickness >= 0,7

Breakthrough time > 480

Use Short-term hand contact

Appropriate Material nitrile

Material thickness >= 0,11

Hand protection must comply with EN 374.

Eye protection

Safety glasses with side protection shield; Eye protection must comply with EN 166.

Trade name: Miele ProCare Dent 11 A

Version: 3 / US

Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Form	solid
Colour	white
Odour	characteristic
Odour threshold	
Remarks	not determined
pH value	
Value	appr. 12
Concentration/H ₂ O	10 %
Temperature	20 °C
Melting point	
Remarks	not determined
Freezing point	
Remarks	not determined
Initial boiling point and boiling range	
Remarks	not determined
Flash point	
Remarks	Not applicable
Evaporation rate (ether = 1) :	
Remarks	not determined
Flammability (solid, gas)	
not determined	
Upper/lower flammability or explosive limits	
Remarks	Not applicable
Vapour pressure	
Remarks	not determined
Vapour density	
Remarks	not determined
Density	
Remarks	not determined
Solubility in water	
Remarks	soluble
Solubility(ies)	
Remarks	not determined
Partition coefficient: n-octanol/water	
Remarks	not determined
Ignition temperature	
Remarks	Not applicable
Decomposition temperature	
Remarks	not determined

Trade name: Miele ProCare Dent 11 A

Version: 3 / US Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

Viscosity

Remarks Not applicable

Explosive properties

evaluation no

Oxidising properties

evaluation None known

9.2. Other information**Bulk density**Value 1015 to 1065 kg/m³**Other information**

None known

SECTION 10: Stability and reactivity**10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

No hazardous reactions known.

Decomposition temperature

Remarks not determined

10.5. Incompatible materials

Strong exothermic reaction with acids.

10.6. Hazardous decomposition products

No hazardous decomposition products known.

SECTION 11: Toxicological information**11.1. Information on toxicological effects****Acute oral toxicity**

Species	rat		
ATE	>	2000	mg/kg
Method	calculated value (Regulation (EC) No. 1272/2008)		
Remarks	Based on available data, the classification criteria are not met.		

Acute oral toxicity (Components)**fatty alcohol, ethoxylated**

Species	rat		
LD50	>	2000	mg/kg

disodium metasilicate pentahydrate

Species	rat		
LD50		1150 to 1350	mg/kg

Trade name: Miele ProCare Dent 11 A

Version: 3 / US

Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

Acute dermal toxicity

Remarks Based on available data, the classification criteria are not met.

Acute inhalational toxicity

Remarks Based on available data, the classification criteria are not met.

Skin corrosion/irritation

evaluation corrosive

Remarks The classification criteria are met.

Serious eye damage/irritation

evaluation corrosive

Remarks The classification criteria are met.

Sensitization

Remarks Based on available data, the classification criteria are not met.

Subacute, subchronic, chronic toxicity

Remarks Based on available data, the classification criteria are not met.

Mutagenicity

Remarks Based on available data, the classification criteria are not met.

Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

National Toxicology Program (NTP)

No component of this product present at levels greater than or equal to 0,1% is identified as a known or reasonably anticipated to be human carcinogen by NTP

International Agency for research on Cancer (IARC)

No component of this product present at levels greater than or equal to 0,1% is identified as probable, possible or confirmed human carcinogen by IARC

Specific Target Organ Toxicity (STOT)**Single exposure**

Remarks The classification criteria are met.

evaluation May cause respiratory irritation.

Repeated exposure

Remarks Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Experience in practice

Inhalation of dusts may irritate the respiratory tract.

Other information

There is no data available on the product apart from the information given in this subsection.

SECTION 12: Ecological information**12.1. Toxicity****General information**

not determined

Trade name: Miele ProCare Dent 11 A

Version: 3 / US

Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

Fish toxicity (Components)**fatty alcohol, ethoxylated**

Species	golden orfe (Leuciscus idus)	
LC50	0,6	mg/l
Method	DIN 38412 / Part 15	

disodium metasilicate pentahydrate

Species	zebra fish (Brachydanio rerio)	
LC50	210	mg/l
Duration of exposure	96	h

Daphnia toxicity (Components)**fatty alcohol, ethoxylated**

LC50	1,2	mg/l
Method	DIN 38412 / Part 11	

disodium metasilicate pentahydrate

Species	Daphnia magna	
EC50	1700	mg/l
Duration of exposure	48	h

12.2. Persistence and degradability**General information**

not determined

Ready degradability (Components)**fatty alcohol, ethoxylated****12.3. Bioaccumulative potential****General information**

not determined

Partition coefficient: n-octanol/water

Remarks not determined

12.4. Mobility in soil**General information**

not determined

12.5. Results of PBT and vPvB assessment**General information**

not determined

Evaluation of persistence and bioaccumulation potential

The product contains no PBT or vPvB substances.

12.6. Other adverse effects**General information**

not determined

General information / ecology

The surfactant(s) contained in this preparation complies(comply) with the biodegradability criteria as laid down in Regulation (EC) No.648/2004 on detergents.

Trade name: Miele ProCare Dent 11 A

Version: 3 / US

Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

SECTION 13: Disposal considerations

13.1. Waste treatment methods




Disposal recommendations for the product

Allocation of a waste code number, according to the European Waste Catalogue (EWC), should be carried out in agreement with the regional waste disposal company.

Disposal recommendations for packaging

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

SECTION 14: Transport information

	Ground transport DOT	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
Tunnel restriction code	E		
IMDG-Code segregation group		18 Alkalis	
14.1. UN number	3253	3253	3253
14.2. UN proper shipping name	DISODIUM TRIOXOSILICATE	DISODIUM TRIOXOSILICATE	DISODIUM TRIOXOSILICATE
14.3. Transport hazard class(es)	8	8	8
Label			
14.4. Packing group	III	III	III
Limited Quantity	5 kg		
Transport category	3		
14.5. Environmental hazards		no	

Information for all modes of transport

14.6. Special precautions for user

See Sections 6 to 8

Other information

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance

Trade name: Miele ProCare Dent 11 A

Version: 3 / US

Replaces Version: 2 / US Date revised: 12.01.2022

Print date: 02.03.22

or mixture**Other information**

The product does not contain substances of very high concern (SVHC).

US. EPA Emergency Planning and Community Right-To-Know Act (EPCRA) SARA Title III Section 302 Extremely Hazardous Substance (40 CFR 355)

The product does not contain any listed components.

US. EPA Emergency Planning and Community Right-To-Know Act (EPCRA) SARA Title III Section 313 Toxic Chemicals (40 CFR 372.65) - Supplier Notification Required

The product does not contain any listed components.

California Safe Drinking Water and Toxic Enforcement Act (Proposition 65)

Warning! This product may contain trace quantities of substance(s) known to the state of California to cause cancer and/or reproductive toxicity – not added as a part of the formulation but remaining as residuals from the manufacturing process of our raw material suppliers.

SECTION 16: Other information**Abbreviations**

ADR: Accord européen relatif au transport international des marchandises Dangereuses par Route

RID: Règlement concernant le transport international ferroviaire de marchandises dangereuses

IMDG: International Maritime Code for Dangerous Goods

ICAO: International Civil Aviation Organization

IATA: International Air Transport Association

CAS: Chemical Abstracts Service

VOC: Volatile Organic Compound

LD: Lethal dose

LC: Lethal concentration

PBT: Persistent, Bioaccumulative and Toxic

vPvB: Very persistent and very bioaccumulative

SVHC: Substances of very high concern

MARPOL 73/78: International Convention for the Prevention of Pollution From Ships, 1973 as modified by the Protocol of 1978 (MARPOL: Marine Pollution)

IBC: Intermediate Bulk Container

OEL: Occupational exposure limit

TSCA: Toxic Substances Control Act (USA)

IMO: International Maritime Organization

GHS: Globally Harmonized System of classification and Labelling of Chemicals

REACH: Registration, Evaluation, Autohorisation and Restriction of Chemicals

UN: United Nations

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: ***

This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.