

PWA 8672 [U DS5 2AC208V/60HZ]

Washer

With liquid dispensing and AutoOpen function.



- Shortest cycle time: 25 min
- Capacity/cycle: up to 24 medium sized cassettes
- Performance/cycle, e.g., 44 handpieces or 44 suction devices
- Maximum user convenience with M Touch Flex full-touch colour display
- Connectivity as standard via Wi-Fi

EAN: 4002516884781 / Material number: 12800830 /
Old Material Number: 62867202NAM

Construction type and design

| | |
|-------------------|--|
| Construction type | Freestanding appliance, can be built-under, slim |
|-------------------|--|

| | |
|--------------|-------------|
| Product line | ExploreLine |
|--------------|-------------|

| | |
|----------------------|-------|
| Control panel colour | Black |
|----------------------|-------|

| | |
|--------------|-----------------|
| Outer casing | Stainless steel |
|--------------|-----------------|

| | |
|--------------------|---|
| Can be built-under | • |
|--------------------|---|

| | |
|-----|-------------|
| Lid | Without lid |
|-----|-------------|

| | |
|-----------------------|---|
| Number of wash levels | 2 |
|-----------------------|---|

| | |
|---------------------------------|---|
| AutoClose – automatic door lock | • |
|---------------------------------|---|

| | |
|---|---|
| Buzzer, acoustic signal at end of program | • |
|---|---|

| | |
|------------------------------|---|
| Easy-to-service construction | • |
|------------------------------|---|

Application

| | |
|-------------------------------|---|
| Suitable for dental practices | • |
|-------------------------------|---|

Capacity

| | |
|---|----|
| Surgical saliva extractors per cycle [number] | 44 |
|---|----|

| | |
|--|----|
| Lumened instruments per cycle [number] | 44 |
|--|----|

| | |
|-------------------------------|----|
| Handpieces per cycle [number] | 44 |
|-------------------------------|----|

| | |
|---|----|
| Mesh trays (176 x 273 x 30 mm) (W x H x D) per cycle [number] | 18 |
|---|----|

| | |
|---|----|
| Mesh trays (176 x 273 x 41 mm) (W x H x D) per cycle [number] | 18 |
|---|----|

Performance data

| | |
|--|-------------|
| Circulation pump, Qmax in gal/min (l/min.) | 132.1 (500) |
|--|-------------|

| | |
|---------------------------------|-----------------------------------|
| Usable wash cabinet volume in l | 5.12 (145.00) ft ³ (l) |
|---------------------------------|-----------------------------------|

| | |
|--|----|
| Shortest program duration (cleaning & disinfection) in minutes | 25 |
|--|----|

| | |
|------------------------|-------|
| Tested operating hours | 15000 |
|------------------------|-------|

| | |
|--|----------|
| Maximum final rinse temperature in °F (°C) | 158 (70) |
|--|----------|

Miele MOVE

| | |
|-----------------------------|---|
| Miele MOVE Connect (option) | • |
|-----------------------------|---|

Controls

| | |
|--------------|--------------|
| Control type | M Touch Flex |
|--------------|--------------|

| | |
|---|---|
| Cycle approval directly via the display | • |
|---|---|

| | |
|-----------------------|----|
| Max. delay start in h | 24 |
|-----------------------|----|

| | |
|-----------------|--------------|
| Programmability | programmable |
|-----------------|--------------|

| | |
|-------------------|---|
| Programs [number] | 5 |
|-------------------|---|

| | |
|-----------------------------|---|
| Free program slots [number] | 5 |
|-----------------------------|---|

| | |
|-------------------|---|
| Time left display | • |
|-------------------|---|

| | |
|----------------------------|---|
| Program sequence indicator | • |
|----------------------------|---|

| | |
|----------------------------|---|
| Choice of display language | • |
|----------------------------|---|

Standard electrical connection

| | |
|-----------------------|----------------|
| Electrical connection | 2 AC 208V 60HZ |
|-----------------------|----------------|

| | |
|------------------------------|------|
| Electric heater rating in kW | 5.40 |
|------------------------------|------|

| | |
|------------------------|------|
| Total rated load in kW | 6.20 |
|------------------------|------|

| | |
|------------------|----|
| Fuse rating in A | 30 |
|------------------|----|

| | |
|---------------------------|----------|
| Length of power cord in m | 6 (1.90) |
|---------------------------|----------|

Possible voltage variant 1

| | |
|-----------------------|----------------|
| Electrical connection | 3 AC 208V 60HZ |
|-----------------------|----------------|

PWA 8672 [U DS5 2AC208V/60HZ]

Washer

With liquid dispensing and AutoOpen function.



EAN: 4002516884781 / Material number: 12800830 /
Old Material Number: 62867202NAM

| | |
|---|---------------------|
| Electric heater rating in kW | 5.40 |
| Total rated load in kW | 6.20 |
| Fuse rating in A | 20 |
| Water connection/drainage | |
| Cold water [number] | 1 |
| Hot water [number] | 1 |
| Required flow pressure in kPa | 200-1000 |
| Maximum drain pump height (delivery head) in cm | 39 3/8 (100) |
| Maximum water hardness (cold water / hot water) in mmol/l | 62.59 (10.70) |
| Drain pump [DN] | 0.86 (22) |
| Integrated water softener | • |
| Water protection system | Waterproof system |
| Dimensions and weight | |
| External dimensions, net height in mm | 32 5/16 (820) |
| External dimensions, net width in mm | 23 9/16 (598) |
| External dimensions, net depth in mm | 23 9/16 (598) |
| Packaging dimensions, gross height in in. (mm) | 36 1/4 (920) |
| Packaging dimensions, gross width in in. (mm) | 26 7/16 (670) |
| Packaging dimensions, gross depth in in. (mm) | 29 3/16 (740) |
| Wash cabinet, usable height in mm | 20 1/2 (520) |
| Wash cabinet, usable width in mm | 20 7/8 (530) |
| Wash cabinet, depth of upper basket in in. (mm) | 18 11/16 (474) |
| Net weight in kg | 149.91 (68) |
| Wash cabinet, depth of lower basket in in. (mm) | 20 1/2 (520) |
| Gross weight in lb (kg) | 169.75 (77) |
| Maximum floor load in N/m ² | 1220 |
| Emission values | |
| Emission sound pressure level in the workplace | ≤70 dB(A) re 20 µPa |
| Heat dissipation into the room in MJ/h | 0.50 |
| Programs | |
| Vario Dental | • |
| Vario Dental + | • |
| Cold water rinse | • |
| Clean machine | • |
| Features | |
| 1 door dispenser for rinse aid | • |
| 1 dispenser pump for liquid cleaning agent | • |
| 1 dispenser pump for neutralizing agent | • |
| Suction lance for media canisters 10–20 L (optional) | • |
| Convenient suction lance DC5 | • |
| Volume flow control | • |
| Spray pressure monitoring | • |
| Spray arm monitoring | • |
| How-to sequences | • |
| Interface for process documentation (optional) | • |
| Multi-component filter system | • |
| Variable-speed heating pump | • |
| Load carrier – direct coupling | • |

PWA 8672 [U DS5 2AC208V/60HZ]

Washer

With liquid dispensing and AutoOpen function.



EAN: 4002516884781 / Material number: 12800830 /
Old Material Number: 62867202NAM

| | |
|--|---|
| AutoOpen | • |
| In-door salt container | • |
| Miele Service optical interface | • |
| Wash cabinet seamlessly welded using laser technology | • |
| Heater elements outside wash cabinet | • |
| Wash cabinet made of high-quality stainless steel (1.4301 / 304) | • |
| Extensive components (option) | • |
| Connection options | |
| Wi-Fi | • |
| Ethernet interface (optional) | • |
| Standards / certs / symbols | |
| IP safety class according to EN 60529: IP 21 | • |
| EN 61010-1 | • |
| EN 61010-2-040 | • |
| EU medical device class IIb | • |
| Health Canada, Medical Device class II | • |