

## PLW 8683 CD [U DC5 LAN]

Laboratory dishwasher with liquid dispensing & hot air drying DryPlus.



- Perf./cycle, e.g. 128 laboratory flasks
- Flexible and simple modular load carriers & EasyLoad system
- Maximum user convenience with M Touch Flex full-touch colour display
- Safety through monitoring wash pressure and spray arm monitoring
- Efficient hot air drying DryPlus
- Connectivity as standard via WiFi or Ethernet

EAN: 4002516884187 / Material number: 12799620 / Old Material Number: 6286C312NWE

Old Material Number: 6286C312NWE	
Construction type and version	
Construction type	Freestanding, can be built- under, wide, with lid
Range	ExpertLine
Control panel colour	Black
Casing	Stainless steel
Glass door	•
Lid	Without lid
Can be built-under	•
Number of wash levels	2
AutoClose – automatic door lock	•
Buzzer, acoustic signal at end of programme	•
Service-friendly design	•
Application	
Suitable for laboratories	•
Capacity	
Narrow-necked glassware per cycle [number]	64
Test tubes per cycle [number]	196
Vials per cycle [number]	196
Laboratory flasks per cycle [number]	128
Pipettes per cycle [number]	98
Performance data	
	500
Circulation pump, Qmax in I/min	145.00
Wash cabinet usable capacity in I	15000
Tested operating hours	<del></del>
Maximum final rinse temperature in °C	93
Miele MOVE	
Miele MOVE Connect (option)	·
Control system	NAT
Control type	M Touch Flex
Max. delay start in h	24
Programmability	programmable
Programmes [number]	20
Free programme positions [number]	10
Time left display	<u> </u>
Programme sequence display	<u> </u>
Selectable display languages	•
Standard electrical connection	
Electrical connection	3N AC 400V 50HZ
Electric heater rating in kW	8.50
Total rated load in kW	9.30
Fuse rating in A	16
Length of mains connection cable in m	1.90
Possible voltage variant 1	
Electrical connection	AC 230V 50HZ
Electric heater rating in kW	2.50
Total rated load in kW	3.30
Fuse rating in A	16



## PLW 8683 CD [U DC5 LAN]

Laboratory dishwasher with liquid dispensing & hot air drying DryPlus.



EAN: 4002516884187 / Material number: 12799620 / Old Material Number: 6286C312NWE

Old Material Number: 6286C312NWE	
Water connection/drainage	
Cold water [number]	_1
Hot water [number]	_ 1
Required flow pressure in kPa	200-1000
DI water [number]	1
Maximum delivery head for drain pump in cm	100
Maximum water hardness (cold water/hot water) in mmol/l	10.70
VEP feed pump for demineralised water (option)	•
Drain pump [DN]	22
Integrated water softener	•
Water protection system	Waterproof system
Dimensions and weight	
External dimensions, net height in mm	820
External dimensions, net width in mm	898
External dimensions, net depth in mm	700
Wash cabinet, usable height in mm	520
Wash cabinet, usable width in mm	530
Net weight in kg	94
Drying unit	
Integrated drying unit	•
Fan in kW	0.30
Heater element in kW	2.20
Air throughput in m³/h	60
Time setting in 1-minute increments	0-120
HEPA filter class	H14
HEPA filter filtration efficiency (in accordance with DIN EN 1822) in %	99.99
HEPA filter service life in h	200
Dimensions and weight	200
External dimensions, gross height in mm	1095
External dimensions, gross width in mm	1200
External dimensions, gross depth in mm	800
Wash cabinet, depth of upper basket in mm	474
Wash cabinet, depth of lower basket in mm	520
<del></del>	127
Gross weight in kg	
Maximum floor load in N/m <sup>2</sup>	1220
Emissions data	-70 dD(A) == 20D=
Emission sound pressure level in the workplace	<70 dB(A) re 20 μPa
Heat dissipation into the room in MJ/h	1.44
Programmes	
Standard	-
Intensive	•
Programmes	
Extra short	•
Universal	•
Inorganic	•
Organic	
Agar	•
Injector Plus	•



**Connection options** WiFi/Ethernet interface

VDE (electrical safety)

RoHS guideline 2011/65/EU

EN 61010-1 EN 61010-2-040

Standards, test certificates and symbols

IP safety class according to EN 60529: IP 21

EAN: 4002516884187 / Material number: 12799620 /

## PLW 8683 CD [U DC5 LAN] Laboratory dishwasher

with liquid dispensing & hot air drying DryPlus.



Old Material Number: 6286C312NWE Pipettes Plastics Oil Hygiene 93/10 Pasteurising Cold water rinse Demineralised water rinse **Features** 1 dispenser pump for liquid cleaning agent 1 dispenser pump for neutralising agent Additional integrated dispensing pump 5–10 I (optional) Suction lance for media canisters 10–20 I (optional) Convenient suction lance DC5 Additional integrated dispenser pumps (optional) Volume flow control Steam condenser Spray pressure monitoring Spray arm monitoring How-to sequences Wash cabinet lighting Multi-component filter system Variable-speed heater pump Load carrier — direct coupling DryPlus Salt container in base of wash cabinet Customer Service optical interface Setting options for process chemicals (3 x 5 l / 2 x 10 l) Wash cabinet seamlessly welded using laser technology Heater elements outside wash cabinet EasyLoad Wash cabinet made from high-quality stainless steel (1.4301/304) Extensive components (option)