

ProCare Dent 30 P - 5 l
Neutralizing agent, acidic, 5 l
for the automatic reprocessing of dental instruments.



- Suitable for dental instruments, not for transmission instruments
- Adapted to the process sequences in Miele washer-disinfectors
- Removes discoloration and soiling which is soluble in acid
- Neutralizes alkali residues that are carried over
- Based on phosphoric acid

EAN: 4029577158231 / Material number: 12796380 / Old Material Number: 66914303V6

Product category	
Thermal disinfectors	•
Large-chamber washer-disinfectors	•
Washer-disinfectors, dental	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing	
Recommended dosage per wash cycle	1 ml/l (0.1%)
Large opening for siphon	•
Automatic dispensing	•
Suitable for all levels of water hardness	•
Main applications	
Dental practices	•
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Phosphoric acid, Citric acid
pH value at 68°F (20°C)	0,7
Suitable for porcelain/ceramics	•
Suitable for standard instruments	•
Not suitable for anodized aluminum	•
Not suitable for light and non-ferrous metals	•
Not suitable for chrome-plated parts	•
Not suitable for nickel-plated parts	•
Not suitable for transmission instruments	•
Does not contain surfactants	•
Free from microplastics	•
Excellent rinsing properties	•
Prevention of limescale formation in wash cabinet	•
Neutralization of alkaline carryover	•
Protects against corrosion	•
Storage information	
Minimum storage temperature in °F/°C	-15
Maximum storage temperature in °F/°C	30
Dimensions and weight	
External dimensions, net height in inches	8 9/16 (217)
External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)
Gross weight in lbs (kg)	17 (7.80)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
CE	•
GHS/CLP compliant	•