

## A safe and perfect fit

### Compatibility overview: plaque removers, ultrasound scalers and air scaler tips

**Adapter A 814/1**  
Mat. no. 11691870  
Art. no. 68.6814.01



Note: Only fit A 814/1 on rails with central filter



- ✓ Ultrasound (PURE NEWTRON)
- ✗ IP2R, IP3L, TK21L, TK21R (not to be used)



Vector scaler



- ✓ Surgery Systems/ Preservative measures / Cavity preparation
- ✓ Scaling & Perio
- ✓ Endodontics



- ✓ SONICflex 2000, 2003/L
- ✓ SONICflex quick 2008/L
- ✓ SONOsoft, PiezoLUX
- ✓ PiezoLED, PiezoSoft



- ✓ Ultrasound (Piezo scaler)
- ✓ PiezoLine, EM1, SI1, KA1
- ✓ Schall (Air Scaler)
- ✓ SonicLine



- ✓ Varios
- ✓ Ti-Max S970, AS2000
- ✓ VarioSurg3



- ✓ SIROSONIC/ L/ TL
- ✓ PerioSonic



- ✓ Proxeo (air scaler)
- ✓ Tigon +, Tigon, Pyon 2 (PiezoScaler)

### Compatibility overview: cannulae from 3F/MF tips

**A 865 adapter**  
Mat. no. 11709860  
Art. no. 69.6865.00



Note:  
Only use in combination with A 800 tubular filter



Should you not find a particular tip in the overview, please contact the Miele hotline on 0800 22 44 644.