


Miele

Move your
Practice.



Miele
MOVE

Miele MOVE. Your digital solution
for efficient machine, data and
process management.

Efficient machine, data and process management in doctors' and dentists' surgeries.

Miele MOVE Digitalisation for Practices

Intuitive, transparent, efficient:

Miele MOVE connects washer-disinfectors/thermal disinfectors and sterilisers to continuously improve instrument reprocessing in your practice. Digital solutions combine convenient user interfaces with useful functions such as process documentation and direct service support.

Make your next move!



Manage
complex machine data
without assistance



Optimise
work processes to spend
more time with patients



Visualise
and approve process data



Enable
greater process security
and independence



Service optimisation

Interfaces & Reports

Miele MOVE will inform you **automatically** of any irregularities. Forward reports to Miele Customer Service or an authorised service agent with ease. This gives service technicians access to all important information and ensures that **requisite spare parts** are on board the service vehicle to ensure **more efficient and faster maintenance and repair processes**. This **reduces service costs, saves time** and reduces downtimes to a minimum.



Short response times

Notifications

More time for your patients: Miele MOVE notifies you proactively and across all products of individually selectable events. You and your team are kept up to date on the status of your machines at all times and from anywhere and are able to respond immediately for example in the event of unforeseen occurrences or when programmes finish. This **streamlines processes** and **frees up capacities**.



Intuitive user interface

Dashboard & Widgets

Convenient, clear and safe in use: Miele MOVE offers an intuitive interface with a **customisable dashboard**, personalised widgets and easy-to-understand symbols. You and your team can see the status of all machines at a glance, have all **relevant data at your fingertips** and can respond fast to any issues that may arise. This promotes independence, increases **process security** and simplifies compliance with **hygiene regulations**.

MOVE your Practice.

Full proof of hygiene

Process Documentation & Protocols

The **highest of hygiene standards** apply to surgical and dental instruments. Miele MOVE provides support in monitoring and documenting programmes with detailed **data capture**. Machine and process data is documented and visualised in easily comprehensible reports and charts. This makes it easy to furnish proof of **compliance with prescribed process parameters** at programme level and collate, approve and export programme data to meet obligations to present proof in a medical environment.



MOVE Ready: Our Compatible Models

The Miele washer-disinfectors PG 8582, PG 8582 CD, PG 8592 and PG 8562 and the thermal disinfectors PG 8581 and PG 8591 are **MOVE ready** thanks to their integrated communication module.

The sterilisers PST 1710, PST 1720, PST 2210 and PST 2220 will follow at the end of 2023. Precise Miele MOVE availability dates depend on the country concerned.




Make your next move with Miele MOVE!



Miele

Miele & Cie. KG, 33332 Gütersloh, Germany
www.miele-professional.com



Miele
MOVE