



Miele

From the Experts in Hygiene.

For the Experts
in Healthcare.

ExpertLine – our new product line
for the reliable reprocessing of your instruments.



Put Safety
Into Practice.
In Every Practice.
With ExpertLine.

Miele




For a wide range of
specialties

Tested to

15,000
operating hours



Glass door with internal lighting*



Infection control with **reliable quality.**

The new ExpertLine machines achieve the highest level of hygiene thanks to their design and versatile programmes.

Spare parts available for up to

15 years


after machine production is discontinued.



Sustainable

thanks to a smart use of resources and custom programmes.

Intuitive and **fast**  **Touch-Control**



**Drying options
suitable for
every application:**

from the AutoOpen function to DryPlus hot-air drying, our Experts offer various options.

Digital connectivity thanks to



**Miele MOVE
Connect**

**Excellent and effective washing and disinfection
performance with high capacity.**

The rinsing pressure is automatically adapted to different cleaning requirements – from the removal of stubborn dirt to the uniform rinsing of reprocessed items.



Discover ExpertLine
[www.miele.co.uk/
pro/pwd86Med](http://www.miele.co.uk/pro/pwd86Med)

* optional

Efficiency for Every Team.

Discover the efficiency of instrument reprocessing with ExpertLine and improve hygiene in your practice.

Intuitive operation

The intuitive colour touch display provides quick access to favorite programmes. It also allows you to adjust individual settings and display the programme status in a clear and user-friendly manner.

Perfect support

If required, how-to sequences on the display provide intuitive support for tasks such as changing filters or cleaning parts. This reduces downtime and ensures smooth operations.

Smart washing and disinfection

The unwelded wash cabinet ensures a hygienic interior, avoiding dirt build-up. Various sensors monitor the proper functioning of spray arms and the water pressure. Multi-stage filters for added safety and much more – ExpertLine offers expertise down to the smallest detail.

Keep an eye on all things

Thanks to the chamber lighting and optional glass door, the cleaning and disinfection process can be monitored even more easily and items can be loaded quickly and conveniently.





Miele

Podology



Flexibility for All Teams.

The modular basket concept offers full flexibility to meet your needs. Inserts and baskets developed specifically for various medical disciplines enable thorough disinfection of all types of instruments. ExpertLine also offers a range of options for drying, from the AutoOpen function to DryPlus hot-air drying.

ENT



Gynaecology



Ophthalmology



Anaesthesia



A Sustainable Addition to Your Team.

For us, durability is not just a word. That is why the ExpertLine is also tested for **15,000 operating hours**. Spare parts availability also ensures repairs **up to 15 years** after the series production is discontinued.

Programmes that are economical and designed for different applications, short running times and resource-saving technologies – all this allows you to save water, electricity and process chemicals and protects the environment.

An International and Digital Asset.

The simple language selection makes the ExpertLine models your international experts in the team. How-to sequences shown on the display ensure smooth processes.

Thanks to the connection with Miele MOVE Connect, machine data can be viewed anytime, from any device, and workflows can be optimised. You also have access to machine updates and receive evaluations of the cleaning process. For direct service support, there is also an interface that transmits important information to customer service so you can save time and effort.



Meet the Hygiene Experts.

All the experts offer maximum efficiency thanks to high capacity and resource-saving cleaning performance with the following highlights:



Large, intuitive
colour touch display




Liquid dispensing
with optional fill level detection



Comprehensive
programme package



Full connectivity: WLAN as standard,
LAN optional



Cabinet lighting for even more
convenient loading



Discover ExpertLine
[www.miele.co.uk/
pro/pwd86Med](http://www.miele.co.uk/pro/pwd86Med)



ExpertLine PWD 8682 with EcoDry drying

Our EcoDry drying has the AutoOpen function, which allows residual moisture to evaporate at the end of the programme, allowing instruments to cool faster and therefore be unloaded promptly.

ExpertLine PWD 8692 with DryPlus hot-air drying for the best drying results

Our DryPlus hot-air drying ensures the best drying results. Here the AutoOpen function ensures better removal of residual moisture and thus faster cooling of medical instruments.



ExpertLine PWD 8682 CD with additional storage space for cleaning agents and glass door option


This machine not only offers DryPlus hot-air drying with the AutoOpen function, but also convenient storage space for reprocessing chemicals agents in the side cabinet – here the glass door option allows you to easily monitor the cleaning and disinfection process.

Miele 360PRO.

The Comprehensive System Solution.

Precision is a responsibility. It is what drives us.

We know that hygiene plays an important role in the medical sector. To provide you with the best possible support, we offer 360PRO. Our comprehensive system with perfectly coordinated solutions from a single source. From planning to sustainability, from components to networking, from reprocessing chemicals to maintenance, including the optimal machines: we are there for you so that you can concentrate fully on your patients.



More about
Miele 360PRO
www.miele.co.uk/pro/med-dent

Miele

Sustainable Approach, Award-Winning Design.

Good for your team. Even better for the planet.

“Forever Better” has been our philosophy since 1899. It drives us to be better tomorrow than we were yesterday. It motivates us every day to develop better, smarter devices. ExpertLine also sets new standards in sustainability. For example by having spare parts available up to 15 years after the machine production is discontinued, with reduced operating costs and ecological footprint, and thanks to a longevity tested to 15,000 operating hours. ExpertLine is also equipped with a digital operating manual as a waste-free alternative.

Miele’s ExpertLine is the proud winner of the iF Design Award 2024.

“We are delighted to receive this award once again, which clearly highlights our premium standards of aesthetics and design.”

Andreas Enslin, Vice President Miele Design Center



So That You Are Fully Protected.

Our service is called PROtect and combines comprehensive services, maximum budget control and a smooth supply of spare parts. We provide you with transparent and industry-specific advice, coordinate appointments and guarantee the availability of suitable spare parts – even 15 years after machine production has been discontinued.



More about PROtect
www.miele.co.uk/pro/protect

Miele
PROtect

Nothing Beats **Personal Contact.**

Please feel free to contact us to find out more about ExpertLine – simply call us or use the contact form on our website.

Telephone: 0330 160 6693

Mon.-Fri. 09:00 - 17:15



More about Miele ExpertLine
and contact form

**[www.miele.co.uk/
pro/pwd86Med](http://www.miele.co.uk/pro/pwd86Med)**