



Miele

Choose Your Expert. Gain a Team Member.

Analytical cleanliness of your laboratory glassware —
with our ExpertLine.



Miele

For every industry
and specialisation

15,000
tested operating hours



Glass door & chamber lighting*



**Reliable quality
for perfect results**

Our ExpertLine laboratory washers ensure analytically pure laboratory glassware cleaning and drying at the highest level for continuously reliable results.

EASYLOAD

Flexible, fast and highly efficient

Spare parts available for up to

15 years

after machine production has been discontinued.



Sustainable

thanks to optimised programmes

**At last, a flexible
team member
who loves cleaning!**

Experience the
ExpertLine highlights
for yourself.

Intuitive and **fast**  **touch control**



**Suitable
drying options
for every application**

From the AutoOpen function to DryPlus hot air drying, our expert offers various options.

Digital networking thanks to



**Miele MOVE
Connect**

Outstanding and effective cleaning performance at high capacities

The wash pressure is adjusted in line with the relevant cleaning requirements – high pressure facilitates the removal of stubborn soiling, while low pressure ensures the laboratory glassware is wetted evenly while conserving resources.



More on
ExpertLine
www.miele.co.uk/pro/expertline-lab

*optional

Your Efficient Team Member.

Intuitive operation thanks to the colour touch display and economical operation, including for high capacities in a compact space.

Intuitive operation

The intuitive colour touch display provides quick access to your favourite programmes. Individual settings can be made and the programme status is clearly visible and easy to understand.

Ready for use quickly

Machine commissioning has been optimised to make it quick and easy. And should interruptions occur, ExpertLine is designed to ensure that services can be carried out quickly and easily thanks to the intelligent and optimised positioning of core components.

Effective cleaning performance

With the variable-speed heater pump, the wash pressure is adjusted in line with the relevant cleaning requirements – high pressure facilitates the removal of stubborn soiling, while low pressure ensures the laboratory glassware is wetted evenly while conserving resources. Various sensors monitor that the spray arms and the wash pressure are working correctly.

Perfect overview

Thanks to the optional glass door and chamber lighting, the cleaning process is even easier to monitor and laboratory glassware can be loaded quickly and conveniently.





Miele

Your Flexible Team Member.



To meet all the requirements of laboratory glassware cleaning, the flexible EasyLoad accessory system makes cleaning considerably faster, simpler and more efficient.

When it comes to drying, our expert offers options such as the AutoOpen function and DryPlus hot air drying.



Your Long-Lasting and Sustainable Team Member.

For us, longevity is not just a word. This is why ExpertLine is also tested for 15,000 operating hours and the availability of spare parts for up to 15 years ensures that repairs can be carried out long after the machine has been discontinued.

Economical programmes designed for the application in question, short running times and resource-saving technologies – all this allows you to conserve water, electricity and detergents and protects the environment.

Your International and Digital Team Member.

The option to easily select a language makes the ExpertLine models the international experts in your team. How-to sequences shown on the display provide assistance to ensure smooth processes.

Thanks to the connection with Miele MOVE Connect, machine data can be viewed from anywhere at any time and work processes can be optimised. You also have access to machine updates and receive analyses of the cleaning process. For direct service support, there is an interface to Customer Service that transmits important information to save you time and effort.



Choose Your Expert.

All experts offer maximum efficiency thanks to high capacities and resource-saving cleaning performance with the following highlights:

- 

Large, intuitive colour touch display
- 

Liquid dispensing with optional fill level detection
- 

Comprehensive programme package
- 

Full connectivity: WiFi as standard, LAN optional
- 

Chamber lighting for even greater convenience when loading
- 

Simply scan the QR code and discover ExpertLine.
www.miele.co.uk/pro/expertline-lab



ExpertLine PLW 8683 with EcoDry technology and glass door option

Our EcoDry technology features the AutoOpen function, which allows residual moisture to escape at the end of the programme, allowing glassware to cool down more quickly and therefore be unloaded promptly. The glass door variant also allows you to easily monitor the cleaning process.

ExpertLine PLW 8693 with DryPlus hot air drying for the best drying results

Our DryPlus hot air drying ensures the best drying results, even with narrow-volume glass. The AutoOpen function also ensures better removal of residual moisture and therefore faster cooling of the laboratory glassware.



ExpertLine PLW 8683 CD with additional storage space for cleaning agents and glass door option

In addition to DryPlus hot air drying with AutoOpen function, this machine also offers convenient storage space for cleaning agents in the side unit – here too, the glass door variant allows you to easily monitor the cleaning process.

Miele 360PRO.

A Tailored Solution for Every Laboratory.

Precision is a responsibility. One that inspires us.

We understand that precision is a responsibility for labs of all types. To support you in the best way, we offer 360PRO. Our holistic system with perfectly matched solutions from a single supplier. From planning to sustainability, accessories to connectivity, detergents and maintenance, with the right laboratory glasswashers: we are in service so that you can focus on delivering excellent results.



More on Miele 360PRO
www.miele.co.uk/pro/lab

Sustainable Aspiration,

Award-Winning Design.

Good for your team. Even better for the planet.

“Immer Besser” has been our philosophy since 1899. It motivates us to be better than we were yesterday. It is the motor which drives us to be better each day and to develop more intelligent machines. Our ExpertLine is setting new standards in terms of sustainability. For example, operating costs are reduced and the ecological footprint is minimised thanks to the 15-year availability of spare parts, while testing for 15,000 operating hours ensures an extremely long service life. Newly equipped with digital operating instructions as a waste-free alternative.

The Miele ExpertLine is the proud winner of the iF Design Award 2024.

“We are delighted to receive these awards once again, as they visibly demonstrate the premium quality of our aesthetics and design to the wider world.”

Andreas Enslin, Vice President Miele Design Centre



For All-Round

Protection.

Our PROtect service combines comprehensive services, maximum cost control and the easy availability of spare parts. We provide transparent advice specific to your industry, coordinate appointments and guarantee the availability of suitable spare parts – even 15 years after machines are discontinued.



More on PROtect
www.miele.co.uk/pro/protect

Nothing Beats **Personal Contact.**

Contact us to find out more about your new team member –
simply call us or use the contact form on our website.

Tel.: 0330 160 6693

Mon.–Fri. 9:00 am – 17:15 pm



More on Miele ExpertLine
and to the contact form
www.miele.co.uk/pro/expertline-lab