


Miele

ProCare Dent

Process Chemicals for
Dental Cleaning and Reprocessing

Miele Professional. Immer Besser.



ProCare Dent

Cleaning Agents



The all-rounder for daily use in practice

ProCare Dent 10 MA

Non-hazardous liquid cleaning agent, mildly alkaline

- For the automatic reprocessing of dental instruments and handpieces
- Perfectly matched to the processes and cleaning cycles of Miele Professional washer and washer-disinfectors
- Very gentle on materials
- Particularly suitable for sensitive materials
- Fast and thorough cleaning with enzymes
- Container size: 5 l
- Mat. nr.: 12680720

The specialist for stubborn cases

ProCare Dent 10 A

Liquid cleaning agent, alkaline

- For the automatic reprocessing of dental instruments
- Perfectly matched to the processes and cleaning cycles of Miele Professional washer and washer-disinfectors
- Very gentle on materials
- Corrosion protection for stainless steel instruments and trays
- Low residue – free from surfactants
- Safe to use – free from oxidative additives
- Container size: 5 l
- Mat.nr.: 12786490



ProCare Dent

Neutralizing Agents



Gentle neutralization for sensitive instruments

ProCare Dent 30 C **Neutralizing agent, acidic**

- For the automatic reprocessing of dental instruments and handpieces
- Perfectly matched to the processes and cleaning cycles of Miele Professional washer and washer-disinfectors
- Good material compatibility
- Environmentally friendly; does not contain phosphates
- Reliably removes residues from previous cleaning steps and neutralizes any remaining alkali residues
- Based on citric acid
- Packaging unit: 5 l
- Mat.nr.: 12796330



Powerful neutralization for standard instruments

ProCare Dent 30 P **Neutralizing agent, acidic**

- For the automatic reprocessing of dental instruments
- Perfectly matched to the processes and cleaning cycles of Miele Professional washer and washer-disinfectors
- Removes discolouration and soiling which is soluble in acid
- Reliably removes residues from previous cleaning steps and neutralizes any remaining alkali residues
- Based on phosphoric acid
- Packaging unit: 5 l
- Mat.nr.: 12796380

ProCare Dent

Rinsing Agent



Fast drying, residue-free results

ProCare Dent 40 **Rinsing agent, acidic**

- For the automatic reprocessing of dental instruments and handpieces
- Perfectly matched to the processes and cleaning cycles of Miele Professional washer and washer-disinfectors
- Improves and accelerates drying
- Reliably removes residues from previous cleaning steps and neutralizes any remaining alkali residues
- Container size: 1 l
- Mat.nr.: 12796420



For a safe finish after cleaning:

ProCare Dent neutralizing and rinsing agents ensure reliable removal of cleaning residues from the instruments. Whether gentle with citric acid or powerful with phosphoric acid, the products are precisely tailored to the requirements of the instruments in question. The rinsing agent rounds off the process and ensures quick and stain-free drying – for hygienically clean results that impress.

Additional products and accessories



ProCare Universal 61 Reactivation salt

- For optimum reactivation of integrated water softener
- Prevents limescale and protects load items
- Container size: 6 kg
- Mat.nr.: 9195780



ProClean Check Cleaning indicators (available in Canada only)

- Routine verification of cleaning cycles in automatic washer and washer-disinfectors
- Test for residues following cleaning or thermal disinfection
- Organic layer disappears upon effective cleaning
- 200 indicators per box
- Works with ProClean Check holder (recommended 2 per machine)
- Mat.nr. ProClean Check indicator: 97200392
- Mat.nr. ProClean Check holder: 97200393



CK/1 DIN 51 Canister key

- Easy opening of ProCare 5 l, 10 l and 20 l canisters
- Enables simple opening of machine-locked canisters
- Stable and long-lasting
- Can be used on both sides
- Mat.nr.: 11851600



DTD 2 Drip tray

- Securely accommodates 5 and 10 l ProCare canisters
- Offers space for 2 canisters
- The angle of incline ensures the residual emptying of canisters
- For storage in an adjacent unit
- Mat.nr.: 10695460



The right machines

Our Miele Professional dental instrument washer PWA 8672, and washer-disinfectors PWD 8682 and PWD 8692 (Canada only), are perfectly matched to ProCare Dent products. Whether with two or three dispensing pumps – we offer the right technology for your requirements.

For the reprocessing of standard instruments, handpieces and sedation hoses in a manner that is reliable and efficient.



PROtect: predictable service

- **Working for you:** Our Miele Professional service team regularly checks your machine – for safe processes and documented quality.
- **Testing made easy:** ProCare Dent and Miele Professional machine provide you with a solution that can be tested easily and reliably.





Miele Limited, Canada
www.mieleprofessional.ca

Miele Inc., USA
www.mieleusa.com/p