

E 130  
Insert  
for the optimum loading of 10 trays.



- equipped with 11 holders, height 170 mm, distance apart 35 mm
- For use in lower basket

EAN: 4002511129375 / Material number: 02648950 / Old Material Number: 69513001D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for trays	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	1 7/16 (35)
Rail height in mm	6 3/4 (170)
Application area	
Reprocessing of trays	•
Capacity	
Trays [quantity]	10
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	180
External dimensions, net width in mm	180
External dimensions, net depth in mm	445
External dimensions, gross height in mm	8 5/16 (210)
External dimensions, gross width in mm	8 5/16 (210)
External dimensions, gross depth in in (mm)	18 9/16 (470)
Net weight in kg	2 (0.76)
Gross weight in kg	2 (0.90)