


E 328

Support rail

for the optimum loading of upright instruments.



- The insert can hold upright instruments

EAN: 4002511117051 / Material number: 02795940 / Old Material Number: 69732801D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
Product category	
Other accessories	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0 5/8 (15) in (mm)
External dimensions, net width in mm	1 3/8 (34) in (mm)
External dimensions, net depth in mm	5 3/4 (145) in (mm)
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 3/16 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Net weight in kg	1 (0.06)
Gross weight in kg	1 (0.07)