

E 338

Insert

for the optimum loading of up to 8 half-shell trays.

EAN: 4002511116160 / Material number: 03080600 / Old Material Number: 69533801D



- Suitable for 8 half-shell trays
- equipped with 10 holders (8 compartments), W 295, D 33 mm
- Maximum tray size 290 x 30 mm
- For use in upper and lower basket

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for trays	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	8
Capacity	
Trays [quantity]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	4 3/8 (110) in (mm)
External dimensions, net width in mm	12 1/16 (305) in (mm)
External dimensions, net depth in mm	17 9/16 (445) in (mm)
External dimensions, gross height in mm	4 15/16 (125)
External dimensions, gross width in mm	12 1/4 (310)
External dimensions, gross depth in in (mm)	19 11/16 (500)
Net weight in kg	2 (0.62)
Gross weight in kg	2 (0.84)