

E 378

Insert

For holding different instruments.



- 0.8 mm wire mesh, 1.7 mm mesh width, 5 mm frames
- equipped with 2 handles
- For use in lower basket

EAN: 4002511116443 / Material number: 03203530 / Old Material Number: 69537801D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	4 3/8 (110) in (mm)
External dimensions, net width in mm	18 1/8 (460) in (mm)
External dimensions, net depth in mm	18 1/8 (460) in (mm)
External dimensions, gross height in mm	6 1/2 (165)
External dimensions, gross width in mm	21 1/8 (535)
External dimensions, gross depth in in (mm)	22 1/16 (560)
Net weight in kg	4 (1.78)
Gross weight in kg	5 (2.25)