

E 146

Insert

For the optimum loading of different instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 8 mm
- equipped with 2 hinged handles
- For use in upper and lower basket

EAN: 4002511116009 / Material number: 03346680 / Old Material Number: 69514601D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	55
External dimensions, net width in mm	150
External dimensions, net depth in mm	225
External dimensions, gross height in mm	3 (75)
External dimensions, gross width in mm	6 1/2 (165)
External dimensions, gross depth in in (mm)	9 3/4 (245)
Net weight in kg	2 (0.61)
Gross weight in kg	2 (0.70)