

SD-B
Injector nozzle
for butyrometers.



- For mobile injector units
- Including thread, 4 x 140 mm*

* Including thread, 4 x 140 mm, plus welded flattened nozzle, 1.5 x 100 mm.

EAN: 4002511116955 / Material number: 03583540 / Old Material Number: 69708001D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	240
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 15/16 (100)
External dimensions, gross depth in in (mm)	5 15/16 (150)
Diameter in in (mm)	0 3/16 (4)
Net weight in kg	1 (0.01)
Gross weight in kg	1 (0.02)