


O 175

Upper basket
for the optimum loading of laboratory glassware.



- Suitable for narrow-necked and laboratory glassware
- Vertical clearance: 6 3/4 (170) in (mm) mm*
- Equipped with connection for hot air drying
- Powder dispensing not possible

* Depending on the laboratory glassware the max. vertical clearance can be lower

EAN: 4002511114876 / Material number: 03807410 / Old Material Number: 69117501D

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Insert / module for narrow necked glassware	•
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	6 3/4 (170)
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	33
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	16 1/4 (412)
External dimensions, net width in in. (mm)	18 3/4 (475)
External dimensions, net depth in in. (mm)	20 15/16 (531)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in lb (kg)	14 (6.30)
Gross weight in lb (kg)	16 (7.20)
Accessory pack	
E 351 injector nozzle, Ø 4 mm, length 160 mm [number]	33
E 353 clip for injector nozzle E 351 [number]	33