


O 176

Upper basket

for the optimum loading of instruments for MIS, arthroscopy and urology.



- Equipped with 14 injector nozzles/funnels
- Left inserts, right hollow instruments (up to 500 mm)
- Vertical clearance left is 360 +/- 20 mm
- Powder dispensing not possible
- No spray arm, with connection for hot air drying

EAN: 4002511114890 / Material number: 03807420 / Old Material Number: 69117601D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Capacity	
Lumened instruments [number]	10
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	16 1/8 (409)
External dimensions, net width in in. (mm)	18 3/4 (475)
External dimensions, net depth in in. (mm)	20 15/16 (531)
External dimensions, gross height in in. (mm)	21 1/2 (545)
External dimensions, gross width in in. (mm)	15 3/4 (400)
External dimensions, gross depth in in (mm)	20 7/8 (530)
Net weight in lb (kg)	11 (4.80)
Gross weight in lb (kg)	15 (6.50)
Accessory pack	
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	2
E 448 silicone hose with female Luer Lock adapter [number]	2
E 452 injector nozzle 2.5 x 60 mm [number]	1
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	1
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	1
Spray jet, Ø 2.5 mm [number]	2
Spray jet, Ø 4 mm [number]	4
Funnel [number]	4