


O 183

Upper basket

for the optimum loading of instruments for MIS, arthroscopy and urology.



- Left side: 10 injector nozzles or free accommodation of inserts
- Right side: 14 injector nozzles/ funnels for hollow instruments
- Vertical clearance 285 +/- 30mm; hollow instruments 370 +/- 30mm
- Powder dispensing not possible

EAN: 4002511114920 / Material number: 03807470 / Old Material Number: 69118301D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Capacity	
Lumened instruments [number]	14
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	10 11/16 (270)
External dimensions, net width in in. (mm)	18 3/4 (475)
External dimensions, net depth in in. (mm)	20 7/8 (530)
External dimensions, gross height in in. (mm)	16 3/16 (410)
External dimensions, gross width in in. (mm)	21 5/16 (540)
External dimensions, gross depth in in (mm)	21 7/8 (555)
Net weight in lb (kg)	9 (4.00)
Gross weight in lb (kg)	11 (4.80)