

O 184

Upper basket

for optimal loading centrifuge tubes, vials or test tubes, for example.



- Equipped with 96 injector nozzles
- Equipped with connection for hot air drying
- Vertical clearance: 6 3/4 (170) in (mm) mm
- Powder dispensing not possible

EAN: 4002511114937 / Material number: 03807480 / Old Material Number: 69118401D

| Product category  |                          |
|---|--------------------------|
| Glassware Washers, laboratory                             | •                        |
| Machine affiliation                                       |                          |
| PG 8535   | •                        |
| PG 8536   | •                        |
| Product category  |                          |
| Basket for vials  | •                        |
| Basket for test tubes                                     | •                        |
| Injector basket   | •                        |
| Upper basket  | •                        |
| Product characteristics                                   |                          |
| Material  | Stainless steel, Plastic |
| Colour  | Stainless Steel, Black   |
| Maximum vertical clearance in mm                          | 6 3/4 (170)              |
| Application area  |                          |
| Reprocessing of pipettes / phials                         | •                        |
| Reprocessing of test tubes                                | •                        |
| Capacity  |                          |
| Test tubes [number]                                       | 96                       |
| Standard electrical connection                            |                          |
| Number of phases  | 0                        |
| Voltage in V  | 0                        |
| Frequency in Hz   | 0                        |
| Total rated load in kW                                    | 0.00                     |
| Fuse rating in A  | 0                        |
| Dimensions and weight                                     |                          |
| External dimensions, net height in in. (mm)               | 18 7/16 (468)            |
| External dimensions, net width in in. (mm)                | 18 3/4 (475)             |
| External dimensions, net depth in in. (mm)                | 20 15/16 (531)           |
| External dimensions, gross height in in. (mm)             | 10 11/16 (270)           |
| External dimensions, gross width in in. (mm)              | 20 11/16 (525)           |
| External dimensions, gross depth in in (mm)               | 21 5/16 (540)            |
| Net weight in lb (kg)                                     | 10 (4.50)                |
| Gross weight in lb (kg)                                   | 13 (5.70)                |
| Accessory pack  |                          |
| ID 110, injector nozzle, Ø 2.5 mm, length 110 mm [number] | 96                       |