

U 167

Lower basket

for the optimum loading of anaesthetic utensils.



- Suitable for 4 soda lime containers, 9 secretion flasks & various utensils
- 20 holders, height 7 7/8 (200) in (mm) mm, spacing approx. 3 3/4 (95) in (mm) mm

EAN: 4002511114784 / Material number: 03807600 / Old Material Number: 69116702D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Lower basket	•
Mobile unit for anaesthesia instruments	•
Product characteristics	
Material	Wire, rilsan-coated
Colour	White
Distance between rails in mm	3 3/4 (95)
Rail height in mm	7 7/8 (200)
Application area	
Reprocessing anaesthetic instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	8 11/16 (220)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	20 3/8 (516)
External dimensions, gross height in in. (mm)	9 1/2 (240)
External dimensions, gross width in in. (mm)	21 11/16 (550)
External dimensions, gross depth in in (mm)	21 11/16 (550)
Net weight in lb (kg)	7 (3.10)
Gross weight in lb (kg)	11 (4.60)