


E 106

Insert

for the optimum loading of wide neck flasks, measuring cylinders, etc.



- Equipped with 10 spring hooks x 175 mm height and 16 x 105 mm height
- The distance between the spring hooks is approx. 60 mm

EAN: 4002511115675 / Material number: 03808310 / Old Material Number: 69510601D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Insert / module for beakers	•
Insert / module for wide-necked glassware	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	7 1/8 (180) in (mm)
External dimensions, net width in mm	7 11/16 (195) in (mm)
External dimensions, net depth in mm	7 3/8 (186) in (mm)
External dimensions, gross height in mm	8 1/8 (205)
External dimensions, gross width in mm	9 11/16 (245)
External dimensions, gross depth in in (mm)	19 5/16 (490)
Net weight in kg	3 (0.93)
Gross weight in kg	4 (1.42)